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

Initiating Instrument Connection...

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About

MicroCal PEAQ-ITC Control Software

Malvern MicroCal PEAQ-ITC Control Software

Run Experiment Maintenance Design Experiment

Start Experiment Start Experiment Load Run Clean

Methods Experiments

Browse... C:\ > Users > Public > Documents > Malvern Instruments > MicroCal PEAQ-ITC > Methods

Name	Modified Date
13 Injections.itcm (MicroCal Method)	
19 Injections.itcm (MicroCal Method)	

Open

Preview - 13 Injections.itcm

This experimental method will make a single 0.4 μ L injection, followed by 12, 3 μ L injections.

Temperature ($^{\circ}$ C) 25.0
Reference Power (μ W) 41.9
Feedback High
Stir Speed (rpm) 750
Initial Delay (s) 60.0
Injection Spacing (s) 150
Injection Duration (s) 6

Instrument Offline

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The screenshot displays the MicroCal PEAQ-ITC Control Software interface. A red box highlights the main navigation bar at the top, which includes the following elements:

- 实验参数设置** (Experiment Parameter Settings): Callout pointing to the 'Run Experiment' button.
- 实验设计** (Experiment Design): Callout pointing to the 'Design Experiment' button.
- 设置** (Settings): Callout pointing to the gear icon in the top right corner.

Below the navigation bar, the interface is divided into several sections:

- Run**: A central control area with buttons for 'Start Experiment', 'Load', 'Run', and 'Clean'. A callout **仪器维护视频** (Instrument Maintenance Video) points to the 'Start Experiment' button.
- Experiment Information**: A sidebar on the left containing fields for [Syr] (M), [Cell] (M), and Comment.
- Instrument Settings**: A table of parameters to be used during the experiment:

Parameter	Value
Reference Power ($\mu\text{cal/s}$)	10.0
Feedback	High
Stir Speed (rpm)	750
Initial Delay (s)	60
Injection Spacing (s)	150
- Process Flow**: A horizontal timeline showing the current state: Idle (active), Setting Temperature, Equilibrating, and Injecting (0/13). A callout **字体大小设置** (Font Size Settings) points to the zoom control in the top right.
- Graph**: A plot of DP ($\mu\text{cal/s}$) versus Time (min) with a y-axis from -0.04 to 0.04 and an x-axis from -0.4 to 0.5.
- Buttons**: 'Start', 'Analyze', and 'Save As Method' buttons are located at the bottom of the graph area.

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The screenshot displays the MicroCal PEAQ-ITC Control Software interface, divided into three panels illustrating the workflow for opening and editing an experimental method.

Left Panel: Method Selection
The 'Methods' tab is selected in the top navigation bar. A list of methods is shown in a table. The 'Open' button at the bottom right of the list is circled in red. A blue callout box with the text '打开一个方法' (Open a method) has an arrow pointing to the 'Open' button.

Name	Modified Date
13 Injections.itcm (MicroCal Method)	
19 Injections.itcm (MicroCal Method)	
guojj.itcm	1/12/2016 3:51:33 PM
H2O.itcm	2/4/2016 10:27:36 AM
husc.itcm	2/17/2016 9:34:37 AM
noise.itcm	8/10/2015 12:05:13 AM
wyichen.itcm	8/31/2015 9:34:04 PM
Min.Yang	11/23/2015 6:11:17 PM

Middle Panel: Run Configuration
The 'Run' configuration window is shown. The 'Experiment Information' section contains fields for [Syr] (M) and [Cell] (M), both set to 0. A blue callout box with the text '修改实验参数' (Modify experimental parameters) has an arrow pointing to the edit icon (a pencil) in the 'Instrument Settings' section, which is circled in red.

Right Panel: Detailed Settings and Graph
The 'Run' configuration window is shown in more detail. The 'Instrument Settings' section is circled in red and contains the following values:
Temperature (°C): 25.0
Reference Power (µcal/s): 10.0
FeedBack: High
Stir Speed (rpm): 750
Initial Delay (s): 60
Injection Settings: # of Injections: 13
A table below shows the injection schedule:

Injection	Volume (µL)	Duration (s)	Spacing (s)
1	0.4	0.8	150
2	3.0	6.0	150
3	3.0	6.0	150
4	3.0	6.0	150
5	3.0	6.0	150
6	3.0	6.0	150
7	3.0	6.0	150
8	3.0	6.0	150

A blue callout box with the text '参数编辑' (Parameter editing) has an arrow pointing to the 'Instrument Settings' section. To the right, a graph shows DP (µcal/s) on the y-axis (ranging from -0.05 to 0.04) and a parameter on the x-axis (ranging from -0.4 to -0.3). The text '样品浓度' (Sample concentration) is written in a blue callout box above the graph.

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4003384toTD2-2.itc - MicroCal PEAQ-ITC Control Software


Malvern MicroCal PEAQ-ITC Control Software

Run Experiment Maintenance Design Experiment

Clean Start Experiment Load Run Clean

0 1 2 3 4 5 6

0 Introduction

To view a video depicting all the cleaning steps, click the Play button.
To enter the Clean Instrument workflow, click Next.

Back Next

- 1 Choose Cleaning Method(s)
- 2 Insert Cell Cleaning Tool
- 3 Attach Fill Port Adapter
- 4 Move Pipette to Clean Location
- 5 Detach Fill Port Adapter
- 6 Remove Cell Cleaning Tool

Show entire video in external player

Instrument Online DP (μcal/s): 4.95 Temperature (°C): 22.7

Stop Idle Temperature (°C): 25 Set Syringe: Load Plunger Down Cell Clean: None Syringe Clean: None Clean

Windows taskbar: ITC 4003384toTD2-2.itc ... Paint 3D 10:36 AM 4/5/2020

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
The screenshot displays the Malvern MicroCal PEAQ-ITC Control Software interface. The main window shows a 'Run' tab with a 'Start Experiment' button and a 'Run' button. A progress bar at the top indicates the current state: Idle, Setting Temperature, Equilibrating, Injecting (0/19), and Ready. The 'Injecting' stage is highlighted with a red box and labeled '实验进度' (Experiment Progress). Below the progress bar is a graph of DP ($\mu\text{cal/s}$) vs. Time (min), showing a series of peaks. A blue box labeled '分析数据' (Analyze Data) points to the 'Analyze' button. Another blue box labeled '另存方法' (Save Method) points to the 'Save As Method' button. A third blue box labeled '实时DP值和温度' (Real-time DP value and temperature) points to the 'DP ($\mu\text{cal/s}$): 4.96' and 'Temperature ($^{\circ}\text{C}$): 21.5' status bar. A fourth blue box labeled '③开始实验' (3 Start Experiment) points to the 'Start' button. A fifth blue box labeled '手动停止实验' (Manually stop experiment) points to the 'Stop' button in the 'Instrument Online' section. A sixth blue box labeled '设置温度' (Set temperature) points to the 'Set' button. A seventh blue box labeled '②滴定针上样' (2 Inject sample into titration needle) points to the 'Load' button. An eighth blue box labeled '①样品池/滴定针清洗' (1 Wash sample cell/titration needle) points to the 'Clean' button. The 'Instrument Online' section also shows 'Plunger Down', 'Cell Clean: None', and 'Syringe Clean: None' buttons. The bottom status bar shows '10:35 AM 4/5/2020' and a logo for Tsinghua University.

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New Analysis - MicroCal PEAQ-ITC Analysis Software

Malvern MicroCal PEAQ-ITC Analysis Software

Analyze Experiment(s) Design Experiment

Assign Controls

Start Analysis Overview Adjust Baseline **Assign Controls** Adjust Fit Presentation

Experiments

Sort by Bin

CaEDTAgGetStart *
Binding >

CaEDTAgGetStart_ctrl *
Control

ΔH (kcal/mol)

Molar Ratio

ΔQ (μ cal)

Injection Number

Control Parameters

Type Single

Titrant -> Buffer (-) CaEDTAgGetStart_ctrl Method Line

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Experiment Information

Filename CaEDTAgGetStart

Temperature ($^{\circ}$ C) 25.0

[Syr] (M) 1.00e-3

[Cell] (M) 100e-6

Ref. Power (μ cal/s) 5.00 (5.35)

Comment

Control Type Single

Control Experiment CaEDTAgGetStart_ctrl

Results

Bin Comment Binding

Competitive Model No

[Syr] (M) 1.00e-3

[Cell] (M) 100e-6

Ligand In Cell No

N (sites) 0.969 \pm 1.5e-3

Reset Show Excluded Injections Show Legends

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Malvern MicroCal PEAQ-ITC Analysis Software

Analyze Experiment(s) Design Experiment

Assign Controls

Start Analysis Overview Adjust Baseline Assign Controls Adjust Fit Presentation

Experiments

Sort by Bin

edta *
Binding

water *
Check data < change to Check data ▾

- Binding
- No binding
- Control
- Check data

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New Analysis - MicroCal PEAQ-ITC Analysis Software

Malvern MicroCal PEAQ-ITC Analysis Software

Analyze Experiment(s) Design Experiment

Adjust Baseline

Start Analysis Overview **Adjust Baseline** Assign Controls Adjust Fit Presentation

Experiments
Sort by Bin
CaEDTAGetStart *
Binding
CaEDTAGetStart_ctrl *
Control

DP vs Time and ΔH vs Molar Ratio

DP ($\mu\text{cal/s}$) vs Time (min)

Area (μcal) = -7.99

Lock Markers Previous Injection # 1 Next

Reset Time Factor 5 Points per Injection 25 Show Excluded Injections

Experiment Information

Filename	CaEDTAGetStart
Temperature ($^{\circ}\text{C}$)	25.0
[Syr] (M)	1.00e-3
[Cell] (M)	100e-6
Ref. Power ($\mu\text{cal/s}$)	5.00 (5.35)
Comment	
Control Type	Single
Control Experiment	CaEDTAGetStart_ctrl
Results	
Bin Comment	Binding
Competitive Model	No
[Syr] (M)	1.00e-3
[Cell] (M)	100e-6
Ligand In Cell	No
N (sites)	0.969 \pm 1.5e-3

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New Analysis - MicroCal PEAQ-ITC Analysis Software

Malvern MicroCal PEAQ-ITC Analysis Software

Analyze Experiment(s) Design Experiment

Overview Start Analysis Overview Adjust Baseline Assign Controls

Fitting Model One Set of Sites

One Set of Sites
Two Sets of Sites
Dissociation
Sequential Binding Sites
One Set of Sites - SIM
Enzyme Kinetics - Multiple Injections
Enzyme Kinetics - Single Injection

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DP (μ)

Time (min)

Display Normalized Heat Raw Heat

ΔH (kJ/mol)

Molar Ratio

Show Excluded Injections Show Legends

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New Analysis - MicroCal PEAQ-ITC Analysis Software

Malvern MicroCal PEAQ-ITC Analysis Software

Analyze Experiment(s) Design Experiment

Adjust Fit

Start Analysis Overview Adjust Baseline Assign Controls **Adjust Fit** Presentation

Experiments

Sort by Bin

CaEDTAgGetStart *
Binding >
CaEDTAgGetStart_ctrl *
Control

CaEDTAgGetStart

Experiment Information

Filename CaEDTAgGetStart

Temperature (°C) 25.0

[Syr] (M) 1.00e-3

[Cell] (M) 100e-6

Ref. Power (μcal/s) 5.00 (5.35)

Comment

Control Type Single

Control Experiment CaEDTAgGetStart_ctrl

Results

Bin Comment Binding

Competitive Model No

[Syr] (M) 1.00e-3

[Cell] (M) 100e-6

Ligand In Cell No

N (sites) 0.969 ± 1.5e-3

Fitting Parameters

Ligand is in Cell

Parameter	Vary	Initial Value	Lower Bound	Upper Bound
N (sites)	<input checked="" type="checkbox"/>	0.969	1.0e-3	10.0
Offset (kcal/mol)	<input type="checkbox"/>	0	-80.0	80.0
[Syr] (M)	<input type="checkbox"/>	1.00e-3	0	1.00

Use Competitive Model

Unknown Binder Strong

[Weak] (M) 100e-6

Known Weak Parameters Enter Manually

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New Analysis - MicroCal PEAQ-ITC Analysis Software

Malvern MicroCal PEAQ-ITC Analysis Software

Analyze Experiment(s) Design Experiment

Presentation

Start Analysis Overview Adjust Baseline Assign Controls Adjust Fit **Presentation**

Result Table **Final Figure** Scatter Plot Injection Table Statistics Plot Signature Plot Raw Plot Integrated Heat Plot

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



29.5 μ M Protein titrated with 1.1mM Compound

11.5 μ M Compound titrated with 179 μ M Protein

Parameter	Ligand in syringe	Ligand in cell
n	0.99	0.97
K_d	104 nM	105 nM
ΔG°	-9.4 kcal/mol	-9.4 kcal/mol
ΔH°	-11.9 kcal/mol	-12.4 kcal/mol
$T\Delta S^\circ$	-2.5 kcal/mol	-3.0 kcal/mol

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

Peptide binding to Batch 1 Protein

N = 1.05
 $K_D = 97 \text{ nM}$

100% of Batch 1
protein Active based
on Stoichiometry

Peptide binding to Batch 2 Protein

50 μM Peptide
10 μM Protein X
N = 0.235
 $K_D = 135 \text{ nM}$


Re-analyzed as
2.3 μM Protein X
N = 1.02
 $K_D = 135 \text{ nM}$

23% of Batch 2
protein Active based
on Stoichiometry

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