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### Fixation of Tissues in 10% Neutral Buffered Formalin

1. *Sacrifice animal by prescribed and approved euthanasia techniques. Tissues to be fixed and processed should be cut to a size no larger than 3 mm thick. Let tissues fix in 10% formalin at room temperature for 8 hours but not to exceed 24 hours. (For small rodent tissue, it is recommended to fix tissues for 4-8 hours.)*
2. *Rinse with running tap water for 1 hour.*
3. *Embedding and sectioning of paraffin blocks*

以上内容来自：[http://www.bdbiosciences.com/in/resources/protocols/paraffin\\_sections.jsp](http://www.bdbiosciences.com/in/resources/protocols/paraffin_sections.jsp)

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**Table 3** Paired T-test of mean cytokeratin scores of stored slides stored for 3 months with nitrogen desiccation and/or paraffin coating

Condition	Automated score	
	% Fresh	P-value
Fresh	100	
Uncoated, room-air	12	<0.0001
Uncoated, nitrogen	54	<0.0001
Paraffin coated, room-air	69	<0.0001
Paraffin coated, nitrogen	96	0.4698

**Table 4** Paired T-test of mean Ki-67 scores of stored slides stored for 3 months with nitrogen desiccation and/or paraffin coating

Condition	Manual score	
	% Fresh	P-value
Fresh	100	
Uncoated, room-air	48	<0.0001
Uncoated, nitrogen	59	<0.0001
Paraffin coated, room-air	48	<0.0001
Paraffin coated, nitrogen	72	<0.0001

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
**Proteins specifically expressed in glomerulus**

The process of urine formation begins in the glomerulus, where an ultrafiltrate of plasma is formed, and the filtered fluid enters the renal tubules. The filter consists of three layers, the fenestrated endothelium, the basement membrane, and the podocyte slit diaphragm. The analysis of the glomerulus elevated proteins is well in line with the function of the glomerulus as a filtration apparatus assembling a slit diaphragm. The list of kidney elevated proteins includes several well-known glomeruli associated genes, such as podocin (NPHS2) and nephrin (NPHS1), well established as proteins creating the filtration diaphragm making up a filter for large molecules. In addition, KIRREL1 is present in the glomerulus as described before. This latter protein is not identified as elevated in kidney, since the placenta shows higher mRNA levels than kidney for this gene. The placenta also acts as a filtration machinery for large and small molecules is therefore interesting. Another kidney elevated gene expressed in the glomerulus is FGF1.

**Proteins specifically expressed in proximal tubule**

Approximately 60% of the filtered  $\text{Na}^+$ ,  $\text{Cl}^-$ ,  $\text{K}^+$ ,  $\text{Ca}^{2+}$ , and  $\text{H}_2\text{O}$  and more than 90% of the filtered  $\text{HCO}_3^-$  are absorbed along the proximal tubule. This is also the segment that normally reabsorbs virtually all the filtered glucose and amino acids. An additional function is the secretion of numerous organic anions and cations. Most of the proteins elevated in the kidney are localized to the proximal tubule, which is in line with the function of proximal

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1. DiVito, K.A., et al., *Long-term preservation of antigenicity on tissue microarrays*. Lab Invest, 2004. **84**(8): p. 1071-8.
2. Economou, M., et al., *Proper paraffin slide storage is crucial for translational research projects involving immunohistochemistry stains*. Clin Transl Med, 2014. **3**(1): p. 4.
3. Gelb, A.B., V.A. Freeman, and S.H. Astrow, *Evaluation of methods for preserving PTEN antigenicity in stored paraffin sections*. Appl Immunohistochem Mol Morphol, 2011. **19**(6): p. 569-73.
4. <https://tissue-microarray-slides.jimdo.com/2017/10/16/storage-options-for-ffpe-tissue-sections/>.
5. Karlsson, C. and M.G. Karlsson, *Effects of long-term storage on the detection of proteins, DNA, and mRNA in tissue microarray slides*. J Histochem Cytochem, 2011. **59**(12): p. 1113-21.
6. Ijsselsteijn, M.E., et al., *A 40-Marker Panel for High Dimensional Characterization of Cancer Immune Microenvironments by Imaging Mass Cytometry*. Front Immunol, 2019. **10**: p. 2534.