




Living Image® Software 4.4 for IVIS® Spectrum

May 2013



PerkinElmer Health Sciences

940 Winter Street, Waltham, Massachusetts 02451 USA

Telephone: 800.762.4000 (US) or +1.203.925.4602

Fax: +1.203.944.4904

For sales: CustomerCareUS@perkinelmer.com

For Technical Support: tech.support@caliperLS.com

877.522.497.3302 (US) or +1.508.435.9500

www.PerkinElmerls.com

Trademarks

Discovery in the Living Organism, IVIS and Living Image are either registered trademarks or trademarks of PerkinElmer Health Sciences, Inc. and/or its parent, affiliates, and/or subsidiary companies (collectively "PerkinElmer"). The names of companies and products mentioned herein may be the trademarks of their respective owners. Apple, Macintosh and QuickTime are registered trademarks of Apple Computer, Inc. Microsoft, PowerPoint and Windows are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries. Adobe and Illustrator are either registered trademarks or trademarks of Adobe Systems Incorporated in the United States and/or other countries. Reproduction of this publication or parts thereof in any form is expressly prohibited without the express written permission of PerkinElmer. Any errors or omission which may have occurred in this publication despite the utmost care taken in its production will be corrected as soon as possible, but not necessarily immediately upon detection. PerkinElmer provides this publication "As Is" without warranty of any kind, either express or implied, including but not limited to the implied warranties of merchantability or fitness for a particular purpose. Some states or jurisdictions do not allow disclaimer of express or implied warranties in certain transactions; therefore, this statement may not apply to you. PerkinElmer reserves the right to revise this publication and to make changes from time to time in the content hereof without obligation of PerkinElmer to notify any person of such revision or changes

Copyright

© 2013 PerkinElmer Health Sciences, Inc. and its parent, affiliated, and subsidiary companies. All rights reserved, including but not limited to those rights to reproduce this publication or parts thereof.

Contents

Chapter 1	Welcome	1
1.1	About This Manual	1
1.2	What's New in Living Image 4.4 Software	2
1.3	Living Image Help	3
1.4	Contact Information	4
Chapter 2	Getting Started	5
2.1	Starting the Living Image Software	5
2.2	Initializing the System and Checking Temperature	7
	Initializing the IVIS Spectrum	8
	CCD Temperature	8
	Stage Temperature	9
2.3	Overview of Image Acquisition	9
	Auto Exposure Feature	10
	Imaging Modes on IVIS Spectrum	10
2.4	Overview of Living Image Tools and Functions	12
2.5	Managing User Accounts	19
	Adding Users	19
	Changing or Adding Passwords	20
	Deleting Users	20
	Locking User Accounts	21
2.6	Tracking System and User Activity	22
	Activity Window	22
Chapter 3	Image Acquisition	23
3.1	Luminescent Imaging	23
	Quick Guide	23
	Acquire a Luminescent Image	24
3.2	Cherenkov Imaging	28
3.3	Fluorescent Imaging With Epi-Illumination	29
	Quick Guide	29
	Acquire a Fluorescent Image With Epi-Illumination	31
3.4	Fluorescent Imaging With Transillumination	35
3.5	Acquire a Sequence Using the Imaging Wizard	39
	Starting the Imaging Wizard	40
	Acquire the Sequence	41
3.6	Acquire Multiple Sequences in Batch Mode	45
3.7	Manually Set Up a Sequence	47
	Editing Image Parameters	49
	Inserting Images in a Sequence	49
	Removing Images From a Sequence	50
3.8	Manually Save Image Data	51

3.9	Exporting Images	51
Chapter 4	Working With Optical Image Data	52
4.1	Loading Optical Image Data	52
	Preview and Load Data Using the Living Image Browser	52
	Load Data From the Menu Bar or Toolbar	56
	About the Image Window and Tool Palette	57
	Organizing Images	60
4.2	Adjusting Image Appearance	60
	Smoothing and Binning	62
	Zooming or Panning	63
4.3	Correcting Optical Image Data	64
4.4	Viewing Intensity Data	65
	X,Y Coordinates and Intensity Data	65
	Line Profile	66
4.5	Measuring Distance	68
	Distance Measurement Tool	68
	Image Crop Box	68
4.6	Managing Image Information	69
	Viewing Image Information	69
	Editing the Image Label	70
	Adding Comments	72
	Tagging an Image	73
4.7	Creating a Transillumination Overview	73
4.8	Overlaying Multiple Images	74
4.9	Rendering Intensity Data in Color	77
4.10	Exporting or Printing Images	79
4.11	Managing Image Sequences	81
	Editing a Sequence	81
	Creating a Sequence From Individual Images	82
Chapter 5	ROI Tools for Optical Data.	83
5.1	About ROIs	83
5.2	Quick Guide: Draw Measurement ROIs on an Optical Image or Sequence	86
5.3	ROI Tools for Optical Images	87
5.4	Measurement ROIs	89
	Automatically Draw ROIs	90
	Manually Draw ROIs	92
	Draw ROIs Using the Free Draw Method	93
5.5	Mirror ROIs	94
5.6	Measuring Background-Corrected Signal	97
5.7	Managing ROI Properties	100
	Viewing ROI Properties	100
	Moving an ROI	103
	Editing ROI Dimensions	104

	Editing the ROI Line	106
	Move or Edit the ROI Label	107
	Save, Load, or Delete ROIs	107
5.8	Managing the ROI Measurements Table	109
	Viewing the ROI Measurements Table	109
	Configuring the ROI Measurements Table	111
	Copying or Exporting the ROI Measurements Table	113
Chapter 6	3D ROI Tools for Volumetric Data	114
6.1	About 3D ROIs	114
6.2	Drawing a 3D ROI	115
	ROI Properties	119
6.3	Managing the 3D ROI Measurements Table	120
	Configuring the 3D ROI Measurements Table	120
	Copying or Exporting the ROI Measurements Table	122
Chapter 7	Image Math	123
7.1	About Image Math	123
7.2	Creating a New Image Using Image Math	123
7.3	Subtracting Tissue Autofluorescence	125
Chapter 8	Spectral Unmixing	128
8.1	About Spectral Unmixing	128
8.2	Acquire an Image Sequence for Spectral Unmixing	128
	Bioluminescence Imaging	129
	Fluorescence Imaging	131
	Cherenkov Imaging	134
8.3	Spectral Unmixing Methods	137
	Guided Method	137
	Library Method	140
	Automatic Method	142
	Manual Method	145
8.4	Correcting Spectra	149
8.5	Spectral Unmixing Results	150
	Spectra Plot	151
	Composite Image	153
	Analyzing Images	154
	Managing Spectral Unmixing Results	155
Chapter 9	DyCE Imaging and Analysis	156
9.1	About DyCE (Dynamic Contrast Enhancement)	156
9.2	Acquire an Image Sequence for DyCE Analysis	157
	Bioluminescence Imaging	157
	Fluorescence Imaging	159
	Cherenkov Imaging	162
9.3	DyCE Analysis	164
	Automatic DyCE Analysis	164

	Manual DyCE Analysis	168
9.4	DyCE Results	171
	Viewing Unmixed Images	171
	Viewing the Composite Image	172
	Correcting Temporal Spectra	174
Chapter 10	Reconstructing a 3D Surface	176
10.1	Generating a Surface	177
	Changing the View Perspective	180
10.2	Managing Surfaces	181
	Export or Import a Surface	182
Chapter 11	3D Reconstruction of Sources	183
11.1	Overview of Reconstructing Sources	183
	General Considerations	185
	Reconstruction Inputs	185
11.2	Acquire a Sequence	186
	Bioluminescence Imaging	186
	Fluorescence Imaging	188
	Manual Sequence Setup	191
11.3	Steps to Reconstruct Luminescent Sources Using DLIT	192
11.4	Steps to Reconstruct Fluorescent Sources Using FLIT	197
11.5	Including or Excluding Data for 3D Reconstruction	200
11.6	3D Reconstruction Results	203
	DLIT or FLIT Results	203
	Managing 3D Reconstruction Results	204
11.7	Checking the Reconstruction Quality	205
	Viewing Photon Density or NTF Efficiency Maps	206
11.8	Measuring Sources	207
	Determining the Source Center of Mass	207
	Measuring Source Depth	209
	Viewing Location Coordinates	209
	Displaying Slices Through a Reconstruction	210
11.9	Viewing Luminescent and Fluorescent Sources in One Surface	211
11.10	Comparing Reconstruction Results	211
	Viewing Results in the Longitudinal Study Window	211
	Measuring Intensity	214
	Viewing Plots	215
11.11	Exporting a 3D Scene as DICOM	216
	Viewing DICOM Data	217
11.12	3D Optical Surface Tools	218
11.13	3D Optical Source Tools	220
11.14	3D Optical Registration Tools	222
11.15	3D Animation	229
	Viewing a Preset Animation	231

Creating a Custom Animation	231
Managing Animation Setups	233
11.16 DLIT/FLIT Troubleshooting	234
Chapter 12 Quantification Database	235
12.1 Preparing and Imaging the Samples	235
12.2 Creating a Quantification Database	236
12.3 Managing Quantification Results	240
Exporting Quantification Results	241
Chapter 13 3D Multi-Modality Tools	242
13.1 About the 3D Multi-Modality Tools	242
3D Multi-Modality Tool Requirements	242
13.2 Classifying 3D Volumetric Data	243
Specifying a Color-Opacity Map	244
13.3 Volume Display Options	247
Adjusting Image Quality	247
Adjusting Volume Opacity	248
Maximum Intensity Projection	249
Gradient Illumination	250
Modifying Volume Resolution	251
13.4 Smoothing a Volume	252
13.5 Viewing and Rendering Slices	253
Viewing Slices	253
Rendering Slices	255
13.6 Volume Information and Results	256
Managing Results	256
13.7 Registering Optical and Volumetric Data	257
Loading Data for Registration	259
Registering Multi-Modal Data	261
13.8 Volume Data Viewer	265
13.9 Viewing RAW Volumetric Data	266
Appendix A IVIS Acquisition Control Panel	269
A.1 Control Panel	269
A.2 Manually Setting the Focus	272
Appendix B Preferences	274
B.1 General Preferences	274
B.2 Options	276
B.3 Acquisition	277
B.4 Theme	278
B.5 Optical Properties	281
Appendix C Menu Commands, Toolbars, and Shortcuts	282
Index	286

1 Welcome

[About This Manual](#)

[What's New in Living Image 4.4 Software](#) on page 2

[Living Image Help](#) on page 3

[Contact Information](#) on page 4

1.1 About This Manual



NOTE: This *Living Image Software 4.4 Software Manual* (PN CLS137703 rev A) is only for use with the IVIS Spectrum instrument.

This manual explains how to acquire optical and volumetric image data on the IVIS Spectrum and analyze the data using the Living Image® software. The manual provides detailed instructions and screenshots. Sometimes the screenshots in the manual may not exactly match those displayed on your screen.

When analyzing data acquired on a different type of IVIS instrument, say for example the IVIS Spectrum CT, please see the Living Image Software Manual specific for the IVIS Spectrum CT.

Table 1.1 Living Image 4.4 Software Manuals


Living Image Software Manual for:	Part Number
IVIS Lumina Series III Note: This manual can also be used with the IVIS Lumina II.	CLS137701 rev A
IVIS Lumina XR Series III Note: This manual can also be used with the IVIS Lumina XR.	CLS137702 rev A
IVIS Lumina K Series III Note: This manual can also be used with the IVIS Lumina Kinetic.	CLS137700 rev A
IVIS 200	CLS137699 rev A
IVIS Spectrum	CLS137703 rev A
IVIS Spectrum CT	CLS137704 rev A

Please see the *IVIS Spectrum Hardware Manual* (PN 121450_Rev00) for information on the IVIS Spectrum instrument.

1.2 What's New in Living Image 4.4 Software


Living Image 4.4 software controls optical image acquisition on the IVIS® Spectrum instrument, and provides tools for optimizing image display and analyzing images. [Table 1.1](#) shows the major new or improved software features (specific to acquisition and analysis of IVIS Spectrum image data).

Table 1.1 Living Image 4.4 – New or Improved Features

New or Improved Features	See Page
Imaging Wizard includes Cherenkov imaging mode.	28
Fluorescent probes <ul style="list-style-type: none"> Expanded probe list includes the entire PerkinElmer catalog of fluorescent probes. A list of recently used fluorescent probes (up to 10 probes) is available in the Probes drop-down menu in the Imaging Wizard and control panel. 	
Image Sequence <ul style="list-style-type: none"> Ability to close and save a sequence data set while another sequence acquisition is in progress. Adjust binning or smoothing of an image sequence. 	
A warning alerts you if DLIT or FLIT sequence acquisition is set up using FOV D, which is not recommended.	
Transillumination Setup Window – Selected transillumination points will include only those points that are located under the subject and at least 5 mm from the edge (to help prevent saturated pixels).	
The  icon in the Living Image Browser indicates images with saturated pixels.	53
A warning appears if loading transillumination image data containing saturated pixels.	
Image information is available in the 3D View window.	
“Jet” color table available for CT volumes.	
Tool Palette <ul style="list-style-type: none"> Reorganized for improved ease of use. Image Information and Filtering tools incorporated into the Image Adjust section. Tool Palette, when docked in the main application window, automatically adjusts to display selected tools (expanded palette section) in full view. 	13
“Show” option in the Surface Topography tools keeps the pink data mask display on, making it easier to see the data selected while adjusting the Threshold slider position.	177
Spectral Unmixing <ul style="list-style-type: none"> Image Cube brightness control. Adjust binning or smoothing for spectral unmixing results. Display a Color Scale for unmixed images. Adjust the Color Scale min and max. 	147 151
DyCE Imaging <ul style="list-style-type: none"> Maximum number of images that can be acquired for DyCE analysis increased to 200 images. Measurement of a DyCE time interval begins at the start of the first luminescent or fluorescent exposure, after the photograph is acquired. Image Cube brightness control. 	
Longitudinal Study – Ability to set minimum and maximum color scale values for voxel display that are outside the range of the image data for improved ability to compare different data sets.	
3D ROIs <ul style="list-style-type: none"> Larger default size for 3D ROI and ROI handles for easier manipulation. Angle parameters in the ROI Properties dialog box provide another way to rotate the 3D ROI in the x, y, or z-plane. 	119

1.3 Living Image Help

There are several ways to obtain help on the software features and related information.

To view:	Do this:
A tooltip about a button function	Put the mouse cursor over the button.
A brief description about an item in the Living Image user interface	Click the  toolbar button, then click the item.
The <i>Living Image Software User's Manual</i>	Press F1 or select Help → User Guide on the menu bar and select the manual specific for your imaging system.
Living Image technical notes (see Table 1.2 on page 3)	Select Help → Tech Notes on the menu bar. Note: Please see the IVIS University download page for the most recent collection of technical notes.

[Table 1.2](#) lists the tech notes that are available under the Help menu. There are three types of tech notes:

- Tech Notes – Quick guides for tasks using the Living Image software tools.
- Biology Tech Notes – Protocols and procedures related to animal subjects.
- Concept Tech Notes – Background information on *in vivo* imaging topics.

Table 1.2 Tech Notes

Tech Notes	Title
1	Adaptive Fluorescence Background Subtraction
2	Auto-Exposure
3	Determine Saturation
4	Bioluminescence Tomography (DLIT) <ul style="list-style-type: none"> ■ 4a – Setup and Sequence Acquisition ■ 4b – Topography ■ 4c – Source Reconstruction and Analysis
5	ROIs (optical data) <ul style="list-style-type: none"> ■ 5a – Drawing ROIs ■ 5b – Subtracting Background ROI from Sequence ■ 5c – Subject ROIs
6	Fluorescence Tomography (FLIT) <ul style="list-style-type: none"> ■ 6a – Setup and Sequence Acquisition ■ 6b – Topography ■ 6c – Source Reconstruction and Analysis
7	High Resolution Images
8	Image Math
9	Image Overlay – 2D
10	Image Overlay – 3D
11	Imaging Wizard
12	Load Groups of Images
13	Spectral Unmixing

Table 1.2 Tech Notes (continued)

14	Transillumination <ul style="list-style-type: none"> ■ 14a – Transillumination Fluorescence ■ 14b – Raster Scan ■ 14c – Normalized Transmission Fluorescence ■ 14d – Well Plates
15	Well Plate Quantification
Biology Tech Notes	Title
1	d-luciferin Prep Sheet
2	Kinetic Analysis of Bioluminescent Sources
3	Imaging Protocol Guide
4	Imaging Procedure
5	Intraperitoneal Injections
Concept Tech Notes	Title
1	Luminescent Background Sources and Corrections
2	Image Display and Measurement
3	Detection Sensitivity
4	Fluorescent Imaging
5	DLIT and FLIT Reconstruction of Sources

1.4 Contact Information

If you have questions regarding this manual or Living Image software, please contact PerkinElmer Technical Support.

PerkinElmer Health Sciences
940 Winter Street, Waltham, Massachusetts 02451 USA
Telephone: 800.762.4000 (US) or +1.203.925.4602
Fax: +1.203.944.4904

For sales: CustomerCareUS@perkinelmer.com

For Technical Support: tech.support@caliperLS.com
877.522.497.3302 (US) or +1.508.435.9500

www.PerkinElmer.com

2 Getting Started

Starting the Living Image Software

Initializing the System and Checking Temperature on page 7

Overview of Image Acquisition on page 9

Overview of Living Image Tools and Functions on page 12

Managing User Accounts on page 19

Tracking System and User Activity on page 22

2.1 Starting the Living Image Software

Living Image software on the PC workstation that controls the IVIS Spectrum includes both the acquisition and analysis features. Living Image software on other workstations includes only the analysis features.

See the Installation Guide on the Living Image CD ROM for software installation instructions. [Table 2.1](#) shows the default software installation locations.

Table 2.1 Living Image Software Installation Locations


Living Image Software	Operating System	Installation Location
32-bit version	32-bit Windows	C:\Program Files\Caliper Life Sciences\Living Image
	64-bit Windows	C:\Program Files(x86)\Caliper Life Sciences\Living Image
64-bit version	64-bit Windows	C:\Program Files\Caliper Life Sciences\Living Image




NOTE: All components of the IVIS Spectrum should be left on at all times due to the long cooling time required to reach operating (demand) temperature. It is also important to leave the system on to enable automatic overnight electronic background measurements. Periodically rebooting the computer is permissible and does not affect the camera operation.

To start the software:

1. **PC Users:** Double-click the Living Image software icon  on the desktop. Alternatively, click the Windows Start button  and select **All Programs → Caliper Life Sciences → Living Image**.


Macintosh Users: Double-click the Living Image icon  on the desktop or run the software from the application folder.

The main window appears ([Figure 2.1](#)).



2. In the dialog box that appears, select a user ID from the drop-down list. If the user ID is password protected, enter the password and click **OK**.


Alternatively, create a new user ID:

- a. In the Select/Add User ID box, click the  button.
- b. Enter a user ID.
- c. Enter and confirm a password. This is optional.
- d. Click **Add** and **OK**.

The control panel appears if the workstation controls the IVIS Spectrum ([Figure 2.2](#)). For more details on the control panel, see [Appendix A on page 269](#).



NOTE: Living Image software has optional password protection for user accounts. See [page 20](#) for more details.

Figure 2.2 Living Image Main Window and IVIS Acquisition Control Panel

NOTE: Living Image software on the PC workstation that controls the IVIS® Imaging System includes both the acquisition and analysis features. Living Image software on other workstations includes only the analysis features. Macintosh users have access to the analysis features only.

2.2 Initializing the System and Checking Temperature

The IVIS Spectrum must be initialized each time Living Image software is started, or if the power has been cycled to the imaging chamber.

The initialization procedure is started from the control panel ([Figure 2.3](#)).




NOTE: The control panel is only available on the PC workstation that controls the IVIS Imaging System.


Initialization moves every motor-driven component in the system (for example, imaging stage and lens) to a home position, resets all electronics and controllers, and restores all software variables to the default settings. Initialization may be useful in error situations.

See the *IVIS® Spectrum Hardware Manual* (PN 133577_Rev A) for further details on instrument operation.

Initializing the IVIS Spectrum

1. Start the Living Image software (double-click the  icon on the desktop).
2. In the control panel that appears, click **Initialize** (Figure 2.3).

After several seconds you will hear the instrument motors move.




CCD Temperature

The IVIS Acquisition Control Panel indicates the temperature status of the charge coupled device (CCD) camera (see Figure 2.4 for a description of the temperature status colors). Immediately after initialization is completed, the temperature box is usually red and will turn green after several minutes. If this is not the case, contact PerkinElmer Technical Support (see page 4).

The demand temperature for the CCD camera is preset and generally should not be changed. Electronic feedback control maintains the CCD camera temperature to within a few degrees of the demand temperature.

The instrument is ready for imaging after the system is initialized and the operating (demand) temperature of the CCD camera is reached (locked).


Figure 2.4 Control Panel – Instrument Temperature Status



Stage Temperature

The stage is temperature-controlled to keep subjects warm during imaging. The temperature control is enabled after the instrument is powered on and initialized from the Living Image software. The default temperature is 37 °C and is self-monitoring after the system is initialized. The imaging stage may be set to a temperature from 20 - 40 °C.

Figure 2.5 Set Imaging Stage Temperature in the Control Panel




2.3 Overview of Image Acquisition

The control panel provides the image acquisition functions (Figure 2.6). See [Appendix A on page 269](#) for details on the imaging parameters in the control panel.



NOTE: The control panel is only available on the PC workstation that controls the IVIS Imaging System. The items available in the control panel depend on the selected imaging mode (luminescent, fluorescent) and acquisition mode (Image Setup or Sequence Setup).

Figure 2.6 IVIS Acquisition Control Panel



Auto Exposure Feature

The Auto exposure setting is useful in situations where the signal strength is unknown or varies widely, for example during a time course study. If Auto exposure is chosen (Figure 2.6), the system acquires an image at maximum sensitivity, then calculates the required settings to achieve, as closely as possible, an image with a user-specified target max count. If the resulting image has too little signal or saturated pixels, the software adjusts the parameters and takes another image.

In most cases, the default auto exposure settings provide a good luminescent or fluorescent image. However, you can modify the auto exposure preferences to meet your needs. See page 277 for more details.

Imaging Modes on IVIS Spectrum

Table 2.1 briefly explains the types of images that can be acquired on the IVIS Spectrum.

Table 2.1 Imaging Modes – IVIS® Spectrum









Imaging Mode	Description	Example
<div>Luminescent optical imaging</div> <div></div>	<div>A longer exposure of the subject taken in darkness to capture low level luminescence emission from the surface of the subject.</div> <div>The optical luminescent image data is displayed in pseudocolor that represents intensity.</div>	<div></div> <div>Luminescent image</div> <div></div> <div>Overlay: Luminescent image on photograph</div>

Table 2.1 Imaging Modes – IVIS® Spectrum (continued)

Imaging Mode	Description	Example
<div>Fluorescent optical imaging</div> <div></div>	<p>An exposure of the subject illuminated by filtered light. The light source is located above the imaging stage (epi-illumination). The target fluorophore emission is captured and focused on the CCD camera.</p> <p>The optical fluorescent image data can be displayed in units of counts or photons (absolute, calibrated), or in terms of efficiency (calibrated, normalized).</p> <p>Note: See the concept tech note Image Display and Measurement for more on quantifying image data (select Help → Tech Notes on the menu bar).</p>	<div></div> <div>Fluorescent image</div> <div></div> <div>Overlay: Fluorescent image on photograph</div>
<div>Photograph</div> <div></div>	<p>A short exposure of the subject illuminated by the lights located in the ceiling of the imaging chamber. The photographic image is displayed as a grayscale image.</p>	<div></div>

2.4 Overview of Living Image Tools and Functions

Living Image tools are organized in the Tool Palette or under "Tools" in the menu bar (Figure 2.7). Some tools are for use with a single image, others require an image sequence.

Table 2.2 provides an overview of the tools available for data acquired on the IVIS® Spectrum. If analyzing data acquired on a different type of IVIS instrument, say for example the IVIS Spectrum CT, please see the Living Image Software Manual specific for that instrument.



NOTE: The tools available in the Tool Palette or menu bar depend on the active image data.




Table 2.2 Living Image Tools Available for Data Acquired on the IVIS Spectrum




Living Image Tools and Functions	See Page
	<p>Image Adjust 65</p> <ul style="list-style-type: none"> Tools for image display (zoom, crop, or pan). Adjust photo display – Tune the photograph brightness, contrast, or opacity. Manage the color table for image display. Apply smoothing or binning to an image. View optical image data (counts or radiance) at an X,Y location. Measure distance in an image. View a line plot of pixel intensities.
	<p>Corrections (Optical Data) 65</p> <p>Apply corrections to the raw data.</p>
	<p>ROI Tools for Optical Data 83</p> <p>Specify a region of interest (ROI) in an optical image and measure the signal intensity within the ROI.</p> <p>3D ROI Tools for Volumetric Data 114</p> <p>Specify a 3-dimensional region of interest (ROI) in a CT volume and measure the signal intensity within the 3D ROI.</p>

Table 2.2 Living Image Tools Available for Data Acquired on the IVIS Spectrum (continued) (continued)

Living Image Tools and Functions

See Page

Spectral Unmixing

128

Use spectral unmixing to:

- Extract the signal of one or more fluorophores from the tissue autofluorescence
- Analyze luminescent or fluorescent images when more than one reporter is used in the same animal model

DyCE (Dynamic Contrast Enhancement)

156

Use DyCE to:

- Determine real-time pharmacokinetic (spatio-temporal biodistribution) of a probe or dye signal
- Extract "temporal spectra" (signal intensity as a function of time) from particular anatomical regions.

Note: DyCE acquisition and analysis tools require a separate license.

Surface Topography

176

Generate 3D reconstruction of the animal surface (*topography*) derived from the CT image.

A surface is a required input for:

- DLIT (*diffuse light tomography*) analysis which generates a 3D reconstruction of luminescent sources
- FLIT (*fluorescence imaging tomography*) analysis which generate a 3D reconstruction of fluorescent sources

Table 2.2 Living Image Tools Available for Data Acquired on the IVIS Spectrum (continued) (continued)




Living Image Tools and Functions		See Page
	<p>3D Multi-Modality Tools</p> <p>Set color and opacity values for different intensity ranges of a CT volume so that the color-opacity map shows the volume regions you are interested in (opaque in the map) and hides unimportant regions.</p> <p>Co-register 3D reconstructions of luminescent or fluorescent sources (biological information) with a CT volume to provide anatomical context for interpreting biological (functional) information.</p> <p>Note: The 3D Multi-Modality tools require a separate license.</p>	242
	<p>3D Optical Tools</p> <p>Surface tools – Adjust the appearance of the reconstructed animal surface and the photon density maps.</p> <p>Source tools – Adjust the appearance of reconstructed sources, make source measurements, export voxel measurements.</p> <p>Registration tools – Display organs on the reconstructed surface, adjust the location or scale of organs on the surface, import an organ atlas.</p>	218
	<p>DLIT 3D Reconstruction</p> <p><i>Diffuse light tomography</i> (DLIT) analysis provides a complete 3D reconstruction of the luminescent source distribution within the subject. The 3D reconstruction is presented as volume elements called <i>voxels</i>.</p> <p>If a luminescent calibration database is available, the number of cells per source can be determined in addition to source intensity (photons/ sec).</p>	186

Table 2.2 Living Image Tools Available for Data Acquired on the IVIS Spectrum (continued) (continued)

Living Image Tools and Functions

See Page

#	ExWL	EmWL	Threshold %
01	745	800	5.1
02	745	800	5.1
03	745	800	5.1
04	745	800	5.0
05	745	800	5.1
06	745	800	5.1
07	745	800	5.1
08	745	800	5.1
09	745	800	5.2
10	745	800	5.1

FLIT 3D Reconstruction

197

Fluorescent imaging tomography (FLIT) analysis provides a complete 3D reconstruction of the fluorescent source distribution within the subject. The 3D reconstruction is presented as volume elements called *voxels*.

If a fluorescent calibration database is available, the number of fluorophore molecules or cells per source can be determined in addition to the total fluorescence yield.

Preset Animations: Spin CW on X-Axis

Frame Factor: 1

Animation Setup

Time Scale %: 0

Key Frame 1
Key Frame 2
Key Frame 3
Key Frame 4
Key Frame 5

Play Record

Frames Per Second: 10

Total Duration (secs): 5

Load Save

3D Animation Tools

229

Select **Tools** → **3D Animation** on the menu bar. Creates an animation from a sequence of 3D views (keyframes). For example, an animation can depict a rotating 3D scene. The animation (series of key frames) can be recorded to a movie file.

Longitudinal Study

211

Select **Tools** → **Longitudinal Study** on the menu bar. Multiple DLIT and/or FLIT reconstruction results can be viewed side-by-side in the Longitudinal Study window. The Longitudinal Study window provides a convenient way to compare different results, for example, results obtained at different time points or results from different types of reporters. Voxel intensity within the entire surface or a user-selected area can be measured in all results in the Longitudinal Study window.

FLIT 3D Reconstruction

Fluorescent imaging tomography (FLIT) analysis provides a complete 3D reconstruction of the fluorescent source distribution within the subject. The 3D reconstruction is presented as volume elements called *voxels*.

If a fluorescent calibration database is available, the number of fluorophore molecules or cells per source can be determined in addition to the total fluorescence yield.

197

| 3DAnimation Preset Animations Presets: Spin CW on X-Axis Frame Factor: 1 Animation Setup Time Scale %: 0 Key Frame 1 Key Frame 2 Key Frame 3 Key Frame 4 Key Frame 5 Play Record Frames Per Second: 10 Total Duration (secs): 5 Load Save |

3D Animation Tools

Select **Tools** → **3D Animation** on the menu bar.

Creates an animation from a sequence of 3D views (keyframes). For example, an animation can depict a rotating 3D scene. The animation (series of key frames) can be recorded to a movie file.

229

| Longitudinal Study Window Select a Sequence Output Files Select Analysis Results Display Results Display Surface Display Fractional Density Map 3D View Angle Of View Intensity 0.0 0.5 1.0 1.5 2.0 |

Longitudinal Study

Select **Tools** → **Longitudinal Study** on the menu bar.

Multiple DLIT and/or FLIT reconstruction results can be viewed side-by-side in the Longitudinal Study window.

The Longitudinal Study window provides a convenient way to compare different results, for example, results obtained at different time points or results from different types of reporters.

Voxel intensity within the entire surface or a user-selected area can be measured in all results in the Longitudinal Study window.

211

Table 2.2 Living Image Tools Available for Data Acquired on the IVIS Spectrum (continued) (continued)


Living Image Tools and Functions	See Page
<div></div>	<div>Well Plate Quantification 235</div> <p>Select Tools → Well Plate Quantification for <sequence name> on the menu bar.</p> <p>Generate a database of luminescence or fluorescence signal intensities by analyzing images of known serial dilutions of luminescent or fluorescent cells or dye molecules.</p> <p>Use the quantification database to extrapolate the number of cells in a DLIT source or the number of dye molecules or cells in a FLIT source.</p>

Table 2.2 Living Image Tools Available for Data Acquired on the IVIS Spectrum (continued) (continued)





Living Image Tools and Functions		See Page
<div></div> <p>Images of Quantum dot nanocrystals (700 or 800 nm) were acquired using different combinations of excitation and emission filters.</p> <div></div> <p>Colorize view of the combined images</p>	<p>Colorize View</p> <p>Select Tools → Colorize for <sequence name> on the menu bar.</p> <p>The colorize tool renders each luminescence or fluorescence image of a sequence in color, and combines them into a single image. This enables you to see both intensity and spectral information in a single view.</p> <p>The tool provides a useful way to visualize multiple probes or scale probe signals that are not in the visible range.</p>	77
<div></div>	<p>Transillumination Overview</p> <p>Select Tools → Transillumination Overview for <sequence name> on the menu bar.</p> <p>The transillumination overview tool combines the images of a FLIT sequence (a fluorescence sequence acquired in transillumination mode) into a single image. All of the individual fluorescent signals are stacked over one photograph and the intensity is summed. One overview is created per filter pair. If two filter pairs were used during acquisition, then two overview images will be created</p>	73

Table 2.2 Living Image Tools Available for Data Acquired on the IVIS Spectrum (continued) (continued)


Living Image Tools and Functions	See Page
	<p>Image Math Window 123</p> <p>Select Tools → Image Math for <sequence name> on the menu bar.</p> <p>Mathematically combine (add, multiply, subtract, or divide) two user-selected images.</p> <p>For example, subtract a blue-shifted background filter image from the primary excitation filter image to remove tissue autofluorescence signal.</p>

2.5 Managing User Accounts

Adding Users

New users can be created in the:


- Main window at startup (see [page 6](#)).
 - User Settings dialog box ([Figure 2.8](#)).
1. Select **Edit → User settings** on the menu bar.
 2. Click the Add User tab in the dialog box that appears.



3. Enter a user ID.
4. Optional: enter and confirm a password.
5. Click **Add**.

Changing or Adding Passwords

1. Select **Edit → User settings** on the menu bar.
2. Click the Change Password tab in the dialog box that appears.




3. Select a User ID, enter and confirm a new password, and click **Submit**.

Deleting Users



NOTE: User accounts can be locked. If this security is applied, a master password is required to delete users from the system. See [page 21](#) for more details on locking user accounts.

1. Select **Edit → User settings** on the menu bar.
2. Click the Delete User tab in the dialog box that appears.



3. Select a User ID.
4. If the accounts are locked, enter the master password.
5. Click **Delete** and **Close**.


Locking User Accounts

If user accounts are locked, a master password is required to change user passwords, delete users, or unlock user accounts.

To lock user accounts:

1. Select **Edit → User settings** on the menu bar.
2. Click the Security tab in the dialog box that appears.
3. Click **Lock User Accounts**.

Figure 2.11 User Settings – Security




4. Enter and confirm a master password. Click **Close**.
The master password will be required to delete users.

To unlock user accounts:

1. In the Security tab, click **Unlock User Accounts**.
2. Enter the master password and click **Unlock**. Click **Close**.

Figure 2.12 User Settings – Security




2.6 Tracking System and User Activity

Activity Window

The Activity window shows the imaging system activities (Figure 2.13). This information may be useful for PerkinElmer field service engineers to understand the imaging system behavior over time or for troubleshooting. The activity log is located at C:\Program Files\Caliper Life Sciences\Living Image.

The software tracks user time on the system (hr/min/sec per user ID) from logon until switching users or system shut down. The software creates a separate record for each month (for example, LI_USAGE_<MONTH>_2011.csv) located at C:\Program Files\Caliper Life Sciences\Living Image\Usage).



3 Image Acquisition

Luminescent Imaging

Cherenkov Imaging on page 28

Fluorescent Imaging With Epi-Illumination on page 29

Fluorescent Imaging With Transillumination on page 35

Acquire a Sequence Using the Imaging Wizard on page 39

Acquire Multiple Sequences in Batch Mode on page 45

Manually Set Up a Sequence on page 47

Manually Save Image Data on page 51

Exporting Images on page 51

3.1 Luminescent Imaging

Luminescent imaging captures signals from luminescent molecular reporters. This section explains how to acquire a single luminescent optical image:

- Quick guide – See below.
- Detailed instructions – See [page 24](#).

See [page 39](#) for information on acquiring a luminescent sequence using the Imaging Wizard.


Quick Guide

1. Start the Living Image® software  and initialize the IVIS Spectrum ([page 7](#)).




NOTE: See the *IVIS Spectrum Hardware Manual* (PN PN121450_Rev00) for more information on the instrument

2. Place the anesthetized subjects in the imaging chamber and close the door.
3. In the control panel ([Figure 3.1](#)):
 - a. Put a check mark next to “Luminescent” and select “Auto” exposure in the control panel.
 - b. Choose “Photograph” (optional). Selecting “Photograph” automatically selects “Overlay”.
 - c. Select “Use subject height” and enter the height in centimeters.
 - d. Click **Acquire**.



4. Select a location for the image data when prompted (optional).
Image data acquired during the session will be automatically saved to this location.
5. Enter experiment and subject information in the dialog box that appears (optional).
The image window and tool palette appear when acquisition is finished ([Figure 3.2](#)).




Acquire a Luminescent Image

This section provides detailed instructions for image acquisition.



NOTE: The IVIS® Spectrum should be initialized and the temperature locked before setting the imaging parameters in the control panel. See [page 7](#) for more details.


1. Put a check mark next to **Luminescent** and select **Auto** exposure (click the  arrows) in the control panel.

The software automatically determines the binning and F/Stop settings.



TIP: See the tech note *Auto-Exposure* for helpful information (select **Help** → **Tech Notes** on the menu bar).

Alternatively, manually set the exposure, binning, and F/Stop. See [Appendix A on page 269](#) for details on these parameters.



- 2. Put a check mark next to Photograph.
- 3. Select a Field of View (size of the stage area to be imaged).



 **TIP:** See the technical note *Detection Sensitivity* for more information about the Field of View (select **Help** → **Tech Notes** on the menu bar).

Table 3.1 Field of View (FOV) Settings – IVIS® Spectrum

FOV Setting	FOV (cm)
A	4
B	6.5
C	13
D	22.5 (19.5)*
E	22.5 (26)*

*Some IVIS Spectrum instruments may have the FOV in parentheses. FOV 19.5 and 26 were replaced by FOV 22.5.

- 4. Select a focus option in the control panel (Figure 3.4).
The focal distance to the camera is set at stage $z = 0$ for each field of view. To focus at the top of the animal, the stage moves down so that the top of the animal is at $z = 0$. You can enter the height of the animal and select the “use subject height” option or use the manual focus option to determine the proper subject height for the area to be imaged. See Appendix A on page 272 for manual focus instructions.

Figure 3.4 Choose a Focus Option in the Control Panel

5. Select **Overlay** to view an overlay image (registered photograph and luminescent image) after acquisition.



NOTE: If you want to check the subject inside the chamber before acquisition, take a photograph—uncheck the Luminescent option, choose the Photograph option, and click **Acquire**. Be sure to select the Luminescent option after taking the photograph.

6. Click **Acquire** when you are ready to capture the image.




NOTE: If necessary click **Image Setup** in the control panel to operate in single image mode. In single image mode, the **Sequence Setup** button appears in the control panel. Use this button to set up sequence acquisition (see [page 39](#) for more details on sequence setup).

7. Enter information about the image in the Edit Image Labels box that appears (optional). Click **OK**.




NOTE: You can enter image label information at any time during or after acquisition. If you do not want to enter image information, click **Cancel**. See [page 70](#) for details on adding information to an image after acquisition.

Figure 3.5 Enter Information About the Image (optional)


Information entered here appears in the image label (see [Figure 3.7 on page 27](#)).

If this is the first image of the session, you are prompted to enable the autosave function (Figure 3.6). When Autosave is enabled, all images acquired during the session are automatically saved to a user-selected location. A different location can be chosen at any time (select **Acquisition** → **Auto-Save** on the menu bar).






8. Click **Yes** in the prompt to enable autosave, then choose a location in the dialog box that appears. Alternatively, click **No** in the prompt and manually save the image data. See page 51 for details. Image acquisition begins and the upper area of the control panel changes to red color. During acquisition, the **Acquire** button in the control panel becomes a **Stop** button. Click **Stop** to cancel acquisition.

The control panel returns to blue color when acquisition is finished and the image window appears (Figure 3.7).



TIP: See the tech note *Determine Saturation* for information on pixel measurements (select **Help** → **Tech Notes** on the menu bar).

Table 3.2 Image Window

Item	Description
Units	Select the measurement units for the image display from this drop-down list. The available units depend on the type of image data. See the concept tech note <i>Image Display and Measurement</i> for more details (select Help → Tech Notes on the menu bar).
Display	A list of image types available for display, for example, overlay. For more details on the different types of image displays, see Table 2.1 on page 10 .
Info	Click to display or hide the image label. The image label includes information you enter in the Edit Image Labels dialog box (Figure 3.7) and other image information automatically recorded by the software.
	Opens a dialog box that enables you to export the active view as a graphic file.
	Creates a preview picture (<i>snapshot</i>) of the image or thumbnails that the Living Image Browser displays when the data are selected in the browser. For more details on the browser, see page 52 .
	Preview picture of the data selected in the browser (blue row)
Color Scale	Provides a reference for the pixel intensities in a luminescent or fluorescent image. Pixels less than the color scale minimum do not appear in the image. Pixels greater than the color scale maximum are displayed in the maximum color.

3.2 Cherenkov Imaging

Cherenkov luminescent imaging captures optical photons produced by Cherenkov radiation from radiotracer probes. Exposure times of several minutes are required because the Cherenkov signal is extremely low.

See [page 39](#) for information on acquiring a Cherenkov sequence using the Imaging Wizard.

To acquire a Cherenkov image:


1. Put a check mark next to **Luminescent** and enter an appropriate exposure time (in minutes) in the control panel ([Figure 3.8](#)).

The software automatically determines the binning and F/Stop settings.



TIP: See the tech note *Auto-Exposure* for helpful information (select **Help → Tech Notes** on the menu bar).

Alternatively, manually set the exposure, binning, and F/Stop. See [Appendix A on page 269](#) for details on these parameters.



2. Perform [step 2 on page 25](#) to [step 8 on page 27](#).

3.3 Fluorescent Imaging With Epi-Illumination

Fluorescent imaging captures signals from fluorescent molecular reporters.

This section explains how to acquire a single fluorescent optical image with epi-illumination (excitation light source located above the stage):

- Quick guide – See [Figure 3.10 on page 30](#).
- Detailed instructions – See [page 24](#).

See [page 35](#) for information on fluorescent imaging with transillumination (excitation light source located below the stage). See [page 39](#) for information on acquiring a fluorescent sequence using the Imaging Wizard.



TIP: See the concept tech note *Fluorescent Imaging* for more about fluorescence imaging theory (select **Help** → **Tech Notes** on the menu bar).


Quick Guide

1. Start Living Image software  and initialize the IVIS Spectrum ([page 7](#)).




NOTE: See the *IVIS Spectrum Hardware Manual* (PN 121450_Rev00) for more information on the instrument

2. Place the anesthetized subjects in the imaging chamber and close the door.
3. In the control panel ([Figure 3.9](#)):
 - a. Put a check mark next to “Fluorescent” and select “Auto” exposure in the control panel.
 - b. Select an excitation and emission filter.
 - c. Choose “Photograph” and “Overlay”.
 - d. Select “Use subject height” and enter the height in centimeters.
 - e. Click **Acquire**.

Figure 3.9 Control Panel

4. When prompted, select a location for the image data (optional).
Image data acquired during the session will be automatically saved to this location.
5. Enter experiment and subject information in the "Edit Image Labels" dialog box that appears (optional).
The image window and tool palette appear when acquisition is finished ([Figure 3.10](#)).

Figure 3.10 Image Window and Tool Palette

Acquire a Fluorescent Image With Epi-Illumination

This section provides detailed instructions for image acquisition.



NOTE: The IVIS Spectrum should be initialized and the temperature locked before setting the imaging parameters in the control panel. See [page 7](#) for more details.


1. Put a check mark next to **Fluorescent** and select **Auto** exposure (click the  arrows) in the control panel.

The software automatically determines the binning and F/Stop settings.




TIP: See the tech note *Auto-Exposure* for helpful information (select **Help** → **Tech Notes** on the menu bar).

Alternatively, manually set the exposure, binning, and F/Stop. See [Appendix A on page 269](#) for details on these parameters.



2. Select an excitation and emission filter from the drop-down lists.

The instrument has 18 narrow band excitation filters that span 490-850nm with a 20nm bandwidth, enabling spectral scanning over the blue to NIR wavelength region ([Figure 3.12](#)).




3. Put a check mark next to Photograph.
4. Select a Field of View (FOV, size of the stage area to be imaged). See [Table 3.1 on page 25](#) for a list of FOV settings.



TIP: See the concept tech note *Detection Sensitivity* for more information about the Field of View (select **Help** → **Tech Notes** on the menu bar).

5. Select a focus option ([Figure 3.13](#)).

The focal distance to the camera is set at stage $z = 0$ for each field of view. To focus at the top of the animal, the stage moves down so that the top of the animal is at $z = 0$. You can enter the height of the animal and select the “use subject height” option or use the manual focus option to determine the proper subject height for the area to be imaged. See [Appendix A on page 272](#) for manual focus instructions.

Figure 3.13 Choose a Focus Option in the Control Panel

6. Select **Overlay** to view an overlay image (registered photograph and fluorescent image) after acquisition.



NOTE: If you want to check the subject inside the chamber before acquisition, take a photograph—uncheck the Fluorescent option, choose the Photograph option, and click **Acquire**. Be sure to check the Fluorescent option after taking the photograph.

7. Click **Acquire** when you are ready to capture the image.




NOTE: If necessary click **Image Setup** in the control panel to operate in single image mode. In single image mode, the **Sequence Setup** button appears in the control panel. Use this button to set up sequence acquisition (see page [page 39](#) for more details on sequence setup).

8. Enter information about the image in the Edit Image Labels box that appears (optional). Click **OK**.




NOTE: You can enter image label information at any time during or after acquisition. If you do not want to enter image information, click **Cancel**. See [page 70](#) for details on adding information to an image after acquisition.

Figure 3.14 Enter Information About the Image (optional)


Information entered here appears in the image label (see [Figure 3.16 on page 34](#)).

If this is the first image of the session, you are prompted to enable the autosave function (Figure 3.15). When Autosave is enabled, all images acquired during the session are automatically saved to a user-selected location. A different location can be chosen at any time (select **Acquisition** → **Auto-Save** on the menu bar).



- Click **Yes** in the prompt to enable autosave, then choose a location in the dialog box that appears. Alternatively, click **No** in the prompt and manually save the image data. See page 51 for details. Image acquisition begins and the upper area of the control panel changes to red color. During acquisition, the **Acquire** button in the control panel becomes a **Stop** button. Click **Stop** to cancel acquisition.

The control panel returns to blue color when acquisition is finished and the image window appears (Figure 3.16). See Table 3.2 on page 28 for details on the image window.



TIP: See the tech note *Determine Saturation* for information on pixel measurements (select **Help** → **Tech Notes** on the menu bar).

3.4 Fluorescent Imaging With Transillumination

Fluorescent imaging captures signals from fluorescent molecular reporters. Transillumination (excitation light source located below the stage) is recommended if the fluorescent source is deep relative to the imaged side of the animal.

Acquisition with transillumination includes a Normalized Transmission Fluorescence (NTF) Efficiency image in which the fluorescent emission image is normalized by the transmission image measured with the same emission filter and open excitation filter (Figure 3.17).




TIP: See these tech notes for helpful information and quick guides (select **Help** → **Tech Notes** on the menu bar):

- *Transmission Fluorescence*
- *Transmission Fluorescence – Raster Scan*
- *Transmission Fluorescence – Normalized Transmission Fluorescence*
- *Transmission Fluorescence – Well Plates*

Figure 3.17 Fluorescent Images Acquired with Transillumination

The NTF Efficiency image in this example highlights the presence of fluorescence in the animal, while the Radiant Efficiency image shows signal ambiguous with autofluorescence.



This section explains how to acquire a single fluorescent optical image with transillumination. See page 39 for information on acquiring a fluorescent sequence.

To acquire a fluorescent image with transillumination:




NOTE: Use only the Single Mouse Anesthesia Manifold when imaging with transillumination. The Dual Mouse or Five Mouse manifolds cannot be used with transillumination.

1. Put a check next to **Fluorescent** and **Transillumination** in the control panel.



NOTE: The Normalization option is selected by default so that NTF Efficiency images can be produced.



2. Select an excitation and emission filter from the drop-down lists.
The instrument has 18 narrow band excitation filters that span 490-850nm with a 20nm bandwidth, enabling spectral scanning over the blue to NIR wavelength region ([Figure 3.12 on page 32](#)).
3. Click Setup. Click **Yes** if prompted to acquire a subject photograph.
4. Choose the location (click a square) for transillumination and image acquisition in the Transillumination Setup box that appears ([Figure 3.19](#)).





Table 3.3 Transillumination Setup Box


Item	Description
Move Motors to Selected Spot	Transillumination motors will move the excitation light source to the grid location selected in the Transillumination Setup dialog box.
Mask Grid points To Subject	When setting up a transillumination sequence, choose this option to automatically select only the grid locations within the subject boundaries. Grid locations outside the subject are masked out. The mask prevents the transillumination excitation source from selecting an uncovered hole. Projecting light through an open hole would saturate the camera.
Raster Scan	If this option is not selected, the software generates one image per transillumination location per filter pair. For example, a sequence setup that includes 20 locations using two filters will generate 20 images. If the raster scan option is selected, the software takes all of the images from the transillumination locations and adds them together into one image. The raster scan option may be helpful when trying to determine the optimal excitation and emission filters for a particular fluorescent probe.
Grid Type	9x19 grid
Update Photograph	Click to acquire a new photographic image. If the chamber door is opened during transillumination setup, you are prompted to acquire a new photograph.
Clear Selections	Clears selected/ highlighted transillumination locations on the grid.

5. Confirm that the Lamp Level is set to High in the control panel.

 **NOTE:** The lamp may be set to Low for certain applications, such as long wavelength data through thin tissue.

6. Select a Field of View (FOV, size of the area to be imaged). See [Table 3.1 on page 25](#) for a list of FOV settings.
7. Select a focus option ([Figure 3.13](#)).

The focal distance to the camera is set a stage $z = 0$ for each field of view. To focus at the top of the animal, the stage moves down so that the top of the animal is at $z = 0$. You can enter the height of the animal and select the “use subject height” option or use the manual focus option to determine the proper subject height for the area to be imaged. See [Appendix A on page 272](#) for manual focus instructions.



8. Select **Overlay** to view an overlay image (registered photograph and fluorescent image) after acquisition.



NOTE: If you want to check the subjects inside the chamber before image acquisition, take a photograph (uncheck the Luminescent option, choose the Photograph and Auto options, and click **Acquire**).

9. Click **Acquire** when you are ready to capture the image.




NOTE: If necessary click **Image Setup** in the control panel to operate in single image mode. In single image mode, the **Sequence Setup** button appears in the control panel. Use this button to set up sequence acquisition (see [page 39](#) for more details on sequence setup).


10. Enter information about the image in the Edit Image Labels box that appears (optional) ([Figure 3.21](#)). Click **OK**.



NOTE: You can enter image label information at any time during or after acquisition. If you do not want to enter image information, click **Cancel**. See [page 70](#) for details on adding information to an image after acquisition.




If this is the first image of the session, you are prompted to enable the autosave function ([Figure 3.22](#)). When Autosave is enabled, all images acquired during the session are automatically saved to a user-selected location. A different location can be chosen at any time (select **Acquisition** → **Auto-Save** on the menu bar).



11. Click **Yes** in the prompt to enable autosave, then choose a location in the dialog box that appears. Alternatively, click **No** in the prompt and manually save the image data. See [page 51](#) for details.

Image acquisition begins and the upper area of the control panel changes to red color. During acquisition, the **Acquire** button in the control panel becomes a **Stop** button. Click **Stop** to cancel acquisition.

The control panel returns to blue color when acquisition is finished and the image window appears ([Figure 3.23](#)). See [Table 3.2 on page 28](#) for details on the image window.



TIP: See the tech note *Identify Saturated Pixels in an Image* for information on pixel measurements (select **Help** → **Tech Notes** on the menu bar).

3.5 Acquire a Sequence Using the Imaging Wizard

The acquisition parameters for each image in a sequence must be specified. The Imaging Wizard ([Figure 3.24](#)) provides a convenient way to set up a sequence for some imaging applications (see [Table 3.4 on page 41](#)). The wizard guides you through a series of steps, prompting you for the information that the software needs to set up the sequence.

This section explains how to start the Imaging Wizard and acquire a sequence of luminescent, fluorescent, or Cherenkov images. A sequence can also be set up manually (see [page 47](#) for details).



TIP: See the *Imaging Wizard* tech note for a quick guide (select **Help** → **Tech Notes** on the menu bar).


Starting the Imaging Wizard






NOTE: The IVIS Spectrum should be initialized and the temperature locked before setting imaging parameters. See [page 7](#) for more details.

1. Click **Imaging Wizard** in the control panel ([Figure 3.24](#)). If necessary, click **Restart** in the Imaging Wizard to show the first screen of the wizard.

Figure 3.24 Open the Imaging Wizard




2. Double-click an imaging mode: Bioluminescence , Fluorescence , or Cherenkov .
3. Double-click an imaging option in the next screen (see [Table 3.4 on page 41](#)).
4. Step through the rest of the wizard.

Each page of the wizard guides you with step-by-step instructions and descriptions. When you finish the wizard, it sets up the sequence to acquire ([Figure 3.25](#)).

Figure 3.25 Control Panel and Sequence Setup

Each row in the sequence table specifies the acquisition parameters for one image in the sequence. See [page 48](#) for details on the sequence table.



5. To clear the sequence, click the **Remove** button  and select **All**.



NOTE: DyCE analysis tools require a separate license.

Table 3.4 Imaging Wizard – Imaging Mode Options


Imaging Mode	Options	See Instructions on Page
Bioluminescence	Open Filter – Acquires a luminescent image at maximum sensitivity.	
	Spectral Unmixing – Acquires an image sequence for analysis using the Spectral Unmixing tools which use a mathematical operation to separate the signals from multiple luminescent probes.	128
	DyCE – Acquires a time series of optical images following a bolus injection of probe (radiotracer, bioluminescent, or fluorescent) to track probe biodistribution.	156
	DLIT – Acquires an image sequence for analysis with the DLIT algorithm that reconstructs the position, geometry, and strength of 3D luminescent sources.	186
Fluorescence	Filter Pair – Selects the best excitation and emission filters for a specific fluorescent probe. Detects the fluorescent signal on the surface of the subject.	
	Spectral Unmixing/Filter Scan – Acquires an image sequence for analysis with the Spectral Unmixing tools to: <ul style="list-style-type: none"> ■ Extract the signal of one or more fluorophores from the tissue autofluorescence. ■ Determine the optimum excitation and emission filter for a probe. 	128
	DyCE – Acquires a time series of optical images following a bolus injection of probe (radiotracer, bioluminescent, or fluorescent) to track probe biodistribution.	156
	FLIT (Fluorescence Imaging Tomography) – Acquires an image sequence for analysis with the FLIT algorithm that reconstructs the position, geometry, and strength of 3D fluorescent sources. This technique is only available in transillumination mode (light source below the stage).	197
Cherenkov	Open Filter – Acquires a Cherenkov image at maximum sensitivity.	
	Spectral Unmixing – Acquires an image sequence for analysis using the Spectral Unmixing tools which use a mathematical operation to separate the signals from multiple luminescent probes.	
	DyCE – Acquires a time series of optical images following a bolus injection of probe (radiotracer, bioluminescent, or fluorescent) to track probe biodistribution.	162

Acquire the Sequence


1. Confirm that the IVIS Spectrum is initialized and the CCD temperature is locked. (See [page 7](#) for details.)
2. Click **Acquire Sequence** in the control panel when ready to begin acquisition.
3. Enter information about the image in the Edit Image Labels box that appears (optional). Click OK ([Figure 3.26](#)).




NOTE: You can enter image label information at any time during or after acquisition. Click **Cancel** if you do not want to enter image information.



If this is the first image of the session, you are prompted to enable the autosave function ([Figure 3.27](#)). When Autosave is enabled, all images acquired during the session are automatically saved to a user-selected location. A different location can be chosen at any time (select **Acquisition** → **Auto-Save** on the menu bar).



4. Click **Yes** in the prompt to enable autosave, then choose a location in the dialog box that appears. Alternatively, click **No** in the prompt and manually save the image data. See [page 51](#) for details. Image acquisition begins and the upper area of the control panel changes to red color. During acquisition, the **Acquire** button in the control panel becomes a **Stop** button. Click **Stop** to cancel acquisition. The image window displays the images as they are acquired. The control panel returns to blue color when acquisition is finished and the Tool Palette appears ([Figure 3.28](#)).



The Image window may include multiple tabs, depending on the type of acquisition:

- Sequence View – Displays the image sequence.
- 3D View – Displays the 3D volume if the acquisition included DLIT or FLIT data.


 **TIP:** See the tech note *Saturated Pixels In an Image* for information on pixel measurements.

Table 3.5 Image Window – Sequence View

Item	Description
Units	Select the measurement units for the image display from this drop-down list. The available units depend on the type of image. See the concept tech note <i>Image Display and Measurement</i> for more details on measurement units.
Use Saved Colors	Choose this option to display the image data using the color table that was specified in the Preferences at the time of acquisition. If this option is not selected, image data are displayed using the color table currently specified in the Preferences.

Table 3.5 Image Window – Sequence View (continued)









Item	Description
Options	<p>Layout - Choose a display option for the images in a sequence (Default, Dynamic, or Film Strip). For example, here is Film Strip mode:</p>  <p>Sort by - Options for ordering images in the sequence window. This option only applies to images that were opened using the “Load as Group” function in the Living Image browser.</p> <p>Default - Order in which the images are stored in the folder.</p> <p>TimeStamp - Ascending order of the image acquisition time.</p> <p>UserID - Ascending alphanumeric order of the user ID</p> <p>Display - Choose the types of information to display with each image.</p>  <p>In this example, exposure time and binning factor are displayed on each image</p>
Info	Click to show or hide the image label information (Figure 3.28).
	Opens all of the images in the sequence.
	Closes all open images.
	Opens the Edit Sequence dialog box that enables you to add or remove images from the sequence.
	Enables you to export the active image as a graphic file (for example, .png, .dcm).

Table 3.5 Image Window – Sequence View (continued)


Item	Description
	Creates a preview picture (<i>snapshot</i>) of the image or thumbnails that the Living Image Browser displays when the data are selected. For more details on the browser, see page 52 .
	Preview picture of the selected data

3.6 Acquire Multiple Sequences in Batch Mode


Use batch mode to set up multiple, separate sequences which will be automatically acquired, one after another, without manual intervention.

To setup and acquire sequences in batch mode:

1. Click **Sequence Setup** in the control panel.
2. Choose the Batch Sequences option ([Figure 3.29](#)).

Figure 3.29 Control Panel

3. To set up the first sequence, do either of the following:
 - Click Imaging Wizard and step through the wizard (see [page 39](#) for details).
 - OR


- Set up the sequence manually (see [page 47](#) for details).
4. To set up the next sequence:
- If using the Imaging Wizard, repeat [step 3](#).
Each sequence is displayed in a separate tab.
 - If setting up the sequence manually, click the button  in the sequence table to add a new tab, then proceed with manual setup in the new tab.




NOTE: Sequence tabs can be renamed. Double-click a tab name to edit it. Alternatively, right-click the selected name to view a shortcut menu of edit commands (for example, Cut, Copy, Paste).

Figure 3.30 Multiple Sequence Tabs

Three sequences are specified in this example.



5. To remove a sequence, click the sequence tab and then click the  button.
6. Click **Acquire Sequence** when you are ready to capture the sequences.
Image acquisition proceeds with no intervening time delay between sequences.



NOTE: If the check mark is removed next to the Batch Sequences option in the control panel ([Figure 3.29](#)), only the sequence in the active tab will be acquired.

To save the batch sequence setup:

1. Click the Save button .
2. Enter a file name (.xsq) and choose a location for the file in the dialog box that appears.

3.7 Manually Set Up a Sequence

This section explains how to set up an image sequence if you do not use the Imaging Wizard. The sequence parameters in the sequence table can be saved as a Living Image Sequence Setup file (.xsq).

For details on image acquisition, see [Acquire the Sequence](#) on page 41.



TIP: It may be convenient to create an image sequence by editing a sequence setup generated with the Imaging Wizard or an existing sequence setup (.xsq). Save the modified sequence setup to a new name.




1. Click Sequence Setup in the control panel (Figure 4.30).
The sequence table appears.
2. If necessary, click the **Remove** button  and select **All** to clear the sequence table.

Figure 3.31 Open the Sequence Table



3. Choose a subject and probe from the drop-down lists (Figure 3.32)

Figure 3.32 Choose Subject and Probe



4. Specify the imaging settings for the first image in the sequence. (See [Appendix A on page 269](#) for details on the imaging parameters in the control panel.)



NOTE: If you selected Photograph and the photograph Reuse option in the control panel ([Figure 3.33](#)), the IVIS Spectrum acquires only one photograph for the entire sequence. If this option is not chosen, the system acquires a photograph for each image in the sequence.

5. Click the **Add** button .

The acquisition parameters appear in the sequence table ([Figure 3.33](#)).


6. Repeat [step 4](#) to [step 5](#) for each image in the sequence.
7. To set a time delay between each acquisition, enter a time (minutes) in the Delay box in the sequence table.
8. To save the sequence setup information (.xsq):
 - a. Click the **Save** button  in the sequence table.
 - b. Select a directory, enter a file name, and click **Save** in the dialog box that appears.

Figure 3.33 Control Panel and Sequence Table with Image Settings




Table 3.6 Sequence Table











Item	Description
	Starts the Imaging Wizard.
	Displays a dialog box that enables you to select and open a sequence setup (.xsq), sequenceinfo.txt, or clickinfo.txt file.
	Displays a dialog box that enables you to save the information in the sequence table to a sequence setup file (.xsq).
Display Photographic Settings	Choose this option to include the photograph exposure time, binning, and F/Stop in the sequence table.
Subject:  Mouse 	If a subject and probe are specified (optional), the software uses the information to automatically set parameters in the Surface Topography, DLIT, FLIT, Spectral Unmixing, and Planar Spectral Imaging tools. If a subject or probe is not selected here, the default parameters appear in the Tool Palette.
Number of Segments	The sequence specified in the sequence table is called a <i>segment</i> . Choose this option to set the number of segments to acquire and the time delay between segments. This is useful for acquiring data for kinetic analysis.
Delay	Specifies a time delay between each segment acquisition.

Table 3.6 Sequence Table (continued)


Item	Description
	Applies the selected cell value to all cells in the same column.
	Remove Selected - Deletes the selected row from the sequence table. Remove All - Removes all rows from the sequence table.
	Updates the selected row in the sequence table with the acquisition parameters in the control panel.
	Inserts a row above the currently selected row using the information from the control panel.
	Adds a new row at the end of the sequence setup list.

Editing Image Parameters

You can edit imaging parameters in the sequence table or in the control panel.


To edit a parameter in the sequence table:

- 1. Double-click the cell that you want to edit (Figure 3.34).




- 2. Enter a new value in the cell or make a selection from the drop-down list. To apply the new value to all of the cells in the same column, click .
- 3. Click outside the cell to lose focus.

To edit a parameter in the control panel:

- 1. Select the row that you want to modify in the sequence table.
- 2. Set new parameter values and/or imaging mode in the control panel.
- 3. Click  in the sequence table.


Inserting Images in a Sequence

Method 1:

- 1. Select the sequence table row that is below where you want to insert a new image (row).
- 2. Set the imaging mode and parameters in the control panel.
- 3. Click  to insert the new image above the selected row,

Method 2:


1. Select the row(s) of interest and right-click the sequence table to view a shortcut menu of edit commands (Figure 4.34 on page 57).

**Table 3.7** Sequence Table – Shortcut Menu Edit Commands

Command	Description
Copy row(s)	Copies the selected row(s) to the system clipboard.
Select All	Selects all rows in the sequence table.
Delete row(s)	Deletes the selected row(s) from the sequence table.
Replace Row(s)	Replaces the row(s) selected in the sequence table with the rows in the system clipboard. Note: The Replace function is only available when the number of rows in the system clipboard is the same as the number of rows selected in the sequence table.
Paste Row(s)	Adds copied rows to end of the sequence.

Removing Images From a Sequence

Method 1:

1. Select the row(s) that you want to delete.
2. Click  and choose **Selected** from the drop-down list.


Method 2:

Select the row(s) of interest and right-click the sequence table to view a shortcut menu of edit commands (Figure 3.35).

3.8 Manually Save Image Data

When you acquire the first image(s) of a session, you are prompted to enable the autosave feature. If autosave is enabled, all images acquired during the session are automatically saved to a user-selected location. You can choose a different location at any time (select **Acquisition** → **Auto-Save** on the menu bar).

This section explains how to manually save data if you do not want to use the autosave feature.

1. Turn off the autosave feature (select **Acquisition** on the menu bar and remove the check mark next to **Auto Save**).
2. After image or sequence acquisition, click the **Save** button . Alternatively, select **File** → **Save** on the menu bar.
3. Select a directory in the dialog box that appears, and click **OK**.




NOTE: The software automatically includes the user ID, and a data and time stamp with the data.

3.9 Exporting Images

The active image view can be saved in different file formats (for example, .bmp, .dcm).

1. Open an image or sequence.
2. Click the **Export Graphics** button  (Figure 3.36).



3. Select a directory in the dialog box that appears and enter a file name.
4. Click **Save**.



NOTE: To export a sequence to DICOM (.dcm) format, select **Export** → **Image/Sequence** as DICOM on the menu bar. This creates a directory that contains the .dcm files and a SequenceInfo.txt.

4 Working With Optical Image Data

Loading Optical Image Data

Adjusting Image Appearance on page 60

Correcting Optical Image Data on page 64

Viewing Intensity Data on page 65

Measuring Distance on page 68

Managing Image Information on page 69

Creating a Transillumination Overview on page 73

Overlaying Multiple Images on page 74

Rendering Intensity Data in Color on page 77

Exporting or Printing Images on page 79

Managing Image Sequences on page 81

4.1 Loading Optical Image Data

You can load (open) optical images from the:

- Living Image Browser (see below).
- Toolbar or menu bar ([page 56](#)).

Alternatively, drag an image file or sequence folder to the Living Image main window.

Multiple data sets can be open at the same time.




NOTE: Select **File** → **Recent Files** on the menu bar to view recently opened files.

Preview and Load Data Using the Living Image Browser

The Living Image Browser provides a convenient way to browse and preview optical data, view information about the data, and load the data.

1. Start the Living Image Browser:

- a. Click the **Browse** button . Alternatively, select **File** → **Browse** on the menu bar.
- b. In the dialog box that appears, select the folder of interest and click **OK**.

The Living Image Browser appears ([Figure 4.1](#)). It displays all Living Image data located in the folder and its subfolders, along with the user ID, label information, and camera configuration information.



NOTE: The next time you start the Living Image software and open the Browse For Folder box, the software automatically returns to the last folder visited.



Figure 4.1 Opening the Living Image Browser

Figure 4.2 Living Image Browser

To expand a sequence, click the arrow next to **SEQ.**

Click a column header to sort the browser contents in ascending alpha numeric order. Click the column header again to sort in descending alpha numeric order.

To view data properties, right-click a row and select **Properties** on the shortcut menu.



To preview data, click a row.

Note: A preview snapshot is automatically taken at the time of image or sequence acquisition. A snapshot can also be captured manually (see [page 59](#) for more details).


2. Load data by doing one of the following:

- Double-click the data row.
- Right-click the data name and select **Load** on the shortcut menu.
- Select the data row and click **Load**.
- Double-click the thumbnail.

The image(s) and Tool Palette are displayed. Green rows in the browser indicate loaded data ([Figure 4.3](#)).

Figure 4.3 Open (“Loaded”) Image Sequences

Multiple data sets can be loaded at the same time.

**Table 4.1** Living Image Browser

Item	Description
Hide Browse View	Closes the browser table.
Close Preview	Closes the image preview box.
Label Set	A drop-down list of the available label sets which specify image information (column headers) that is displayed in the Living Image Browser.
Add to List	If this option is chosen, the data selected in the Browse for Folder box is added to the Living Image Browser. If this option is not chosen, the data selected in the Browse for Folder box replaces the contents of the Living Image Browser, except for the loaded data.
Browse	Opens the Browse For Folder box.
View	The name of the Living Image Browser configuration (the column headers and their order in the browser).
Configure	<p>Opens a dialog box that enables you create and save custom Living Image Browser configurations.</p> <p>Note: To reorder a column in the browser, click the column header, then press the mouse key while you drag the header left or right. Release the mouse key to set the new position.</p>

Table 4.1 Living Image Browser (continued)

Item	Description
Load as Group	<p>Enables you to select particular images that you want to view as a sequence. The images may be acquired during different sessions.</p> <p>To select adjacent images in the browser, press and hold the Shift key while you click the first and last file in the selection.</p> <p>To select non-adjacent images in the browser:</p> <p>PC users: Press and hold the Ctrl key while you click the images in the browser</p> <p>Macintosh users: Press and hold the Cmd key (apple key) while you click the images in the browser.</p> <p>Note: The Load as Group option is only available when two or more images (non-kinetic) are selected in the browser.</p> <p>Tip: See the tech note Loading Groups of Images for a quick guide (select Help → Tech Notes on the menu bar).</p>
Load	Opens the selected image or image sequence.
Remove	Removes a user-selected image sequence(s) from the browser.
Close	Closes the Living Image Browser.

Load Data From the Menu Bar or Toolbar



NOTE: To open a recently viewed file, select **File** → **Recent Files** on the menu bar.

1. Click the **Open** button  on the toolbar. Alternatively, select **File** → **Open** on the menu bar.
2. In the box that appears, choose a file type filter from the drop-down list (Figure 4.4).


Figure 4.4 Opening Data From the Toolbar or Menu Bar

Table 4.2 File Filters

File Type Filter	Shows:
Living Image files	Click*.txt – an image (Living Image file format). Sequence*.txt – an image sequence (Living Image file format). *.dcm – kinetic data or an image that was exported to a DICOM file.
TIFF Image Files	Graphic files (*.tif, *.tiff).
All Files (*.*)	All file types.

3. Navigate to the file and click double-click it. Alternatively, select the data and click **Open**.

About the Image Window and Tool Palette

An image or image sequence is displayed in an image window (Figure 4.5). Multiple image windows can be open at the same time.




Table 4.3 Image Window

Item	Description
Units	Select the measurement units for the image display from this drop-down list. The available units depend on the type of image data. See the concept tech note <i>Image Display and Measurement</i> for more details on measurement units (select Help → Tech Notes on the menu bar).
Use Saved Colors (image sequence)	Choose this option to display the image data using the color table that was specified in the Preferences at the time of acquisition. If this option is not selected, image data are displayed using the color table currently specified in the Preferences.

Table 4.3 Image Window (continued)










Item	Description
Options (image sequence)	<p>Layout – Choose a display option for the images in a sequence (Default, Dynamic, or Film Strip). For example, here is Film Strip mode:</p>  <p>Sort by – Options for ordering images in the sequence window:</p> <ul style="list-style-type: none">■ Default – Order in which the images are stored in the folder.■ TimeStamp – Ascending order of the image acquisition time.■ UserID – Ascending alphanumeric order of the user ID. <p>Display - Choose the types of information to display with each image.</p>  <p>In this example, exposure time and binning factor are displayed on each image</p>
Labels	<p>Labels – Enables you to select the information to include in the image label.</p>  <p>Image label</p>
Info	<p>Click to show or hide the image label. The image label includes information you enter in the Edit Image Labels dialog box (see page 26) and other information automatically recorded by the software.</p>


Table 4.3 Image Window (continued)

Item	Description
	Opens all of the images in a sequence.
	Closes all open images of a sequence.
	Opens the Edit Sequence dialog box that enables you to add or remove images from the sequence display.
	Opens a dialog box that enables you to export the active view as a graphic file.
	Takes a “snapshot” that is displayed with the data in the Living Image Browser. See page 52 for more details on the browser.




Snapshots of an image sequence

The Tool Palette appears when you open an image or sequence ([Figure 4.6](#)). The options available in the Tool Palette depend on the type of active image data. A tool is only available if the data set includes the components that the tool requires to perform the analysis.



NOTE: The 3D Multi-Modality tools and DyCE tools require a separate license.




Organizing Images

When multiple image windows are open, you can organize them in a cascade or tile arrangement.

Choose **Window** → **Cascade** or **Window** → **Tile** on the menu bar.

Figure 4.7 Image Windows



4.2 Adjusting Image Appearance

Use the Image Adjust tools to adjust image display (Figure 4.8).



NOTE: Not all tools are available for all image display modes. Some tools are available for single images, but not image sequence and vice versa. For example, Corrections tools are available for an image, but not for an image sequence.

Figure 4.8 Tool Palette – Image Adjust Tools




Table 4.4 Image Adjust Tools














Item	Description
	Click this button to incrementally zoom out on the image (reduces the image dimensions in the image window). Note: The zoom tools are also available in the shortcut menu when you right-click the image (Cmd -click for Macintosh users).
	Click this button to incrementally zoom in on the image (incrementally magnifies the image in the image window).
	Draw a box on the image, then click this button to magnify the area inside the box. (Sets the dimensions of the magnified area equal to image window dimensions.) See Zooming or Panning on page 63 for more information on zooming.
	Click this button to draw a box on an image that can be used to: <ul style="list-style-type: none"> ■ Make measurements (see page 68) ■ Select an area of the image to copy to the system clipboard.
	Click this button to return the image to the default display magnification.
	Click this button to move a magnified image (<i>pan</i>) in the image window. For more details, see page 64 .
	Click this button to hide or display the image min/max information in the image window
	Click this button to hide or show the color scale in the image window.
Photo Adjustment	<p>Brightness – Click and move the slider left or right to adjust the brightness of an image displayed in overlay or photograph mode. Alternatively, enter a brightness value.</p> <p>Contrast – Click and move the slider left or right to adjust the <i>gamma</i> of an image displayed in overlay mode. Alternatively, enter a gamma value. (Gamma is related to image contrast.)</p> <p>Opacity – Click and move the slider left or right to adjust the opacity of the pseudocolor luminescent data of an image displayed in overlay mode. Alternatively, enter an opacity value.</p>
Color Scale	<p>Min – The minimum pixel intensity associated with the color scale for an image. Pixels less than the minimum value are not displayed.</p> <p>Max – The maximum pixel intensity associated with the color scale for an image. Pixels greater than the maximum value are displayed in the maximum color.</p>
Color Scale Limits	<p>Auto – If this option is chosen, the software sets the Min and Max values to optimize image display and suppress background noise. The Min and Max settings can be manually adjusted to further optimize the image display for your needs.</p> <p>Full – Choose this option to set the Max and Min values to the maximum and minimum data values in the image.</p> <p>Manual – Choose this option to enter Max and Min values for the image display.</p> <p>Individual – Applies a separate color scale limits to each image in a sequence. Note: This option is only available when an image sequence is active.</p>
Color Table	<div>  </div> <p>Click the drop-down arrow to select a color table for the image data. (For more details on color tables, see the concept tech note <i>Image Display and Measurement</i>.)</p> <p>Reverse – Choose this option to reverse the selected color table.</p> <p>Logarithmic Scale – Choose this option to apply a log scale to the relationship between numerical data and the color range in the color table. A log scale improves the visibility of dark areas in an image.</p>

Table 4.4 Image Adjust Tools (continued)


Item	Description
Smoothing	<p>Computes the average signal of a specified number of pixels (for example, 5 x 5) and replaces the original signal with the average signal (see Figure 4.9). Smoothing removes signal noise without changing pixel size. Smoothing can be applied to an image or a sequence.</p> <p> Click this button to return smoothing to the previous setting and update the image.</p>
Binning	<p>Specifies the number of pixels in the image data that are grouped together to form a larger pixel (called <i>soft binning</i>). Binning changes the pixel size in the image (see Figure 4.9). Binning can be applied to an image or a sequence. See the tech note <i>Detection Sensitivity</i> for more details on binning (select Help → Tech Notes on the menu bar).</p> <p> Click this button to return binning to the previous setting and update the image.</p>
Image X,Y	<p>The x,y pixel coordinates of the mouse pointer location in an image and the intensity (counts or photons) at that location. Note: This tool is only available when an image is active.</p>
	<p>Click this button to display a line profile (see page 66.) Note: This tool is only available when an image is active.</p>
	<p>Click this button to display the distance measurement cursor in the image window (see page 68). Note: This tool is only available when an image is active.</p>

Smoothing and Binning

Smoothing computes the average signal of a specified number of pixels (for example, 5 x 5) and replaces the original signal with the average signal ([Figure 4.9](#)). Smoothing removes signal noise without changing pixel size

Binning specifies the number of pixels in the image data that are grouped together to form a larger pixel (called soft binning) ([Figure 4.9](#)). Binning changes the pixel size in the image.

Smoothing and binning can be applied to a single image or all of the images in a sequence.



To set smoothing and/or binning:

1. Load an image or a sequence.
2. Make a selection from the Smoothing and/or Binning drop-down lists in the Image Adjust tools (Figure 4.10).

The image or all images in the sequence will be updated.



3. Click the  button to return the smoothing or binning to the previous setting and update the image or sequence images.

Figure 4.10 Select Smoothing and/or Binning Levels




Zooming or Panning


To incrementally zoom in or out on an image:

Click the  or  button. Alternatively, right-click the image and select **Zoom In** or **Zoom Out** on the shortcut menu.

To magnify a selected area in an image:

1. Click the  button. Alternatively, right-click the image and select **Area Zoom** on the shortcut menu.
2. When the pointer becomes a +, draw a rectangle around the area that you want to magnify.
The selected area is magnified when you release the mouse button.

To reset the magnification (remove magnification):

Click the  button. Alternatively, right-click the image and select **Reset Zoom** on the shortcut menu.

To pan the image window:

NOTE: Panning helps you view different areas of a magnified image. Panning is only available if the image has been magnified.

1. Click the button.
2. When the pointer becomes a , click and hold the pointer while you move the mouse.

4.3 Correcting Optical Image Data

Use the Corrections tools to subtract background or apply corrections to the optical image data.



TIP: See these technical notes for helpful information (select **Help** → **Tech Notes** on the menu bar).

- *Detection Sensitivity* (includes information about binning and smoothing).
- *Luminescent Background Sources and Corrections*.
- *Fluorescent Imaging* for more about fluorescent background.

Figure 4.11 Tool Palette – Corrections Tools




Table 4.5 Optical Image Data Corrections

Item	Description
Lens Distortion Correction	Select this option to correct for distortion at the perimeter of an image due to curvature of the CCD lens. Lens distortion correction is available for data acquired by Living Image® software version 4.3 and higher. The correction is particularly important for IVIS® Spectrum CT data acquired for DLIT or FLIT.
Read Bias Subtraction/Dark Charge Subtraction	<p>Select this check box to subtract dark background from the image data. If a dark charge image is available for the imaging conditions, the dark background image, including read bias noise, will be subtracted. Otherwise, only read bias noise will be subtracted.</p> <p>Note: In Radiance (Photons) mode, dark background or read bias subtraction is a mandatory default. In counts mode, the check box can be cleared.</p> <p>Tip: See the tech note <i>Luminescent Background Sources and Corrections</i> (select Help → Tech Notes on the menu bar).</p>
Flat Field Correction	<p>Select this check box to apply flat field correction to the image data.</p> <p>Note: In photons mode, flat field correction is a mandatory default. In counts mode, the check box can be cleared.</p>

Table 4.5 Optical Image Data Corrections (continued)

Item	Description
Cosmic Correction	Select this check box to correct image data for cosmic rays or other ionizing radiation that interact with the CCD. See the tech note <i>Image Data Display and Measurement</i> for more about cosmic correction (select Help → Tech Notes on the menu bar).
Adaptive FL Background Subtraction	Opens the Photo Mask Setup box that enables you to set the photo mask for adaptive fluorescent background subtraction. Tip: See the tech note <i>Adaptive Fluorescence Background Subtraction</i> (select Help → Tech Notes on the menu bar).
Normalization Threshold (Counts)	The minimum number of counts required to perform normalization.


4.4 Viewing Intensity Data

You can view intensity data:

- At a particular x,y location.
- Along a line drawn on the image.
- Within a user-selected region of interest (ROI). See [Chapter 5 on page 83](#) for more information on ROIs for optical data.

X,Y Coordinates and Intensity Data

1. Open an image and choose Cm or Pixels from the Units drop-down list in the Image Adjust tools.
2. Put the mouse pointer over the image to view the:
 - x,y pixel coordinates of the mouse pointer location in the image ([Figure 4.12](#)).
 - Intensity at the pixel location of the mouse pointer. The intensity is represented in the units currently selected for the image.

Figure 4.12 X,Y Coordinates and Intensity Data at the Mouse Pointer Location

Line Profile


The line profile plots intensity (y-axis) at each pixel (x-axis) along a user-specified line in the image. It is particularly useful for inspecting the detailed character of the image data.



NOTE: In the Overlay display mode, the line profile plots the luminescent data. To obtain a histogram of the photograph, select Photograph from the Display drop-down list.

To display a line profile:

1. Open an image and click the **Line Profile** button  in the Image Adjust tools (Figure 4.13). A line appears on the image and the Line Profile window opens. See Table 4.6 on page 67 for details on the Line Profile window.











2. To view the line profile at another location in the image, put the mouse pointer over the line. When the pointer becomes a , drag the line over the image. The line segment colored blue indicates the pixel intensities that are plotted in the line profile graph. The line profile is updated as you move the line move over the image.

Table 4.6 Line Profile Window

Item	Description
Line Orientation	Choose Vertical, Horizontal, or Free Hand from the drop-down list to set the orientation of the line in the image window. The Free Hand orientation enables you to drag each line segment endpoint to a user- selected position.
Width	Sets the line width. The Line Profile window displays the average of the pixel values included in the line width.
Position	Line position (pixels).
	Enables you to choose the grid line pattern to display in the line profile window.
	Exports the line profile data to a .csv or .txt file.
	Copies the line profile graph to the system clipboard.
	Opens the Print dialog box.
X Min X Max	Displays the minimum and maximum value of the x-axis. Use the  arrows to change the x-axis min or max. If a calibrated unit such as “radiance” is selected in the image window, the x-axis units = cm. If “counts” is selected in the image window, the x-axis units = pixels. To display the range available for the Min or Max, place the mouse pointer over the Min or Max edit box.
Y Min Y Max	Displays the minimum and maximum value of the y-axis. Use the  arrows to change the y-axis min or max. To display the range available for the Y Min or Y Max, place the mouse pointer over the Min or Max edit box.
	Click to reset the X and Y Min and Max values to the defaults.
Full Scale	Select this option to display the full X and Y-axis scales.
Logarithmic Scale	Select this option to apply a log scale to the y-axis.

4.5 Measuring Distance

Measure distance in an image using the distance measurement tool or image crop box.

Distance Measurement Tool

1. Open an image and in the Image Adjust tools, select Cm or Pixels from the Units drop-down list, then click the **Distance Measurement** button .



A measurement cursor () appears on the image (Figure 4.14). The Tool Palette displays the length of the cursor.

Figure 4.14 Measuring Distance Using the Measurement Cursor
The Tool Palette displays the measurement cursor position and length.




2. Change the cursor position or size by dragging the A or B end of the cursor to a new location on the image.
The measurement information in the Tool Palette is updated.
3. Click the  button to hide the cursor.

Image Crop Box






1. Open an image and select Cm or Pixels from the Units drop-down list in the Image Adjust tools.
2. Click the **Image Crop** button  in the Image Adjust tools (Figure 4.15).
The mouse pointer changes a "+".
3. Draw a rectangle on the image using the mouse pointer.
The length of the crop box diagonal is displayed.
4. Change the size or position of the crop box by dragging a  handle on the crop box.
5. Click the  button to hide the crop box.

Figure 4.15 Measuring Distance Using a Crop Box

4.6 Managing Image Information

Viewing Image Information


At acquisition, the software captures image information such as camera parameters and any image label information you entered at acquisition time (Figure 4.16).

Figure 4.16 Image Label Information

Detailed information about images is available in the View menu.

1. Open an image or sequence.
2. Select **View** → **Image Information** on the menu bar.
The Image Information window appears.
3. Choose an image by making a selection from the Sequences drop-down list and the Images drop-down list (Figure 4.17).

Figure 4.17 Viewing Image Information




4. To view particular information, select a category in the upper box to show the associated information in the lower box. For example, select luminescent image in the upper box to show the luminescent image acquisition parameters.

Editing the Image Label

You can edit image label information or add information to the label after acquisition.


To edit the image information:

1. Open an image or sequence.
2. Click **Info** to display the image label (Figure 4.18).




3. Edit the label information.

To add information to the image label:

1. Click the  toolbar button. Alternatively, select **Edit** → **Image Labels** on the menu bar.
2. In the Edit Image Labels box that appears, select information and/or enter a comment (Figure 4.19).



NOTE: If a single image is active, changes are applied to that image only. If a sequence is active, changes are applied to each image of the sequence.



3. Click **OK** when finished.
The image information is updated.
4. Save the image to save the updated image label (select **File** → **Save** or **File** → **Save As** on the menu bar).

Adding Comments

Comments can be added to an image and saved with the image.

To add comments:



1. Open an image.
2. Right-click the image and select **Insert Comment** on the shortcut menu. Enter comments in the yellow box that appears (Figure 4.20).
3. To move a comment in an image:
 - a. Position the mouse pointer over the comment.
 - b. When the hand tool appears , drag the comment box, then click the mouse to set the location.
4. Remove comments by doing either of the following:
 - Right-click a comment and select **Remove Comment** on the shortcut menu.
 - To remove all comments, right-click the image and select **Remove All Comments** on the shortcut menu.

Figure 4.20 Adding Comments




Tagging an Image


An image tag displays the x,y pixel coordinates and the pixel intensity (z, counts or photons) at a user-selected location.

To apply a tag:

1. Right-click a location in the image.
2. Select **Insert Tag** on the short cut menu.

Figure 4.21 Tag an Image (left), Move the Tag Label (right)



3. To move a tag:
 - a. Position the mouse pointer over the tag.
 - b. When the hand tool appears , drag the tag, then click the mouse to set the tag location. A line between the pixel and the tag identifies the location associated with the tag.

4.7 Creating a Transillumination Overview

The transillumination overview tool combines the images of a FLIT sequence (a fluorescence sequence acquired in transillumination mode) into a single image. All of the individual fluorescent signals are stacked over one photograph and the intensity is summed. One overview is created per filter pair. If two filter pairs were used during acquisition, then two overview images will be created.

All transillumination locations are displayed simultaneously; a tool tip displays the transillumination position when you mouse over a transillumination point. An overview image is displayed by default in radiant efficiency, and if transmission images are available, in normalized transmission fluorescence efficiency.


Transillumination overview images can be analyzed using the tools in the Tool Palette.



NOTE: If you choose the Raster Scan option in the Transillumination Setup box, the overview image is automatically generated (see [Figure 3.19 on page 36](#)).

1. Load a sequence that was acquired in fluorescence transillumination mode.
2. Click the **Overview** button. Alternatively, select **Tools** → **Transillumination Overview for <name>_SEQ** on the menu bar.
The overview appears.

Figure 4.22 Transillumination Overview



4.8 Overlaying Multiple Images


The image overlay tool provides a convenient way to view multiple reporters in one image. You can use the image overlay tool to display multiple luminescence or fluorescence images on one photographic image.



TIP: See the technical note *Image Overlay – 2D* for a quick guide (select **Help** → **Tech Notes** on the menu bar).

To coregister multiple images:

1. Acquire an image sequence using the appropriate filters for each reporter. Alternatively, create a sequence from images acquired during different sessions. (For more details, see [page 82](#).)
2. Load the image sequence.




3. Open one of the images and optimize the image display using the color scale Min and Max sliders in the Image Adjust tools.

To view all images in the sequence, click the **Display All** button  to open each image (overlay mode) in a separate image window.


4. Select **Tools**→ **Image Overlay for <sequence name>_SEQ** on the menu bar.

The image overlay window appears and shows the first photograph in the sequence. To view a different photograph, make a selection from the photograph drop-down list.



5. To overlay all images, click the  button.

The overlay appears. The photograph is at the bottom of the stack and the last fluorescent or luminescent image in the list is at the top of the stack.

**Table 4.7** Image Overlay Window





Item	Description
Units	Choose the type of units for displaying the fluorescent or luminescent image. See the concept tech note <i>Image Display and Measurement</i> for more details on measurement units.
Photograph	A drop-down list of the photographs in the image sequence.
Fluorescent or Luminescent Images	The sequence images.
	Copies the overlay to the system clipboard.
	Click to export the overlay to a graphic file.
	Click to include all fluorescent or luminescent images in the overlay.
	Click to remove all fluorescent or luminescent images from the photograph.
Image Adjust	<p>Tools for adjusting the appearance of the highlighted fluorescent or luminescent image. Adjustments can only be made on one image at a time.</p> <p>Min – The minimum pixel intensity associated with the color scale for an image. Pixels less than the minimum value are not displayed.</p> <p>Max – The maximum pixel intensity associated with the color scale for an image. Pixels greater than the maximum value are displayed in the maximum color.</p> <p>Opacity – Controls the opacity of the fluorescent or luminescent image.</p>

Table 4.7 Image Overlay Window (continued)


Item	Description
Color Table	<p>Tools for selecting and modifying the color scale associated with an image.</p> <p>Color Scale Type – Choose BlackLevel to show black at the low end of the color scale. Choose WhiteLevel to show white at the low end of the color scale.</p> <p>Red Click the drop-down arrow to select a color table for the image data. See the concept tech note <i>Image Display and Measurement</i> for more details on color tables (select Help → Tech Notes on the menu bar).</p> <p>Reverse – Choose this option to reverse the selected color table.</p> <p>Logarithmic – Choose this option to apply a log scale to the relationship between numerical data and the color range in the color table. A log scale improves the visibility of dark areas in an image.</p>
Palette label	<p>To include a brief line of text next to the color scale, enter text in the palette label box, then press the Enter key. To remove the text from the image window, delete the text in the palette label box and press Enter.</p>
Scales per Column	<p>Sets the number of color scales to display in a column.</p>

4.9 Rendering Intensity Data in Color

The colorize tool renders luminescence or fluorescence data in color, enabling you to see both intensity and spectral information in a single view. The tool provides a useful way to visualize multiple probes or scale probe signals that are not in the visible range.


To view colorized intensity data:

1. Load an image sequence.



2. Select **Tools** → **Colorize** on the menu bar.

The software renders each luminescent or fluorescent image in color and combines them into a single image (Figure 4.27).

**Table 4.8** Colorize Tools




Item	Description
Colorize View	
Color Map	NIR – A special camera setup that extends the color response into the near infrared range. Near infrared fluorophores appear red to purple using the NIR camera setup. VIS – Regular camera setup that mainly renders color in the visible range. It is similar to the color response of a commercial digital camera. NIR fluorophores appear dark red to invisible using the VIS camera setup.
Color Range	The color map indicates the color range of the selected camera setup from short to long wavelength. The two sliders determine the lower and upper limits of the color range that is used to render color. The parts of the color map outside the selected range are not used in the color rendering process. By default, the entire color range is selected.
Filter Range	The wavelength range of the luminescent images in the sequence. The two sliders determine the lower and upper end of the filter range. Only the parts of the image that are within the selected wavelength range are colorized. By default, the entire filter range is selected.
Log Scale	If this option is chosen, the dynamic range of the brightness in the image is compressed using a log scale. This improves the visibility of dark areas in the image.
Real Color	If this option is chosen, the colors are rendered using the wavelengths that directly correspond to the camera setup. For example, GFP appears green using real color rendering. If this option is not chosen, the original wavelength range of the image is modified to include the entire visible wavelength range of the camera setup. This helps improve the color contrast.
	Click this button to copy the colorize view to the system clipboard.
	Click this button to export the colorize view as a graphic file (for example, .jpg).



Table 4.8 Colorize Tools (continued)


Item	Description
	Click this button to print the colorize view.

4.10 Exporting or Printing Images

The Image Layout window (Figure 4.28) provides an alternative way to:

- Annotate and export an image (for example, .bmp)
- Print an image
- Copy an image to the system clipboard

1. Select **View → Image Layout Window** on the menu bar to open the Image Layout window.
2. Click the  button to paste the active image into the Image Layout window.
3. Drag a handle  at a corner of the image to resize the image.
4. Drag the image to reposition it in the window.

**Table 4.9** Image Layout Window


















Item	Description
	Clears the Image Layout window. Note: If you do not clear the layout (click the button) before you close the Image Layout window, the same window contents are displayed the next time the window is opened
	Opens a dialog box that enables you to save the Image Layout window contents to a graphic file.
	Pastes an image of the active data in the Image Layout window.
	Copies the contents of the Image Layout window to the system clipboard.

Table 4.9 Image Layout Window (continued)

Item	Description
	Pastes the contents of the system clipboard to the Image Layout window.
	Rectangle drawing tool
	Ellipse drawing tool
	Pointer tool
	Arrow and line drawing tool
	Select an item in the Image Layout window. To move the item to the front or back in the window, choose an option from the  drop-down list.
	Deletes the selected image.
	A drop-down list of formatting options for the Image Layout window. For example, the 2x2 layout style provides 4 separate layout areas in the window. A different image can be pasted into each layout area.
	To apply notes to an image, enter text in the annotation box and press Enter . Drag the text to the location of interest in the image.
	Opens a dialog box that enables you to select a font or edit the font style and size.
	Opens a color palette that enables you to select a font color or specify a custom font color.
	Opens a text editor that enables you to edit the selected text.

4.11 Managing Image Sequences

Editing a Sequence

You can add or remove individual images from a sequence. Only individual images, not an image sequence, can be added to a sequence.



1. Open the image sequence that you want to edit.
2. If you plan to add images to the sequence, browse for the images to add in the Living Image® browser. (See [page 52](#) for more details on browsing.)
3. Click the **Edit** button  in the image window ([Figure 4.29](#)).

Figure 4.29 Editing a Sequence



4. Choose the image(s) to add or remove (*retire*) from the sequence in the Edit Sequence box that appears ([Figure 4.29](#)).
To add an image to the sequence, select an image from the “Browser Images” and click **Copy**.
To remove an image from the sequence, choose an image from “Sequence Clicks” and click **Retire**.
5. To restore a retired image to the sequence, select the retired image and click **Reactivate**.
6. To reorder the sequence, select an image and click **Move Up** or **Move Down**.



NOTE: The **Move Up** and **Move Down** buttons are only available when the sequence view window displays images in the default sort order. If the TimeStamp or UserID sort order is selected, the images cannot be reordered.

7. Click **Close** when you are finished editing the sequence.
The updated image sequence is displayed.

Creating a Sequence From Individual Images

This section explains how to create a sequence from images acquired during different sessions.




TIP: Also see the tech note *Loading Groups of Images* for helpful information (select Help → Tech Notes on the menu bar).

1. In the Living Image Browser, browse for the images of interest. (See [page 52](#) for more details on browsing.)



NOTE: Browse for individual images (which may or may not be part of a sequence), not image sequences.

Figure 4.30 Living Image Browser



2. In the browser, select the images that you want to group together.
To select adjacent images in the browser, press and hold the **Shift** key while you click the first and last file in the selection.
To select non-adjacent images in the browser:
 - PC users - Press and hold the **Ctrl** key while you click the images of interest in the browser.
 - Macintosh users - Press and hold the **Cmd** key (apple key) while you click the images of interest in the browser.
3. Click **Load as Group**.
The image thumbnails are displayed together in an image window.
4. Save the images as a sequence:
 - a. Click the Save button . Alternatively, select **File** → **Save** on the menu bar.
 - b. In the dialog box that appears, select a folder and click **OK**.

5 ROI Tools for Optical Data

- About ROIs
- Quick Guide: Draw Measurement ROIs on an Optical Image or Sequence on page 86
- ROI Tools for Optical Images on page 87
- Measurement ROIs on page 89
- Mirror ROIs on page 94
- Measuring Background-Corrected Signal on page 97
- Managing ROI Properties on page 100
- Managing the ROI Measurements Table on page 109

5.1 About ROIs

This chapter explains how to draw and measure signal within a *region of interest* (ROI) on an optical image. Four types of ROIs are available for optical data (Table 5.2).

Table 5.1 Types of ROIs for Optical Images

ROI Name	Description	Shape	See Page
Measurement ROI for optical data	Measures the signal intensity in an area of an optical image.	Circle, square, grid, or contour	86 (Quick Guide) 89 (detailed steps)

Table 5.1 Types of ROIs for Optical Images (continued)




ROI Name	Description	Shape	See Page
Mirror ROI for left or right views of optical data obtained using the Side Imager	 <p>Images acquired using the Side Imager have three views: left, right, and center.</p> <p>Note: Use mirror ROIs to measure signal in the left or right views which are reflected from the mirrors. Use measurement ROIs to measure signal in the direct, non-reflected center view only.</p>	Circle or square	94
Average Background ROI for optical data	<p>Measures the average signal intensity in a user-specified area of an optical image that is considered background.</p> <p>Note: Using this type of ROI is optional. If the animal has significant autoluminescence or autofluorescence, you can determine a background-corrected signal in a measurement ROI by subtracting an average background ROI from a measurement ROI.</p> 	Circle or square	97


Table 5.1 Types of ROIs for Optical Images (continued)

ROI Name	Description	Shape	See Page
Subject ROI for optical data	Identifies a subject animal in an optical image. Note: Using this type of ROI is optional. It provides a convenient way to automatically associate (link) a measurement and average background ROI for background-corrected ROI measurements when there is significant autoluminescence or autofluorescence.	Square	97




Living Image software records information about the ROIs you create during a session and computes statistical data for the ROI measurements. The ROI Measurements table displays the data and provides a convenient way to review or export ROI information (Figure 5.1).

Figure 5.1 Example Measurement ROIs on an Optical Image and the ROI Measurements Table



5.2 Quick Guide: Draw Measurement ROIs on an Optical Image or Sequence


These steps provide a quick guide on how to apply a measurement ROI on an optical image or image sequence. See [page 89](#) for details on measurement ROIs.

1. Open an image or sequence and click ROI Tools in the Tool Palette.
2. Select Measurement ROI from the Type drop-down list.
3. Click the  button and select Auto All on the drop-down list.

The software automatically draws measurement ROIs on all images. The ROI label shows the total intensity in the ROI and the Threshold % ([Figure 5.2](#)).



NOTE: Auto ROIs are created and numbered in order from highest to lowest maximum signal within the ROI (ROI 1 contains the highest maximum signal). You may want to arrange the ROIs in a known order for easier comparison between images. To renumber the ROIs (ascending order from right to left), right-click the image and select Sort ROIs on the shortcut menu. If the "Apply to Sequence" option is selected in the ROI tools, choose "Sort ROIs in Sequence" to sort all of the ROIs in the sequence. The sort options are only available if the ROIs have not been sorted.



4. Use the Threshold % slider or  arrows to adjust the ROI boundaries.



NOTE: After the ROIs have been created, right-click an ROI to view a shortcut menu of ROI commands (Ctrl-click for Macintosh users). The shortcut menu provides easy access to many functions for managing ROIs and viewing ROI properties.

5. Click the **Measure** button  in the ROI tools to show the ROI Measurements table.

Figure 5.3 ROI Measurements Table

ROI Measurements

Refresh

Image Number	ROI	Image Layer	Total Count	Avg Count	Stdev Count	Min Count	Max Count
TLT20050624145507_001	ROI 1	Overlay	6.536e+03	2.334e+02	4.395e+01	1.723e+02	3.301e+02
TLT20050624145507_002	ROI 2	Overlay	2.653e+04	6.982e+02	1.402e+02	4.730e+02	9.393e+02
TLT20050624145507_003	ROI 3	Overlay	8.290e+05	2.961e+04	5.807e+03	2.209e+04	4.280e+04
TLT20050624145507_003	ROI 4	Overlay	2.687e+05	4.405e+03	9.302e+02	3.042e+03	6.048e+03
TLT20050624145507_004	ROI 5	Overlay	1.621e+06	4.053e+04	8.599e+03	2.976e+04	5.899e+04
TLT20050624145507_004	ROI 6	Overlay	7.273e+05	8.456e+03	1.784e+03	6.007e+03	1.199e+04

Customized Selections

Measurements Types: Counts

Image Attributes: _none_

ROI Dimensions: _none_

Copy

Select All

Configure...

Export...

Close

The ROI Measurements table displays data for all ROIs created in images or sequences during a session (one ROI per row). The table provides a convenient way to review and export ROI data. For more details on the table, see “ROI Measurements Table,” page 110.

6. Click **Yes** in the save prompt when closing a data set, to save the ROIs with the data.

5.3 ROI Tools for Optical Images

This section provides an overview of the ROI tools for optical images (Table 5.2). The ROI tools that appear in the Tool Palette depend on the type of ROI selected from the ROI Type drop-down list, and whether an image or sequence is active. Some ROI parameters are only available if “Show Advanced Options” is selected in the General Preferences (Figure 5.4).

Figure 5.4 ROI tools for Optical Images

Tool Palette

Image Adjust

ROI Tools

Measure ROIs

Apply to Sequence

Type: Measurement ROI

Save ROIs

Name: ROI_1_WLA

Delete

Load

Save

Auto ROI Parameters

Threshold %: 23

Lower Limit: 1.0

Minimum Size: 20

☒ Preview

☐ Use Bkg Offset

☐ Replace ROIs

Restore Defaults

Save / Load

Corrections

These Auto ROI parameters are available if “Show Advanced Options” is selected in the General Preferences. For more details on setting Preferences, see Appendix B on page 274.

Table 5.2 ROI Tools for Optical Images


Item	Description
	Click to select the number of circle ROIs to add to the active image.

Table 5.2 ROI Tools for Optical Images (continued)






Item	Description
	Click to select the number of square ROIs to add to the active image.
	Click to specify the grid pattern for a measurement ROI that you want to add to the active image. This tool is useful for an image of a multi-well culture plate or microplate.
	Click and select Auto All to automatically draw ROIs in the image using the auto ROI parameters. Click and select Auto 1 to automatically draw one ROI at a user-selected location using the auto ROI parameters. For more details on using the auto ROI features, see page 90 .
	Click to display the ROI Measurements table or compute intensity signal in an ROI.
	Click to display a drop-down list of options to delete an ROI(s) in the active image. For more details, see page 108 . Note: These commands do not delete the ROIs that are saved to the system (listed in the Menu Name drop-down list).
Apply to Sequence	Choose this option to apply the selected ROI to all images in a sequence.
Type	Choose the ROI type from the drop-down list: Measurement – Measures the signal intensity in an area of an image. Average Bkg – Measures the average signal intensity in a user-specified area of the image that is considered background. Subject ROI – Identifies a subject animal in an image. The software automatically associates a measurement and an average bkg ROI that are included in the same subject ROI. Using this type of ROI is optional. Mirror ROI – Measures the signal intensity in an area of an image acquired using the Side Imager, taking mirror reflection effects into account.
Save ROIs	Creates a file that includes the ROI parameters (for example, the X,Y coordinates, type of ROI, color, shape, width/height). ROIs that have been saved to file can be recalled and applied for another image at any time. Name – The name of the selected ROI set or the default name for a new ROI set. Delete – Deletes the selected ROI set from the system. Note: This permanently removes the ROI from the system. Load – Applies the ROI set selected from the Name drop-down list to the active image. Save – Saves the ROI set in the active image. Note: This is a global save (the ROI is saved to the system) and the ROI set can be loaded onto any image. If you use the File → Save commands to save an image that includes an ROI, the ROI is saved with the image only (not a global save) and is not available for loading onto other images. For more details, see Save, Load, or Delete ROIs, page 107 .
Auto ROI Parameters	Parameters that specify how the auto ROI tool draws an ROI. Threshold % – If the Auto All or Auto 1 method is selected, the Threshold % specifies the minimum percent of peak pixel intensity that a pixel must have to be included in an ROI identified by the software. After ROIs are drawn on an image, if you modify the Threshold% (move the slider or enter a new value), the software automatically updates the ROIs.

Table 5.2 ROI Tools for Optical Images (continued)

Item	Description
	<p>Note: The following Auto ROI parameters are only available if “Show Advanced Options” is selected in the general preferences. For more details on setting Preferences, see Appendix B on page 274.</p> <p>Lower Limit – Specifies a multiple (1 to 10) of the color scale minimum that sets the lower threshold for identifying an ROI. For example, if the lower limit = 2 and the color scale minimum = 1000 counts, then the auto ROI tool will only draw an ROI on areas of 2000 counts or greater. This helps create ROIs only within pixels visible on the image.</p> <p>Minimum Size – Sets the minimum size of an ROI (measured in pixels). For example if the minimum size is set at 50, then ROIs created on the image must be greater than 50 pixels in size.</p> <p>Preview – If this option is chosen, the software draws the ROI each time a parameter is changed. ROI parameters can be saved without drawing the ROI.</p> <p>Use Bkg Offset – Choose this option to measure background-corrected signal. This is typically used to remove natural animal background luminescence, and should not be confused with the dark-charge and read-bias 'background' corrections that are applied (by default) to the raw CCD data to remove electronic noise before any measurements. For more details, see page 97.</p> <p>Replace ROIs – If this option is chosen, all auto ROIs are replaced when new ROI(s) are created.</p> <p>Restore Defaults – Restores the factory-set defaults for the auto ROI parameters.</p>
Save/Load	Click to display or hide the tools that enable you to save, load, or delete auto ROIs in the active data. Note: The save function saves parameters, the not actual ROIs. This means that when you load saved auto ROI parameters, the software draws a new ROI using the saved values (Threshold%, Lower Limit, Minimum Size).

5.4 Measurement ROIs

This section explains in detail how to draw a measurement ROI on an optical image to obtain the intensity signal in a user-specified area. [Table 5.3](#) lists the three methods for drawing measurement ROIs on an image.



NOTE: See [page 86](#) for a quick guide to drawing measurement ROIs on an optical image or sequence.




Table 5.3 Methods for Drawing Measurement ROIs

Drawing Method	Description	See Page
Automatic	The software automatically locates and draws an ROI(s) on the image. To do this, the software locates the peak pixel intensities in the image and searches the neighborhood around a peak pixel. A pixel is included in the ROI if the pixel intensity is greater than the threshold%, a user-specified percentage of the peak pixel intensity.	90
Manual	Places one or more ROIs (circular, square, or grid shape) on the image.	86
Free draw	Draw line segments that define the ROI.	93

Automatically Draw ROIs

Living Image® software can automatically identify all of the ROIs in an image or image sequence that meet the auto ROI parameter thresholds or draw one ROI at a user-specified location.

To automatically identify and draw all ROIs:

1. Open an image or image sequence, and in the ROI tools, select Measurement ROI from the Type drop-down list (Figure 5.5).
2. Click an ROI shape button (Circle , Square , or Contour ) and select **Auto All** from the drop-down list.

The ROIs appear on the image or sequence thumbnails. The ROI label includes the ROI intensity threshold (Threshold%) and intensity measurement.



NOTE: Auto ROIs are created and numbered in order from highest to lowest maximum signal within the ROI (ROI 1 contains the highest maximum signal). You may want to arrange the ROIs in a known order for easier comparison between images. To renumber the ROIs (ascending order from right to left), right-click the image and select Sort ROIs on the shortcut menu. If the "Apply to Sequence" option is selected in the ROI tools, choose "Sort ROIs in Sequence" to sort all of the ROIs in the sequence. The sort options are only available if the ROIs have not been sorted.

3. Click the **Measure** button  in the ROI tools to show the ROI Measurements table.

Figure 5.5 Automatically Detect and Draw Measurement ROIs





Figure 5.6 ROI Measurements Table




The screenshot shows the 'ROI Measurements' window with a table of data for six ROIs. The table columns are: Image Number, ROI, Image Layer, Total Count, Avg Count, Stdev Count, Min Count, and Max Count. Below the table are dropdown menus for 'Measurements Types' (set to 'Counts'), 'Image Attributes' (set to '_none_'), and 'ROI Dimensions' (set to '_none_'). There are also buttons for 'Refresh', 'Copy', 'Select All', 'Close', 'Configure...', and 'Export...'.


Image Number	ROI	Image Layer	Total Count	Avg Count	Stdev Count	Min Count	Max Count
TLT20050624145507_001	ROI 1	Overlay	6.536e+03	2.334e+02	4.395e+01	1.723e+02	3.301e+02
TLT20050624145507_002	ROI 2	Overlay	2.653e+04	6.982e+02	1.402e+02	4.730e+02	9.393e+02
TLT20050624145507_003	ROI 3	Overlay	8.290e+05	2.961e+04	5.807e+03	2.209e+04	4.280e+04
TLT20050624145507_003	ROI 4	Overlay	2.687e+05	4.405e+03	9.302e+02	3.042e+03	6.048e+03
TLT20050624145507_004	ROI 5	Overlay	1.621e+06	4.053e+04	8.599e+03	2.976e+04	5.899e+04
TLT20050624145507_004	ROI 6	Overlav	7.273e+05	8.456e+03	1.784e+03	6.007e+03	1.199e+04


The ROI Measurements table displays data for all ROIs created in images or sequences during a session (one ROI per row). The table provides a convenient way to review and export ROI data. For more details on the table, see [“Managing the ROI Measurements Table,” page 109](#)

To automatically draw an ROI at a user-specified location:

1. Open an image.
2. Click an ROI shape button (Circle , Square , or Contour ) and select **Auto 1** from the drop-down list.

The create tool appears on the image.






3. Move the create tool using the ring  to the area where you want to draw the ROI, then click **Create**.

The ROI appears on the image and the ROI label displays the intensity signal.

4. To draw another ROI on the image, repeat step 2. to step 3.


For information on how to save ROIs, see [page 107](#).

Manually Draw ROIs

1. Open an image or image sequence, and in the ROI tools, select Measurement ROI from the Type drop-down list.
2. Select the ROI shape:
 - a. Click the **Circle** , **Square** , or **Grid**  button.
The grid shape is useful for drawing a grid of ROIs on an image of a well plate.
 - b. On the drop-down list that appears, select the number of ROIs that you want to add to the image or the grid ROI dimensions.
The ROIs and intensity measurements appear on the image.




NOTE: Manual ROIs are numbered in the order they are created. You may want to arrange the ROIs in a known order for easier comparison between images. To renumber the ROIs (ascending order from right to left), right-click the image and select Sort ROIs on the shortcut menu. If the "Apply to Sequence" option is selected in the ROI tools, choose "Sort ROIs in Sequence" to sort all of the ROIs in the sequence. The sort options are only available if the ROIs have not been sorted.




3. Adjust the ROI position:
 - a. Place the mouse pointer over the ROI. When the pointer becomes a , click the ROI.
 - b. Drag ROI(s).



NOTE: To move multiple ROIs at the same time, press and hold the Shift key while you click the ROIs, and then drag them to a new location. Contour ROIs () cannot be moved using this method.

Figure 5.8 Two Circle ROIs Placed on Image



4. Adjust the ROI dimensions:
 - a. Place the mouse pointer over the ROI. When the pointer becomes a , click the ROI.
 - b. Place the mouse pointer over an ROI handle  so that it becomes a . Drag the handle to resize the ROI.



NOTE: You can also change the ROI position or size using the adjustment controls in the ROI Properties box (see [Moving an ROI](#), page 103 and [Editing ROI Dimensions](#), page 104).

5. Click the **Measure** button .

The ROI measurements and table appear. For more details on the table, see “[Managing the ROI Measurements Table](#),” page 109.

For information on how to save ROIs, see page [page 107](#).

Draw ROIs Using the Free Draw Method









1. Open an image, and in the ROI tools, select the type of ROI that you want to draw from the Type drop-down list.
2. Click an ROI shape button (**Circle** , **Square** , or **Contour** ) and select **Free Draw** from the drop-down list. In this example, the Contour shape  was selected for the free draw method. The ROI shapes that are available depend on the type of ROI selected.
3. If you selected:
 -  or  - Use the pointer (+) to draw the ROI.
 -  - Use the pointer (+) to click around the area of interest and draw line segments that define the ROI. Right-click when the last point is near the first point in the ROI.

Figure 5.9 Drawing an ROI – Free Draw Method



4. Click the **Measure** button .

The ROI measurements and table appear. For more details on the table, see “[Managing the ROI Measurements Table](#),” page 109.

For information on how to save ROIs, see page [page 107](#).

5.5 Mirror ROIs

Use a mirror ROI to measure bioluminescence or fluorescence in the right or left mirror-reflected view of images acquired using the Side Imager. Measure signals in the center view using a measurement ROI. See [page 89](#) for more details on drawing a measurement ROI.




NOTE: Do not apply mirror ROIs on the center view or measurement ROIs on the left or right mirror-reflected views. Placing an ROI on the wrong view will result in incorrect ROI measurements.

1. Open an image or image sequence acquired with the Side Imager.



NOTE: Fluorescent image data acquired in reflectance/epi-illumination mode must include a photograph.

2. Select "Mirror ROI" from the Type drop-down list in the ROI tools. If analyzing a fluorescent image, choose the Photo Mask option.





3. Select the ROI shape:
 - a. Click the **Circle** , or **Square** , button.
 - b. Select the number of ROIs to add to the image on the drop-down list that appears.
 - If analyzing a reflectance/epi-illumination fluorescent image, go to [step 4](#); otherwise, go to [step 5](#).
4. For reflectance/epi-illumination fluorescent images only:
 - a. Confirm the purple data mask in the dialog box that appears ([Figure 5.11](#)).
The data mask includes the entire subject by default and defines the area of excitation light projection onto the animal. If you do not want to analyze the entire subject, select the Data Mask option and mask a particular area using the data mask options ([Table 5.4](#)).
 - b. Click **OK**.
The mirror ROIs and intensity measurements appear on the image ([Figure 5.12](#)).

Figure 5.11 Excitation Projection Setup Dialog Box
For fluorescent images only.




Table 5.4 Data Mask Options



Option	Description
Photograph	If this option is chosen, the software automatically draws the data mask by using higher intensities in the photograph. The mask selects high-valued photograph image pixels which are located continuously and centrally in the photograph image. The photograph mask works best with light-colored subjects.
Threshold	If necessary use the threshold slider or  arrows to adjust the mask so that it matches the underlying subject photograph as closely as possible without including any area outside the subject image.
Draw Mask	Choose this option to manually draw a data mask on an area of the photograph.
Rectangle	Specifies a rectangular shape for the manual data mask.
Ellipse	Specifies an elliptical shape for the manual data mask.




Figure 5.12 Mirror ROIs on Fluorescent Image Acquired with Side Imager

NOTE: The ROIs are numbered in the order they are created. You may want to arrange the ROIs in a known order for easier comparison between images. To renumber the ROIs (ascending order from right to left), right-click the image and select Sort ROIs on the shortcut menu. If the "Apply to Sequence" option is selected in the ROI tools, choose "Sort ROIs in Sequence" to sort all of the ROIs in the sequence. The sort options are only available if the ROIs have not been sorted.

5. Adjust the ROI position:
 - a. Place the mouse pointer over the ROI. Click the ROI when the pointer becomes a .
 - b. Drag ROI(s).



NOTE: To move multiple ROIs at the same time, press and hold the Shift key while you click the ROIs, and then drag them to a new location. Contour ROIs () cannot be moved using this method.

6. Adjust the ROI dimensions:
 - a. Place the mouse pointer over the ROI. Click the ROI when the pointer becomes a .
 - b. Place the mouse pointer over an ROI handle  so that it becomes a . Drag the handle to resize the ROI.



NOTE: You can also change the ROI position or size using the adjustment controls in the ROI Properties box (see [Moving an ROI, page 103](#) and [Editing ROI Dimensions, page 104](#)).

7. Click the **Measure** button .

The ROI table appears. For more details on the table, see [“Managing the ROI Measurements Table,” page 109](#).

5.6 Measuring Background-Corrected Signal

If a subject has significant autoluminescence or autofluorescence, you can obtain a background-corrected ROI measurement by subtracting an average background ROI from a measurement ROI. The software computes:

Background-corrected intensity signal = Signal in the measurement ROI - Average signal in the average background ROI




NOTE: This is an optional "background" correction that is applied in addition to the electronic dark-charge and read-bias corrections that are applied to the raw CCD data.

The Image Adjust tools and zoom feature are helpful for selecting an appropriate area for an ROI. By setting the image minimum close to zero and zooming in on a background area in the image, you can determine where naturally occurring background luminescence or autofluorescence is present. For more details on the Image Adjust tools and the zoom feature, see [Viewing Intensity Data, page 65](#) and [Zooming or Panning, page 63](#).

Draw a Subject ROI

A subject ROI identifies a subject animal in an image. It provides a convenient way to automatically associate (link) a measurement and average background ROI for background-corrected ROI measurements when there is significant autoluminescence or autofluorescence. Using a subject ROI is optional.


To draw a subject ROI using the auto ROI feature:

1. Select Subject ROI from the Type drop-down list.
2. Click the  button.
3. Select **Auto All**.



To manually draw a subject ROI:



NOTE: If the image was acquired using the Side Imager, draw three subject ROIs, one for each view.

1. Select Subject ROI from the Type drop-down list.
2. Click the  button, and select **1**.
3. Position the subject ROI so that it includes the measurement ROI(s) and the associated average background ROI.


Measure Background-Corrected Signal

1. Draw one or more measurement ROIs on the subject (see [page 89](#) for more details).
2. Draw an average background ROI on the subject.
 - a. Select Average Bkg ROI from the Type drop-down list.
 - b. Click the **Square**  or **Circle**  button and select **1**.

The ROI is added to the image. For more details on adjusting the ROI position or dimensions, see [page 103](#) and [page 104](#).




NOTE: The average background ROI and measurement ROI do not need to be the same shape or size because the software computes the average intensity signal in each ROI.



NOTE: If the image was acquired using the Side Imager, draw a background ROI on each view

Figure 5.13 Draw a Background ROI on Each Image View Acquired Using the Side Imager



3. Associate each background ROI with a measurement ROI(s) or mirror ROI(s) using one of the methods in Table 5.5.

Table 5.5 Methods for Associating Measurement or Mirror ROIs With a Background ROI


Method	Description
Draw a subject ROI See page 97 for details.	

Table 5.5 Methods for Associating Measurement or Mirror ROIs With a Background ROI (continued)





Method	Description
Right-click a measurement ROI and select an average background ROI from the shortcut menu.	
<div>1. Right-click a background ROI and select Properties on the shortcut menu.</div> <div>2. In the ROI Properties box that appears, click the Background ROI tab and put a check mark next to Use as BKG for future ROIs in.</div> <div>3. Choose the image name or the Entire sequence option.</div>	<div></div> <div></div>

Figure 5.14 Measurement ROIs Successfully Associated with a Background ROI



5.7 Managing ROI Properties

In the ROI Properties box, you can view information about an ROI, change the position of the ROI on the image, and edit the ROI label or line characteristics.

Viewing ROI Properties

- To view ROI properties, do one of the following:
 - Double-click an ROI in the image.
 - Right-click the ROI and select **Properties** from shortcut menu that appears.
 - Select the ROI, then select **View** → **Properties** on the menu bar.

The ROI Properties box appears (for more details see [Figure 5.17](#)).
- To view properties for another ROI, click the ROI in the image. Alternatively, select an ROI from the ROI drop-down list in the ROI Properties dialog box ([Figure 5.15](#)).





Figure 5.16 ROI Properties – Background ROI Tab

The items in the ROI Properties box depend on the type of ROI selected in the image. For more details see [Table 5.6](#), [page 102](#).








Table 5.6 ROI Properties

Item	Description
ROI	A drop-down list of ROIs in the active image or image sequence. To select an ROI, double-click the ROI in the image or make a selection from the drop-down list. Shape – The shape of the ROI (circle, square, grid, or contour) selected in the image. Type – Indicates the method that was used to draw the selected ROI (automatic, manual, or free draw).
ROI Label	Click to edit the selected ROI label name.
Image Number	A drop-down list of open images.
Background ROI tab	The Background ROI tab shows a drop-down list shows all average background ROIs in active image that can be linked to a user-specified measurement ROI or subject ROI (selected from the drop-down list at the top of the dialog box).

Table 5.6 ROI Properties (continued)

Item	Description
Subj ROI	The Subject ROI tab shows a drop-down list of all subject ROIs in the image number selected above that can be linked to a user-specified measurement ROI or average background ROI (selected from the drop-down list at the top of the dialog box). The Background ROI tab shows a drop-down list of all average background ROIs in the click number selected above that can be linked to a user-specified measurement ROI or subject ROI (selected from the drop-down list at the top of the dialog box).
ID	User-entered information about a subject ROI.
Label	Label name of the selected subject ROI.
Lock Position	Choose this option to lock the position of the ROI selected in the image.
Xc	x-axis coordinate at the center of the ROI selected in the image.
Yc	y-axis coordinate at the center of the ROI selected in the image.
Lock Size	Choose this option to lock the dimensions of the ROI selected in the image.
Width	Width (pixels or cm) of the ROI selected in the image (for more details on setting the units, see ROI Dimensions, page 110).
Height	Height (pixels or cm) of the ROI selected in the image.
Line Size	Specifies the ROI line thickness. To change the line thickness, enter a new value or click the up/down arrows  .
Line Color	Specifies the color of the ROI line. To select a line color, click the Browse button  .
Done	Click to close the ROI Properties box and apply any new settings, including: <ul style="list-style-type: none"> ■ Linkage between a measurement ROI and subject ROI (for more details, see Draw ROIs Using the Free Draw Method, page 93). ■ ROI size dimensions or position ■ Subject ROI ID information

Moving an ROI


To move an ROI on an image, select the ROI and do one of the following:

- Press a keyboard arrow key
- Drag the ROI
- Edit the settings in the ROI Properties box




NOTE: An ROI cannot be moved if it was created using the auto ROI tool or if the ROI position is locked.

To drag an ROI:

1. Put the mouse pointer over the ROI so that it becomes a  arrow.
2. Drag the ROI.
3. Release the mouse button when the ROI is properly positioned.

To move an ROI using the ROI Properties dialog box:

1. Double-click the ROI in the image.
The ROI Properties box appears and displays the position and dimensions of the selected ROI.



2. To set ROI position, enter new coordinates for the center of the ROI (Xc (pix or cm) and Yc (pix or cm values)) in the ROI Properties box.
3. To rotate the ROI clockwise, enter the degrees in the Angle (deg) box and click outside the box.
4. To lock the current ROI position, choose the Lock Position option.



NOTE: The ROI position cannot be changed until the Lock Position option is cleared.

Editing ROI Dimensions

There are two ways to resize a circle or square ROI:

- Drag a handle on the ROI
- Edit the settings in the ROI Properties box



NOTE: You cannot change the size of an ROI that was created using the auto ROI or free draw tool.


To resize an ROI using a handle:

1. Select the ROI and put the mouse pointer over a handle (■) on the ROI.
2. When the pointer becomes a ↖ arrow, drag the handle.

To resize an ROI using the ROI Properties box:

1. Double-click the ROI in the image.

The ROI Properties box appears and displays the positions and dimensions of the selected ROI.



Dimensions of the ROI
selected in the image


2. Enter a new width or height value in the ROI Properties box.
3. To lock the current ROI size, choose the Lock Size option.





NOTE: The ROI size cannot be changed until the Lock Size option is cleared.

Editing the ROI Line


1. Double-click the ROI that you want to edit.
The ROI Properties box appears (Figure 5.20).




2. To edit the ROI line thickness, enter a new value in the Line Size box. Alternatively, click the .
3. To change the ROI line color:
 - a. Click the **Browse** button .
The Select Color box appears.
 - b. To select a basic color for the ROI line, click a basic color swatch, and click **OK**.
 - c. To define a custom color, drag the crosshairs in the custom color field, adjust the brightness slider, and click **Add to Custom Colors**.
 - d. To select a custom color for the ROI line, click a custom color swatch, and click **OK**.

Move or Edit the ROI Label

To move the ROI label:

1. Put the mouse pointer over the ROI label.
2. When the pointer becomes a , drag the label, and then click to release the label at the new location (Figure 5.21).



To edit the ROI label:


1. Double-click the ROI of interest. Alternatively, right-click the ROI (**Ctrl**-click for Macintosh users) and select Properties on the shortcut menu.
2. In the ROI Properties box that appears, edit the name in the ROI Label box and click **Done** (Figure 5.21).

Save, Load, or Delete ROIs

The software automatically saves ROIs with an image. The ROI measurements are saved in the AnalyzedClickInfo.txt file associated with the image. ROIs are saved per user and can be applied to other sequences. Additionally, ROI parameters can be saved per user and applied to other sequences.

To save ROIs to the system:

1. In the Name drop-down list, confirm the default name or enter a new name for the ROI(s).

Figure 5.22 Name and Save ROIs to the System

2. Click **Save**.

The ROI(s) from the image are saved to the system and can be selected from the Name drop-down list.

To load ROIs on an image:

1. Open an image.
2. In the ROI tools, make a selection from the Name drop-down list and click **Load**.



NOTE: If you load ROI(s) onto an image, then draw additional ROIs, the Save button changes to Overwrite. If you want to save this collection of ROIs using the existing name, click Overwrite.

To delete ROIs from an image:



NOTE: This does not delete ROIs saved to the system (global save).

- Select the ROI and press the Delete key.

OR




- Click the  button in the ROI tools and select a delete command from the drop-down list.

Figure 5.23 Delete ROIs From an Image

To permanently delete ROIs from the system:

1. Select the ROI(s) that you want to delete from the drop-down list of saved ROIs.
2. Click **Delete**.

Figure 5.24 Delete ROIs From the System

5.8 Managing the ROI Measurements Table

The ROI Measurements table shows information and data for the ROIs created during a session. The ROI measurements can be displayed in units of counts, radiance, Radiant Efficiency, Efficiency, or NTF Efficiency, depending on the type of image data. See the technical note *Quantifying Image Data* for more details (select **Help** → **Tech Notes** on the menu bar).

Viewing the ROI Measurements Table


Click the  **Measure ROIs** button to display the ROI measurement table. Alternatively, select **View** → **ROI Measurements** on the menu bar.

Figure 5.25 ROI Measurements Table

Image Number	ROI	Image Layer	Total Count	Avg Count	Stdev Count	Min Count	Max Count	ROI Pixels	Area [ccd Pixels]	Xc [pixels]	Yc [pixels]	Width [pixels]
HX20070420121444_001	ROI1	Overlay	1.425e+06	4.167e+02	6.726e+02	-8.462e+00	2.545e+03	3420	2.189e+05	4.141e+02	1.020e+03	2.876e+02
HX20070420121444_001	ROI2	Overlay	1.489e+06	3.309e+02	6.061e+02	-8.642e+00	2.545e+03	4500	2.880e+05	4.168e+02	9.001e+02	2.822e+02
HX20070420121444_001	ROI3	Overlay	1.676e+06	3.637e+02	5.988e+02	-1.056e+01	2.545e+03	4608	2.949e+05	4.168e+02	1.152e+03	2.845e+02
HX20070420121444_002	ROI1	Overlay	7.261e+06	2.123e+03	3.445e+03	-2.181e+01	1.031e+04	3420	2.189e+05	4.141e+02	1.020e+03	2.876e+02
HX20070420121444_002	ROI2	Overlay	7.572e+06	1.683e+03	3.107e+03	-2.181e+01	1.031e+04	4500	2.880e+05	4.168e+02	9.001e+02	2.822e+02

Table 5.7 ROI Measurements Table

Item	Description
Measurement Types	Make a selection from the this drop-down list to select the type of image unit for the ROI measurements in the table.
None	Excludes ROI measurements from the table.
Counts (luminescence)	<p>Includes Total Counts, Avg Counts, Stdev Counts, Min Counts, and Max Counts in the table.</p> <p>Total Counts = the sum of all counts for all pixels inside the ROI.</p> <p>Avg Counts = Total Counts/Number of pixels or super pixels.</p> <p>Stdev Counts = standard deviation of the pixel counts inside the ROI</p> <p>Min Counts = lowest number of counts in a pixel inside the ROI.</p> <p>Max counts = highest number of counts in a pixel inside the ROI.</p> <p>Note: These numbers are displayed if the units selected in the ROI Measurements table and the image are the same. Otherwise, N/A appears in each column.</p> <p>Tip: See the tech note <i>Image Display and Measurement</i> for more details on count units (select Help → Tech Notes on the menu bar).</p>
Radiance (Photons) (fluorescence)	<p>Total Flux (photons/sec) = the radiance (photons/sec/cm²/steradian) in each pixel summed or integrated over the ROI area (cm²) x 4π.</p> <p>Average Radiance = the sum of the radiance from each pixel inside the ROI/ number of pixels or super pixels (photons/sec/cm²/sr).</p> <p>Stdev Radiance = standard deviation of the pixel radiance inside the ROI</p> <p>Min Radiance = lowest radiance for a pixel inside the ROI.</p> <p>Max Radiance = highest radiance for a pixel inside the ROI.</p> <p>Tip: See the tech note <i>Image Display and Measurement</i> for more details on photon units (select Help → Tech Notes on the menu bar).</p>
Radiant Efficiency (fluorescence)	<p>Epi-fluorescence - Fluorescence emission radiance per incident excitation intensity: p/sec/cm²/sr/μW/cm²</p> <p>Transillumination fluorescence - Fluorescence emission radiance per incident excitation power: p/sec/cm²/sr/mW</p>
Efficiency (epi-fluorescence)	Fluorescent emission yield normalized to the incident excitation intensity (radiance of the subject/illumination intensity)
NTF Efficiency (transillumination fluorescence)	Fluorescent emission image normalized by the transmission image which is measured with the same emission filter and open excitation filter.
Image Attributes	Make a selection from the drop-down list to specify the click number (image file) information to include in the table. Click attributes include label name settings and camera settings.
None	Excludes image attributes from the table.
All Possible Values	Includes all of the image attributes (for example, label name settings and camera settings) in the table.
All Populated Values	Includes only the image attributes with values in the table.
Living Image Universal	Includes all Living Image Universal label name settings in the table.
ROI Dimensions	Make a selection from the drop-down list to specify the ROI dimensions to include in the table.
None	Excludes the ROI area, x,y-coordinates, and dimensions from the table.
Pixels	Includes ROI area, x,y-coordinates, and dimensions (in pixels) in the table.


Table 5.7 ROI Measurements Table (continued)

Item	Description
cm	Includes ROI area, x,y-coordinates, and dimensions (in cm) in the table.
Copy	Copies the selected row(s) in the table to the system clipboard.
Select All	Copies all rows in the table to the system clipboard.
Refresh	Updates the ROI Measurements table (for example, after you draw new ROIs, move an ROI, and close or open image data).
Configure	Displays the Configure Measurements box that enables you to specify and organize the data categories (column headers) for the table.
Export	Displays the Save Measurements box so that the data can be saved to a .txt or .csv file. Note: Grid ROI measurements exported to a .csv file can be opened in a spreadsheet application like Microsoft® Excel®.
Close	Closes the ROI Measurements table.

Configuring the ROI Measurements Table

You can customize the data and information (column headers) in the ROI Measurements table. Several preset categories are available in the Measurement Types, Click Attributes, and ROI Dimensions drop-down lists.

1. Drag a column header (left or right) in the table to reorder the columns.
2. Make a selection from the Measurement Types drop-down list to change the measurement units.



3. Make a selection from the Image Attributes drop-down list to include image information in the ROI table.
4. Select units (Pixels or cm) from the ROI Dimensions drop-down list to include ROI dimensions in the table.

Creating a Custom ROI Table Configuration


A table configuration specifies the column headers in the ROI table. Several preset configurations are available (selected from the Measurements Types drop-down list in the ROI table, [Figure 5.26](#)). You can also create a custom table configuration.



NOTE: Preset table configurations cannot be edited. You can modify a preset configuration and save it to a new name.

1. In the ROI Measurements table, click **Configure**.
The Configure Measurements box appears.

Figure 5.27 Configure Measurements Dialog Box



2. Select a configuration from the User Lists drop-down list and click **Customize**.
3. To add column header to the ROI table, make a selection from the “Available Item” list and click **Add**.
4. To remove column header from the ROI table, select the item that you want to remove in the Selected Items list, and click **Remove**.
5. To reorder an item in the Selected Items list, select the item and click **Move Up** or **Move Down**.
The columns in the ROI Measurements table are updated.
6. Enter a name for the custom configuration in the Name box and click **Save**.

To delete a custom table configuration:

Select the configuration from the User Lists drop-down list and click **Delete**.




NOTE: Preset table configurations cannot be deleted.

Copying or Exporting the ROI Measurements Table

To copy the table to the system clipboard:

- Copy selected rows – Select the rows of interest and click **Copy**. Alternatively, select the rows, then right-click the table and choose Copy on the shortcut menu.
- Copy all rows – Click Select **All** and click **Copy**. Alternatively, press **Ctrl+A**, then right-click the table and choose **Copy** on the shortcut menu.

Figure 5.28 Copy All Rows in the ROI Measurements Table to the System Clipboard



To export the table:

1. Click **Export** in the ROI Measurements table.
2. In the dialog box that appears:
 - a. Select a folder and enter a name for the file.
 - b. Select a file type (.txt or ,csv) and click **Save**.

6 3D ROI Tools for Volumetric Data

About 3D ROIs

Drawing a 3D ROI on page 115

Managing the 3D ROI Measurements Table on page 120

6.1 About 3D ROIs


A 3D *region of interest* (ROI) can be drawn on a:

- DLIT reconstruction of a luminescent source
- FLIT reconstruction of a fluorescent source
- CT volume



NOTE: The 3D Multi-Modality tools (see [page 242](#)) are required to load IVIS Spectrum CT volumetric data or import volumetric data (PET, MRI, or CT data from instruments other than the IVIS Spectrum CT).


A 3D ROI measures the signal intensity within a user-specified bounding box.




Living Image software records information about the ROIs you create during a session and computes statistical data for the ROI measurements. The ROI Measurements table displays the data and provides a convenient way to review or export ROI information ([Figure 6.2](#)).

If a data set includes ROIs on both 2D optical and 3D data, the measurements for the two types of ROIs are displayed in separate tabs of the ROI table ([Figure 6.2](#)):

- ROI Measurements tab shows ROI measurements for 2D optical data.
- 3D ROI Measurements tab shows measurements of structural data (for example, CT) or measurements of the 3D optical reconstruction.

Figure 6.2 3D ROI Measurements Table


6.2 Drawing a 3D ROI

1. Load DLIT or FLIT results.
2. Click the 3D ROI button  in the ROI tools (Figure 6.3).
A red bounding box appears in the 3D View.




If you do not see the red bounding box in the 3D View, do either of the following:

- Select the "Maximum Intensity Projection (MIP)" option in the 3D Multi-Modality tools
- Reduce the volume opacity by adjusting the position of the Air/Noise Boundary in the 3D Multi-Modality tools.

Figure 6.3 3D ROI

3. Adjust the position of the 3D ROI using the transform tools:



NOTE: It may be helpful to view the surface and/or reconstruction results from different perspectives to check the 3D ROI position and size. To turn and rotate the surface, press and hold the left mouse key, then drag the mouse when the hand  appears.





- a. Click the 3D ROI Transform button  and select the ROI from the drop-down list.
The first 3D ROI created during a session is named “ROI 1” by default. A tooltip shows the ROI name when you put the mouse pointer over an ROI.
- b. Click the 3D ROI to begin using the transform tools.
[Figure 6.4](#) explains the tool functions. The ROI position is updated in the slice windowpanes (coronal, sagittal, and transaxial views) after each adjustment.
- c. Press the Tab key to switch between the transformations tools.
- d. Turn off the transform tool when you finish positioning the ROI (click the 3D ROI Transform button ).

Figure 6.4 3D ROI Transformation Tools




Click and drag the 3D ROI when the yellow “+” appears.



Click and drag a handle to scale (increase or decrease) the ROI size.

Red ■ – Scales on the z-axis.
Blue ■ – Scales on the x-axis.
Green ■ – Scales on the y-axis.



To rotate the 3D ROI on the x,y, or z-axis, click the blue, green, or red circle and drag the mouse arrow in the direction of interest.



NOTE: The 3D ROI location (x, y, or z-coordinates) and dimensions (width, height, or depth) can be viewed and modified in the 3D ROI Properties dialog box. See [page 119](#) for details.



4. Click the 3D ROI Measurement button  in the tool palette to view the intensity measurements ([Figure 6.5](#)).

Figure 6.5 3D ROI Measurements**Table 6.1** 3D ROI Measurements Table

Item	Description
Data Types	<p>3D Volumetric Data – Select this data type to measure the grayscale values of 3D volumetric data such as CT or MRI.</p> <p>Source Voxels – Choose this option to measure the source intensity of the voxels of a 3D optical image.</p>
Measurement Types	<p>3D Volumetric Data:</p> <p>Counts – A measurement of a voxel value. The scale is image specific and may not be consistent between images.</p> <p>Absorption – A measurement of the amount of X-rays absorbed by the voxels.</p> <p>Hounsfield – A measurement of voxel grayscale value in Hounsfield units.</p> <p>Note: Absorption and Hounsfield units are only available for IVIS® Spectrum CT data.</p> <p>Source Voxels:</p> <p>photons/sec – The total flux of a luminescent source.</p> <p>cells – The number of cells for calibrated sources integrated over the 3D ROI.</p> <p>pmol M⁻¹ cm⁻¹ – Fluorescence yield for uncalibrated sources integrated over the 3D ROI.</p> <p>pmol – The number of picomoles for calibrated sources integrated over the 3D ROI.</p>
Sequence Number	The identifier of the active image data.
ROI	Name of the 3D ROI.
Voxels	The number of voxels within the 3D ROI.
3D Volumetric Data: Counts measurements (16-bit scale with values that change from image to image)	<p>Total Counts – the sum of all counts for all voxels inside the 3D ROI.</p> <p>Average Counts – Total Counts/Number of voxels in the 3D ROI</p> <p>Min Counts – The smallest number of counts in a voxel within the 3D ROI.</p> <p>Max Counts –The largest number of counts in a voxel within the 3D ROI.</p>

Table 6.1 3D ROI Measurements Table (continued)

Item	Description
3D Volumetric Data: Absorption Measurements (Fixed 32-bit scale with values that are consistent between images.) Note: These measurements are only available for IVIS Spectrum CT data.	Total Value – The sum of the absorption measurements of all voxels in the 3D ROI. Average Value – Total Value/Number of voxels in the 3D ROI. Stdev Value – Standard deviation of the absorption values for all voxels inside the ROI. Min Value – The smallest absorption value for any single voxel in the 3D ROI. Max Value – The largest absorption value for any single voxel in the 3D ROI.
3D Volumetric Data: Hounsfield measurements (Calibrated CT scale. Fixed from image to image.) Note: These measurements are only available for IVIS Spectrum CT data.	Total Hounsfield – The sum of the Hounsfield unit values for all of the voxels in the 3D ROI. Average Hounsfield – Total Hounsfield unit value/Number of voxels in the 3D ROI. Stdev Hounsfield – Standard deviation of the Hounsfield unit values for all voxels inside the ROI. Min Hounsfield – The minimum Hounsfield unit value for any single voxel in the 3D ROI. Max Hounsfield – The maximum Hounsfield unit value for any single voxel in the 3D ROI.
Source Voxels: photons/sec measurements	Total Flux [ph/s] – The flux in each voxel summed or integrated over the 3D ROI. Average Flux [ph/sec] – Total flux/Number of voxels in the 3D ROI. Stdev Flux – Standard deviation of the flux of the voxels inside the ROI. Min Flux – The smallest flux value of a voxel. Max Flux – The largest flux value of a voxel.
Source Voxels: cells Note: This measurement type requires a quantification database. See Chapter 12 on page 235 for more details.	Total Cells – The number of cells in the 3D ROI. Average Cells – Total number of cells/Number of voxels in the 3D ROI. Stdev Cells – Standard deviation of the number of cells in the 3D ROI. Min Cell – The smallest number of cells in a voxel included in the 3D ROI. Max Cell – The largest number of cells in a voxel included in the 3D ROI.
Source Voxels: pmol M ⁻¹ cm ⁻¹ measurements	Total pmol M ⁻¹ cm ⁻¹ – The fluorescence yield summed or integrated over the 3D ROI. Average pmol M ⁻¹ cm ⁻¹ – Total fluorescence yield/Number of voxels in the 3D ROI. Stdev pmol M ⁻¹ cm ⁻¹ – Standard deviation of the fluorescence yield of the voxels in the 3D ROI. Min pmol M ⁻¹ cm ⁻¹ – The smallest fluorescence yield in the 3D ROI. Max pmol M ⁻¹ cm ⁻¹ – The largest fluorescence yield in the 3D ROI.
Source Voxels: pmol measurements Note: This measurement type requires a quantification database. See Chapter 12 on page 235 for more details.	Total pmol – Total picomoles of fluorescent probe within the 3D ROI. Average pmol – Total picomoles/Number of voxels. Stdev pmol – Standard deviation of the picomole values in the 3D ROI. Min pmol – Smallest picomole value in the 3D ROI. Max pmol – Largest picomole value in the 3D ROI.
Refresh	Updates the ROI Measurements table (for example, after you draw new ROIs, move an ROI, and close or open image data).
Copy	Copies the selected row(s) in the table to the system clipboard.
Select All	Copies all rows in the table to the system clipboard.


Table 6.1 3D ROI Measurements Table (continued)

Item	Description
Configure	Displays the Configure Measurements box that enables you to specify and organize the data categories (column headers) for the table. See page 120 for more details.
Export	Opens a dialog box that enables you to export the ROI measurements (.txt or .csv).
Close	Closes the ROI Measurements table.


ROI Properties

You can view information about the location and dimensions of a 3D ROI and edit these properties.

1.

Click the 3D ROI Transform button  and select an ROI from the drop-down list.
2.

Double-click the 3D ROI.
The 3D ROI Properties dialog box appears.




3.

Enter new values or use the arrows in the dialog box to modify the location or dimensions of the 3D ROI in the x, y, or z-planes.
4.


Enter new values or use the arrows in the dialog box to rotate the 3D ROI in the x, y, or z-planes.
5.

To change the color of the 3D ROI:

a.

Click the **Browse** button .

The Select Color box appears ([Figure 6.7](#)).


Figure 6.7 Select 3D ROI Color

- b. To select a basic color for the ROI line, click a basic color swatch, and click **OK**.
- c. To define a custom color, drag the crosshairs in the custom color field, adjust the brightness slider, and click **Add to Custom Colors**.
- d. To select a custom color for the ROI line, click a custom color swatch, and click **OK**.

6.3 Managing the 3D ROI Measurements Table

Configuring the 3D ROI Measurements Table

You can customize the data and information (column headers) in the 3D ROI Measurements table. Several preset categories are available in the Measurement Types drop-down list.

Figure 6.8 3D ROI Measurements Table

1. Drag a column header (left or right) in the table to reorder the columns.
2. Click a column header to sort the table in ascending or descending alphanumeric order.
3. Make a selection from the Data Types and Measurement Types drop-down lists to change the data and measurements displayed in the table.

Creating a Custom 3D ROI Table Configuration


A table configuration specifies the column headers in the 3D ROI table. Several preset configurations are available (selected from the Measurements Types drop-down list in the ROI table (Figure 6.8). You can also create a custom table configuration.




NOTE: Preset table configurations cannot be edited. You can modify a preset configuration and save it to a new name.

1. In the ROI Measurements table, click **Configure**.
The Configure Measurements box appears.

Figure 6.9 Configure Measurements Dialog Box



2. Select a configuration from the User Lists drop-down list and click **Customize**.
3. To add column header to the ROI table, make a selection from the “Available Item” list and click **Add**.
4. To remove column header from the ROI table, select the item that you want to remove in the Selected Items list, and click **Remove**.
5. To reorder an item in the Selected Items list, select the item and click **Move Up** or **Move Down**.
The columns in the ROI Measurements table are updated.
6. Enter a name for the custom configuration in the Name box and click **Save**.
7. Select the custom configuration from the Measurements Unit drop-down list (Figure 6.10).

Figure 6.10 Select a Custom Configuration for 3D ROI Measurements Table**To delete a custom table configuration:**

Select the configuration from the User Lists drop-down list and click **Delete**.



NOTE: Preset table configurations cannot be deleted.


Copying or Exporting the ROI Measurements Table

To export the table:

1. In the ROI Measurements table, click **Export**.
2. In the dialog box that appears:
 - a. Select a folder and file type (.txt or .csv).
 - b. Enter a name file and click **Save**.

To copy the table to the system clipboard:

- Copy selected rows - Select the rows of interest and click **Copy**. Alternatively, select the rows, then right-click the table and choose **Copy** on the shortcut menu.
- Copy all rows - Click **Select All** and click **Copy**. Alternatively, press **Ctrl+A**, then right-click the table and choose **Copy** on the shortcut menu.

Figure 6.11 3D ROI Table Shortcut Menu

7 Image Math

About Image Math

Creating a New Image Using Image Math

Subtracting Tissue Autofluorescence on page 125

7.1 About Image Math

The Image Math tool is used to mathematically combine two images to create a new image. Image math is primarily for subtracting tissue autofluorescence background from signal.

To perform image math, open an image sequence or a group of images. For more details on creating a sequence from individual images, see [page 82](#).




TIP: See the tech note *Image Math* for a quick guide (select **Help** → **Tech Notes** on the Help menu).


7.2 Creating a New Image Using Image Math


1. Load an image sequence.
2. Select **Tools** → **Image Math for <name>_SEQ** on the menu bar.


Figure 7.1 Opening the Image Math Window



3. In the Image Math window that appears, select an image from box A and from box B. The Image Math window shows a thumbnail of image A, image B, and the new image.



 **NOTE:** For more details on items in the Image Math window, see [Table 7.1, page 124](#)

4. Select a mathematical function from the Result drop-down list.
5. To include a scaling factor (k) in the function, enter a value for k.
6. To view the new image, click **Display Result for Measuring**.
7. To save the new image:
 - a. Click the **Save** button . Alternatively, select **File** → **Save** on the menu bar.
 - b. In the dialog box that appears, select a directory, and click **Save**.


A folder of data is saved to the selected location (AnalyzedClickInfo.txt, ClickInfo.txt, luminescent and photographic TIF images).
8. To export the image to a graphic file:
 - a. Click the **Export** button  ([Figure 7.2](#)).
 - b. Select a directory in the dialog box that appears, enter a file name, and select the file type from the “Save as type” drop-down list.
 - c. Click **Save**.

Table 7.1 Image Math Window

Item	Description
Color Ranges for A and B	<p>Full - Choose this option to set the Max and Min values to the maximum and minimum data values in the image.</p> <p>Auto - When this option is chosen, the software sets the Min and Max values to optimize image display and suppress background noise. The Min and Max settings can be manually adjusted to further optimize the image display for your needs.</p> <p>Note: The color scale does not affect the image math result.</p>

Table 7.1 Image Math Window (continued)

Item	Description
Color Ranges for Result Image	Full - See above. Auto - See above. Min = 0 - Choose this option to set the minimum data value to zero.
Results	Drop-down list of mathematical functions that can be used to generate the new image, including: A - B*k A + B*k A * B*k A/B if Counts(B)>k (Useful for fluorescence tomography.) (A/B)*k
k, Image Math window	A user-specified scaling factor applied in the results function.
Compute 'k' from ROI	This option is useful for subtracting fluorescence background. Draw one ROI in an image on an area considered background. In the "Compute 'k' from ROI" drop-down list, select the this ROI.
with Photo from	Choose this option to display the new image in overlay mode using the selected photographic image. (This option is only available if one of the selected images is an overlay.
Display Result for Measuring	Opens the image generated by image math in an image window.

7.3 Subtracting Tissue Autofluorescence

To remove tissue autofluorescence from image data, you can use a subtraction method that uses a second excitation filter which is blue-shifted (a background filter) from the primary excitation filter.

The objective of using a background filter is to excite the tissue autofluorescence without exciting the fluorophore. To reduce autofluorescence signal in the primary image data, use the image math tool to subtract the background filter image from the primary excitation filter image.

The software computes the signal corrected for background: $(A - B) \times k$, where:

A = primary image (acquired using the excitation filter)

B = background image (acquired using the background filter)


k = (primary signal/background signal)

The background signal is obtained from a measurement ROI that is located in an area where no fluorophore signal is present. The scale factor k accounts for different levels of tissue autofluorescence due to different excitation wavelengths and filter transmission characteristics.

After you acquire an image sequence that includes a primary and background image, use the image math tool to subtract tissue autofluorescence. (For more details on acquiring an image sequence, see [Chapter 3 on page 39](#).)

To subtract tissue autofluorescence:


1. Load the image sequence that includes the primary and background fluorescent images.

Figure 7.3 Image Sequence

2. Open either the primary or background image and:
 - a. Optimize the image display using the color scale Min and Max sliders in the Image Adjust tools.
 - b. Draw a measurement ROI on an area of the animal that represents background signal (area where no fluorophore signal is present).



NOTE: You only need to draw the ROI on one of the images. The software copies the ROI to the other image.


Figure 7.4 Draw Measurement ROI on an Area Representative of Background Signal

3. Select **Tools** → **Image Math** for <name>_SEQ on the menu bar.

4. In the Image Math window that appears, select the primary image in box A. Select the background image in box B.

For more details on items in the Image Math window, see [Table 7.1, page 124](#).

5. Select the math function 'A-B*k' in the Result drop-down list.



6. Click **Compute 'k' from ROI** and select the ROI (created in step 2) from the drop-down list.
The background-corrected signal is displayed.
7. To view the mathematical result (overlay mode) in a separate image window, click **Display Result For Measuring**.
If necessary, use the Color Scale Min and Max sliders in the Image Adjust tools to adjust the image display.
8. To save the new image:
 - a. Click the **Save** button . Alternatively, select **File** → **Save** on the menu bar.
 - b. Select a directory in the dialog box that appears and click **Save**.
A folder of data is saved to the selected location (AnalyzedClickInfo.txt, ClickInfo.txt, luminescent and photographic TIF images).
9. To export the new image to a graphic file:
 - a. Click the **Export** button .
 - b. Select a directory in the dialog box that appears, enter a file name, and select the file type from the "Save as type" drop-down list.
 - c. Click **Save**.

8 Spectral Unmixing

About Spectral Unmixing

Acquire an Image Sequence for Spectral Unmixing on page 128

Spectral Unmixing Methods on page 137

Correcting Spectra on page 149

Spectral Unmixing Results on page 150

8.1 About Spectral Unmixing

Living Image software applies spectral unmixing to distinguish the spectral signatures of different fluorescent or luminescent reporters and calculate the respective contribution of each on every pixel of an image. Use spectral unmixing to:

- Extract the signal of one or more fluorophores from the tissue autofluorescence. Images are acquired using epi-illumination (excitation light above the stage) or transillumination (excitation light below the stage).
- Analyze luminescent or fluorescent images when more than one reporter is used in the same animal model.

8.2 Acquire an Image Sequence for Spectral Unmixing

Set up an image sequence for spectral unmixing using the Imaging Wizard.



TIP: See the *Imaging Wizard* tech note for a quick guide on sequence acquisition (select **Help** → **Tech Notes** on the menu bar).

Choose an imaging mode in the wizard based on the type of probes.

Probe Type	Follow the Instructions for:
Luminescent	<i>Bioluminescence Imaging on page 129</i>
Fluorescent	<i>Fluorescence Imaging on page 131</i>
Radio-isotope	<i>Cherenkov Imaging on page 134</i>

If you are not using the Imaging Wizard to set up the image sequence, it is recommended that the sequence include images acquired using several filters that sample the emission or excitation spectra at multiple points across the entire range. Make sure that the band gap between the excitation and emission filters is sufficiently large so that the excitation light does not leak through the emission filter where it can be detected by the CCD.

If a data set includes multiple excitation and emission filter scans, the software automatically unmixes signal according to the filter type with the most entries. For example, a data set acquired using three excitation filters and four emission filters will be unmixed by emission wavelength.


Bioluminescence Imaging



NOTE: The IVIS Spectrum should be initialized and the temperature locked before setting the imaging parameters. See [page 7](#) for more details.


1. Start the Imaging Wizard. See [Starting the Imaging Wizard on page 40](#) for instructions.
2. Double-click the Bioluminescence option. Double-click the Spectral Unmixing option in the next screen ([Figure 8.1](#)).

Figure 8.1 Imaging Wizard – Choose Bioluminescence and Spectral Unmixing Options




3. Select a probe from the Name drop-down list ([Figure 8.2](#)).
4. Click **Add** and select another probe from the Name drop-down list. Repeat until all of the probes are added.

Figure 8.2 Select Probes




5. Optional: If you generated a spectrum library (a set of reference spectra), you can select it in the Imaging Wizard. Click **Filter Config**, then click **Select by Spectrum Library** in the dialog box that appears (Figure 8.3).

See [Guided Method on page 137](#) for instructions on creating a spectrum library.




6. Click **Next** and in the screen that appears (Figure 8.4):
 - a. Select the type of subject.
 - b. Select a field of view.
 - c. Set the focus options.



7. To acquire a time series of images:

- a. Choose the Time Series Study option (Figure 8.4).
 - b. Enter the number of segments and the delay between segments.
 8. Click **Next**.
- The specified sequence appears in the sequence table.

Figure 8.5 Sequence Setup Complete




9. Acquire the sequence following the instructions on page 41.
- The image window appears when acquisition is completed (Figure 8.17). See Table 3.2 on page 28 for more details on the Image window.


Fluorescence Imaging

1. Start the Imaging Wizard. See *Starting the Imaging Wizard* on page 40 for instructions.
2. Double-click the Fluorescence option.
3. Select Spectral Unmixing in the next screen (Figure 8.6).
4. Select the type of illumination and click **Next**.
 - Epi-Illumination – Light source above the stage.
 - Trans-Illumination – Light source below the stage.

Figure 8.6 Imaging Wizard – Choose Fluorescence and Spectral Unmixing Options




5. Select a probe from the Name drop-down list in the next screen (Figure 8.7).
6. Click **Add** and select another probe from the Name drop-down list. Repeat until all of the probes are added.




7. Optional: If you generated a spectrum library (a set of reference spectra), you can select it in the Imaging Wizard. Click **Filter Config**, then click **Select by Spectrum Library** in the dialog box that appears (Figure 8.8).


See [Guided Method on page 137](#) for instructions on creating a spectrum library.



8. Click **Next** and in the screen that appears (Figure 8.9):
 - a. Select the type of subject.
 - b. Select a field of view.
 - c. Set the focus options.




9. To acquire a time series of images:
 - a. Choose the Time Series Study option (Figure 8.9).
 - b. Enter the number of segments and the delay between segments.
10. If using transillumination, select the transillumination locations.
 - a. Click **Transillumination Setup**.
 - b. Choose the transillumination locations in the Transillumination Setup box that appears (Figure 8.10).
See Table 3.3 on page 37 for more details on Transillumination Setup.



11. Click Next.

The specified sequence appears in the sequence table (Figure 8.11).

Figure 8.11 Sequence Setup Complete**12. Acquire the sequence following the instructions on page 41.**


The image window appears when acquisition is completed (Figure 8.17). See Table 3.2 on page 28 for more details on the Image window.

Cherenkov Imaging




NOTE: The IVIS Spectrum should be initialized and the temperature locked before setting the imaging parameters. See page 7 for more details.

1. Start the Imaging Wizard. See *Starting the Imaging Wizard* on page 40 for instructions.
2. Double-click the Cherenkov option. Double-click the Spectral Unmixing option in the next screen (Figure 8.12).


Figure 8.12 Imaging Wizard – Choose Cherenkov and Spectral Unmixing Options

3. Select a probe from the Name drop-down list (Figure 8.13).
4. Click **Add** and select another probe from the Name drop-down list. Repeat until all of the probes are added.



5. Optional: If you generated a spectrum library (a set of reference spectra), you can select it in the Imaging Wizard. Click **Filter Config**, then click **Select by Spectrum Library** in the dialog box that appears (Figure 8.14).


See [Guided Method on page 137](#) for instructions on creating a spectrum library.



6. Click **Next** and in the screen that appears (Figure 8.15):
 - c. Select the type of subject.
 - d. Select a field of view.


- e. Set the focus options.

Figure 8.15 Select Subject and Set Acquisition Parameters



7. To acquire a time series of images:
 - a. Choose the Time Series Study option (Figure 8.15).
 - b. Enter the number of segments and the delay between segments.
8. Click **Next**.
The specified sequence appears in the sequence table (Figure 8.16).

Figure 8.16 Sequence Setup Complete



9. Acquire the sequence following the instructions on page 41.
The image window appears when acquisition is completed (Figure 8.17). See Table 3.2 on page 28 for more details on the Image window.

8.3 Spectral Unmixing Methods

Living Image software provides four spectral unmixing methods ([Table 8.1](#)).

Table 8.1 Spectral Unmixing Methods


Method	Description	See Page
Guided	<p>Use this method when:</p> <ul style="list-style-type: none"> ■ Probe locations are known. ■ Probe signals are mixed with background signal, but not other probe signals. <p>Note: This method is not recommended if probe signals are overlapping.</p> <p>Use this method to generate a spectrum library (a set of reference spectra) for probes with known spectra and known locations.</p>	137
Library	<p>This method requires a user-generated spectrum library. The library method identifies pixels in the data with spectral characteristics that match the spectrum library.</p> <p>Note: The data being analyzed must be acquired using the same, or a subset of, the excitation/emission filter pairs of the spectrum library. The probe depth in the data being analyzed and the spectrum library data set should be similar for optimum analysis results. For example, do not use a spectrum library generated from <i>in vivo</i> data to analyze <i>in vitro</i> data and vice versa.</p>	137
Automatic	Use this method when the probe locations are unknown.	142
Manual	If necessary, perform a manual analysis after an automatic analysis to identify additional probe locations.	145

Guided Method


Use the guided method:

- When the probe locations are known and probe signals do not overlap.
 - To generate a spectrum library for probes with known spectra and known locations
1. Load the image sequence.

In [Figure 8.17](#), the fluorophores are Alexa Fluor 680 and Alexa Fluor 750. Images were acquired using a 605 nm excitation filter and emission filters from 660 to 800 nm in 20 nm increments.

Figure 8.17 Sequence for Spectral Unmixing

2. Click the Analyze tab of the Spectral Unmixing and DyCE tools.
By default, all wavelengths are included in the analysis. Remove the check mark next to wavelengths that you want to exclude from the analysis.
3. Select "Guided" from the Methods drop-down list and click **Start Unmixing**.
The Unmixing window appears ([Figure 8.18](#)).

Figure 8.18 Unmixing Window

See [Table 8.3](#) on page 152 for more details on the toolbar buttons.


List of the spectral components to unmix. If the Imaging wizard was used to set up the sequence, the list includes the probes selected in the wizard. The default list also includes a tissue autofluorescence background component.

The image cube represents a “stack” of the sequence images (sorted according to the spectral axis). When the Overview option is selected, the image cube shows a pseudo color image that is a composite of the stack images which have been colorized to encode spectral information.

The entire image cube is calibrated and visualized on the same scale. To view a particular image, remove the check mark next to the Overview option and move the slider or enter an image number.



NOTE: In the Guided method, the Tissue AF component is preset as background. After you define the Tissue AF component (mark a region of tissue autofluorescence only on the image cube), the spectra of the other components that you mark on the image cube will be background-subtracted, not raw spectra from the data.


4. Move the mouse pointer over the image cube to see the spectrum at a particular location.
The raw spectrum at the pointer location is updated as you move the pointer.
5. To specify a probe location for unmixing:
 - a. Click the  button for a spectrum.
 - b. Using the mouse, draw a mark on an area of the image cube which represents the probe signal.
The software plots a background-subtracted spectrum of the signal ([Figure 8.19](#)).



NOTE: For “Tissue AF”, draw a mark on an area of the image cube where no probe signal is present.

- c. If necessary, right-click the image cube to erase the mark.
6. Repeat step [step 5](#) to specify other probe locations.

Figure 8.19 Mark Probe Locations for Unmixing on the Image Cube



7. Click **Next** after you finish marking the probe locations.

The Unmixing window shows the analysis results which include unmixed spectra corrected for tissue autofluorescence, unmixed images, and a composite of the unmixed images ([Figure 8.20](#)).

See [Spectral Unmixing Results](#) on page 150 for information about the results.



8. To save the results as a spectrum library:
 - a. Click the  button in the Spectrum List toolbar (Figure 8.20).
 - b. Enter a file name in the dialog box that appears and click **Save**.

Figure 8.20 Spectral Unmixing Results



Library Method

The library method uses a user-generated spectrum library to analyze a data set. If you plan to analyze data by this method, the data must be acquired using the same, or a subset of, the excitation/emission filter pairs of the spectrum library.


The probe depth in the data set being analyzed and the spectrum library data set should be similar for optimum analysis results. For example, do not use a spectrum library generated from in vivo data to analyze in vitro data.




NOTE: Use the "guided" method to generate a spectrum library of known probes with known locations (see [page 137](#) for more details on the guided method).

1. Load the image sequence.

In [Figure 8.21](#), the fluorophores are Alexa Fluor 680 and Alexa Fluor 750. Images were acquired using a 605 nm excitation filter and emission filters from 660 to 800 nm in 20 nm increments.



- Click the Analyze tab of the Spectral Unmixing and DyCE tools.
By default, all wavelengths are included in the analysis. Remove the check mark next to wavelengths that you want to exclude from the analysis.
- Select "Library" from the Methods drop-down list and click **Start Unmixing**.
- Select a reference spectral library in the dialog box that appears and click **Apply** (Figure 8.22).
The software identifies pixels with spectral characteristics that match the spectrum library. The Unmixing window shows the analysis results which include unmixed spectra, unmixed images, and a composite of the unmixed images (Figure 8.20 on page 140).
See *Spectral Unmixing Results* on page 150 for information about the results.




Automatic Method

Use the automatic method to analyze data when the probe locations are unknown.

1. Load the image sequence.

In [Figure 8.23](#), the fluorophores are Alexa Fluor 680 and Alexa Fluor 750. Images were acquired using a 605 nm excitation filter and emission filters from 660 to 800 nm in 20 nm increments.




2. Click the Analyze tab of the Spectral Unmixing and DyCE tools.

By default, all wavelengths are included in the analysis. Remove the check mark next to wavelengths that you want to exclude from the analysis.

3. Select "Automatic" from the Methods drop-down list and click **Start Unmixing**.


The Auto Unmix window appears. The purple data mask shows the data that will be included in the analysis (the entire subject is included by default).




4. If you do not want to analyze the entire subject, draw a mask on a particular area (Figure 8.25).



Table 8.2 Data Mask Options

Option	Description
Photograph	If this option is chosen, the software automatically draws the data mask so that it includes the entire photograph.
Threshold	If necessary use the threshold slider or  arrows to adjust the mask so that it matches the underlying subject photograph as closely as possible without including any area outside the subject image.
Draw Mask	Choose this option to manually draw a data mask on an area of the photograph.
Rectangle	Specifies a rectangular shape for the manual data mask.
Ellipse	Specifies an elliptical shape for the manual data mask.


5. Choose an imaging subject and background signal(s).



6. Click the **PCA** button.

The Principle Component Analysis window shows the amount of signal explained by the suggested components (Figure 8.27). The three components in this example (tissue autofluorescence, probe AF680, and probe AF750) explain more than 99.5% of the signal. The small residual is due to noise.

If the explained variance is low, add more components (probes) to unmix using the  button.



7. Click Finish.


The Unmixing window shows the analysis results which include unmixed spectra, unmixed images, and a composite of the unmixed images (Figure 8.20 on page 140).

See *Spectral Unmixing Results* on page 150 for information about the results.

Manual Method

Sometimes you may want to manually analyze results, for example, if the explained variance of the principle component analysis of an automatic analysis seems low. The example in this section shows how to manually analyze results from a previous analysis.

1. Open the image sequence.
2. Select the results and click **Load**.




3. Click the Analyze tab of the Spectral Unmixing and DyCE tools.

All wavelengths are selected by default. Remove the check mark next to wavelengths that you want to exclude from the analysis.


In [Figure 8.29](#), the fluorophores are Alexa Fluor 680 and Alexa Fluor 750. Images were acquired using a 605 nm excitation filter and emission filters from 660 to 800 nm in 20 nm increments.

Figure 8.29 Unmixing Window – Results Loaded




4. Select "Manual" from the Methods drop-down list and click **Start Unmixing**.
The Unmixing window appears ([Figure 8.18](#)).

Figure 8.30 Unmixing Window



The image cube represents a “stack” of the sequence images (sorted according to the spectral axis). When the Overview option is selected, the image cube shows a pseudo color image that is a composite of the stack images which have been colorized to encode spectral information.

The entire image cube is calibrated and visualized on the same scale. To view a particular image, remove the check mark next to the Overview option and move the slider or enter an image number (Figure 8.31).




5. Move the mouse pointer over the image cube to see the spectrum at a particular location. The spectrum at the pointer location is updated as you move the pointer.
6. To specify a probe location for unmixing:
 - a. Click the button for a spectrum.
 - b. Using the mouse, draw a mark on an area of the image cube which represents the probe location. The software plots a normalized spectrum of the signal (Figure 8.32).



NOTE: Mark a region of tissue autofluorescence only (where no probe signal is present) on the image cube for the Tissue AF component. The spectra of components that you mark on the image cube are raw spectra from the data when using the manual method.

- c. If necessary, right-click the image cube to erase the mark.
7. Repeat step [step 6](#) to specify other probe locations.
8. Manually subtract autofluorescence background. See [Correcting Spectra on page 149](#) for instructions.

Figure 8.32 Mark Probe Locations for Unmixing on the Image Cube


9. Click **Unmix** after you finish marking the probe locations and correct spectra for tissue autofluorescence.

The Unmixing window shows the analysis results which include unmixed spectra, unmixed images, and a composite of the unmixed images (Figure 8.20 on page 140).

See [Spectral Unmixing Results](#) on page 150 for information about the results.

8.4 Correcting Spectra

Spectra can be corrected for overlapping signal by subtracting one spectrum from another.

1. Click the  button in the Unmix window.
2. Choose the spectra to subtract in the dialog box that appears. (Figure 8.33).
3. Click **Apply** to add the computed spectrum to the spectrum plot and list in the Unmixing window. Alternatively, select an existing spectrum from the Name drop-down list and click **Apply** to overwrite the results.





Table 8.1 Computed Spectrum

Item	Description
Normalized	Choose this option to display spectra normalized on a scale from zero to one.
Result: $C = A - x*B$	The subtraction performed by the software where "x" is a factor that ensures the residual signal is positive.
Autoscaling	Choose this option to display computed results on a normalized scale starting a zero.
Fit Offset	If this option is chosen, the software computes and removes an intensity baseline from the spectra.


Table 8.1 Computed Spectrum (continued)

Item	Description
Error Tolerance	The software computes a default error tolerance (the factor "x" for A - x*B) such that signal B is maximally removed from signal A with no negative result. Moving the slider adjusts the error tolerance and automatically updates the computed spectrum.
	Choose "New" to save computed spectrum with the specified name and color. Click Apply to add the computed spectrum to the spectrum plot and list in the Unmixing window. Choose a spectrum number from the drop-down list to overwrite that spectrum with the computed spectrum when you click Apply .

8.5 Spectral Unmixing Results


The results include a signal distribution map of each unmixed result and a composite image of all signals, each displayed in a different color.

Figure 8.34 Spectral Unmixing Results



Remove the check mark next to “Individual” to view a signal color scale (Figure 8.35).

Figure 8.35 Spectral Unmixing Results With Color Scale



Spectra Plot

Spectra plots show the unmixed spectra.












Table 8.3 Spectra Window

Item	Description
Normalized	Choose this option to display signals normalized on a scale from zero to one.
Legend	Choose this option to display a key for the spectra plot.
	Opens a dialog box that enables you to export the spectra plot data to a .csv file.
	Opens a dialog box that enables you to select and load a spectrum library.
	Opens a dialog box that enables you to save spectral unmixing results as a reference spectrum library for use with the "library" method of spectral unmixing. See page 140 for more details on the library method.
	Enables you to view and save the unmixed images as a sequence data set which can be analyzed using the tool palette.
	Opens a dialog box that enables you to correct a spectrum for overlapping signal by subtracting one spectrum from another (see page 149).
	Adds a component to the spectrum list.
	Deletes the last spectrum in the spectrum list.

Adding Spectra to the Plot

To Add:	Do This:
A spectrum library	Click the  button and select a spectrum library in the dialog box that appears. Note: A spectrum library is a user-created set of reference spectra generated by analyzing probes with known spectra and known locations.
A spectrum from a user-defined region	Add a new spectrum to the list in the Unmix window and identify the region by drawing a mark on the image cube. See page 146 for more details.

Composite Image

The composite image includes all of the signals, each displayed in a different color. Double-click the composite image to view it in a separate window (Figure 8.37).

Figure 8.37 Composite Window









Table 8.2 Composite Window

Item	Description
Units	The type of data displayed in the composite image.
Image list	A list of the images that comprise the composite (background component(s), probe(s), and a photograph).
Min/Max	Sets the minimum and maximum count to display in the image.
Brightness	Adjusts the brightness of the component signals.
Logarithmic Scale	Choose this option to display signals using a logarithmic scale. This may be useful when probe signal strengths differ significantly, for example, a bright source and a dim source.
Color	Shows the color of the figure legend for the image selected in the image list. Click the color swatch to open a color palette that enables you to select a new color for the figure legend.
Label	The name of the image selected in the image list. To edit the name, double-click the name in this box. Right-click the label name to show a short-cut menu of edit commands (for example, Cut, Copy, Paste).
	Sends the composite image to the "top" of the image cube. This helps improve the pseudo color visualization of the image cube.
	Copies the composite image to the system clipboard.
	Click to export the composite image to a graphic file (for example, .jpg).
	Opens the Print dialog box.


Analyzing Images

Do either of the following:

- Click the  button toolbar button to view all images as a sequence.
- Double-click a particular unmixed image.

The image(s) appear in a separate window and the tool palette is available for image analysis. When closing the window, the software prompts you to save the sequence or image.

Figure 8.38 View Unmixed Images as a Sequence



Managing Spectral Unmixing Results

Figure 8.39 Spectral Unmixing Results




Table 8.4 Spectral Unmixing Results Tab

Item	Description
Name	The name for the active spectral unmixing results. Select results from this drop-down list.
Delete	Deletes the selected results.
Load	Opens the selected results in the Unmixing window.
Save	Saves the active results using the selected name. The results are saved to the sequence click number folder and are available in the Name drop-down list.
Overwrite	If you reanalyze results, saves the new results and overwrites the previous results.

9 DyCE Imaging and Analysis

About DyCE (Dynamic Contrast Enhancement)

Acquire an Image Sequence for DyCE Analysis on page 157

DyCE Analysis on page 164

DyCE Results on page 171

9.1 About DyCE (Dynamic Contrast Enhancement)



NOTE: The DyCE acquisition and analysis features of Living Image software require a separate license.


DyCE imaging and analysis is intended for biodistribution studies. DyCE imaging captures a time series of optical images immediately following a bolus injection of a probe or dye. Living Image software temporally unmixes the data on a pixel-by-pixel basis for each image of the time series and determines real-time spatio-temporal distribution of the probe or dye signal.

Living Image software presents the spatio-temporal information as:

- Temporal spectra – Line plots of signal intensity as a function of time. Each line plot represents the signal time course within a particular anatomical region.
- An unmixed image – An image representing the peak signal time point for a particular temporal spectrum.
- A composite image – An overlay of the unmixed images.

Figure 9.1 Example DyCE Results

Images were obtained using the Mouse Side Imaging Kit.



9.2 Acquire an Image Sequence for DyCE Analysis

A DyCE sequence is set up using the Imaging Wizard and includes a user-specified time delay between exposures. An acquisition can include up to three different time intervals where each interval is defined by duration and the delay between exposures.



NOTE: For optimum DyCE analysis results, acquire images using the Side Imager accessory (PN CLS135111).

Choose an imaging mode in the wizard based on your probe type.

Probe Type	Follow the Instructions for:
Luminescent	Bioluminescence Imaging (below)
Fluorescent or near infrared	Fluorescence Imaging on page 159
Radiotracer	Cherenkov Imaging on page 162


Bioluminescence Imaging



NOTE: The IVIS Spectrum should be initialized and the temperature locked before setting the imaging parameters. See [page 7](#) for more details.

1. Start the Imaging Wizard. See [Starting the Imaging Wizard on page 40](#) for instructions.
2. Double-click the Bioluminescence option. Double-click the DyCE option in the next screen ([Figure 9.2](#)).

Figure 9.2 Imaging Wizard – Choose Bioluminescence and DyCE Options



If this screen does not appear when the wizard starts, click [Restart Wizard](#) (at the lower left wizard screen).

3. Select the type of imaging subject in the next screen ([Figure 9.3](#)).

Figure 9.3 Imaging Wizard – Bioluminescence DyCE

Click the button to view definitions of the time series parameters.

Interval – Defined by a Duration and Delay Time.

Duration – The amount of time that the interval lasts.

Delay – The amount of time between imaging time points in the interval.

Note: Each image exposure time must be less than the Delay Time.

4. Choose “Manual Settings” and set appropriate exposure parameter values for your probe.



NOTE: Typical exposures are in the range of minutes because Cherenkov light emission is very low. Beta decays with higher energies can afford shorter exposure times than lower energy beta decays.

5. Select a field of view from the drop-down list.
6. Set the focus by doing either of the following:
 - Enter a subject height and choose the “use subject height” focus option.
 OR
 - Choose the “manual focus” option from the Focus drop-down list and set the focus parameters in the Manual Focus Window that appears.



NOTE: If using the Side Imaging accessory for bioluminescence DyCE, set the subject height = 0.0 cm and FStop = 2 (or larger). If using the Side Imaging accessory for fluorescence DyCE, choose the Manual Settings options and set the subject height = 0.0 cm and FStop = 4 (or larger).

7. Specify the time series.



A time series can include up to three intervals. Each interval is defined by duration (minutes) and delay between images (seconds) (Figure 9.3).

Maximum number of images = Duration/Exposure if exposure is greater than delay.

Maximum number of images = Duration/Delay if exposure is less than delay.

A time series can include a maximum of 200 images.

- a. Enter the number of intervals.
- b. Enter the duration and the delay between images for each interval.
The software computes the number of images to acquire during the interval.




NOTE: The software alerts you if the number of images in the time series exceeds 200. If necessary, adjust the duration or delay between images of one or more intervals to reduce the number of images.

c. Click **Next**.

The specified sequence appears in the sequence table (Figure 9.4).

Figure 9.4 Sequence Setup Complete



8. Acquire the sequence following the instructions on page 41.

The image window appears when acquisition is completed (Figure 9.13 on page 164). See Table 3.2 on page 28 for more details on the Image window.

Fluorescence Imaging




NOTE: The IVIS Spectrum should be initialized and the temperature locked before setting the imaging parameters. See page 7 for more details.

1. Start the Imaging Wizard (see *Starting the Imaging Wizard* on page 40).
2. Double-click the Fluorescence option (Figure 9.5).
3. In the next screen, select DyCE and the type of illumination (Figure 9.5):
 - Epi-Illumination – Excitation light source above the stage.
 - Trans-Illumination – Excitation light source below the stage. If this option is selected, NTF Efficiency images will be produced in which the fluorescent emission image is normalized by the transmission image measured with the same emission filter and open excitation filter.



TIP: See these tech notes for helpful information and quick guides (select **Help → Tech Notes** on the menu bar):

- *Transmission Fluorescence*
- *Transmission Fluorescence – Normalized Transmission Fluorescence*


Figure 9.5 Choose the Fluorescence and DyCE Options

If this screen does not appear when the wizard starts, click **Restart Wizard** (at the lower left wizard screen).

4. Click Next.


5. Select a probe from the Name drop-down list in the next screen (Figure 9.6).

If your fluorescent probe is not in the list, select “Input” and enter the fluorescence excitation and emission peak wavelengths. Click **Next**.

Figure 9.6 Select a Probe

6. Select the type of imaging subject in the next screen (Figure 9.7).


7. Choose the Auto Settings option.



8. If using transillumination, select a transillumination location. Only one excitation point location is allowed.

- a. Click **Transillumination Setup**.
- b. Choose the location for transillumination by clicking a square of the grid in the Transillumination Setup box that appears (Figure 9.8).


See Table 3.3 on page 37 for more details on the Transillumination Setup.



9. Perform [step 5 on page 158](#) to [step 7 on page 158](#).

The specified sequence appears in the sequence table ([Figure 9.9](#)).

Figure 9.9 Sequence Setup Complete



10. Acquire the sequence following the instructions on [page 41](#).

The image window appears when acquisition is completed ([Figure 9.13 on page 164](#)). See [Table 3.2 on page 28](#) for more details on the Image window.


Cherenkov Imaging



NOTE: The IVIS Spectrum should be initialized and the temperature locked before setting the imaging parameters. See [page 7](#) for more details.

1. Start the Imaging Wizard. See [Starting the Imaging Wizard on page 40](#).
2. Double-click the Cherenkov option. Double-click the DyCE option in the next screen ([Figure 9.10](#)).

Figure 9.10 Choose Cherenkov and DyCE Options



3. Select the subject type and radio-isotope from the drop-down lists ([Figure 9.11](#)). If your radio-isotope is not available in the list, choose "Undefined".

- Choose the Manual Settings option and set exposure parameter values that are appropriate for your radiotracer probe.



NOTE: Typical exposures are in the range of minutes because Cherenkov light emission is very low. Beta decays with higher energies allow shorter exposure times than lower energy beta decays.

Figure 9.11 Imaging Wizard – Cherenkov DyCE

Imaging Subject: Mouse Isotope: F-18

Exposure Parameters

☐ Auto Settings ☒ Manual Settings

Exposure: 5.00 min Binning: Medium F/Stop: 1

☒ Luminescent ☒ Photograph

Field of View: C - 12.8 cm

Focus: Subject Height: 1.50 cm Focus: use subject height

Options

Time Series

Total number of Intervals: 3

Interval	Duration (min)	Delay (sec)	Max Images
1	10	60	2
2	15	180	3
3	30	300	6

Restart Wizard Cancel Back Next

- Perform [step 5 on page 158](#) to [step 7 on page 158](#).
The specified sequence appears in the sequence table ([Figure 9.12](#)).

Figure 9.12 Sequence Setup Complete

Imaging Wizard - Cherenkov - DyCE

Imaging Subject: Mouse Isotope: F-18

Exposure Parameters

☐ Auto Settings ☒ Manual Settings

Exposure: 5.00 min Binning: Medium F/Stop: 1

☒ Luminescent ☒ Photograph

Field of View: C - 12.8 cm

Focus: Subject Height: 1.50 cm Focus: use subject height

Options

Time Series

Total number of Intervals: 3

Interval	Duration (min)	Delay (sec)	Max Images
1	10	60	2
2	15	180	3
3	30	300	6

Restart Wizard Cancel Back Next

IVIS Acquisition Control Panel

Imaging Mode: Exposure Time: 5.00 min Binning: Medium F/Stop: 1 Excitation Filter: Block Emission Filter: Open

☒ Luminescent ☒ Photograph

Structure: Overlay Lights Alignment Grid Match Sequenced

Field of View: C

System Status: Idle

Acquire Sequence Imaging Wizard Image Setup Initialize

Number of Intervals: 3 Delay: 0.0 min Apply to All Remove Update Insert Add

Seq-1	Mode	Exposure	Binning	F/Stop	Excitation	Emission	Structure	FOV	Height	DyCE	Duration	DyCE Delay
1	Lum	300	Medium	1	Block	Open	No	C	1.50	300min	5sec	
2	Lum	300	Medium	1	Block	Open	No	C	1.50	300min	5sec	
3	Lum	300	Medium	1	Block	Open	No	C	1.50	300min	5sec	

- Acquire the sequence following the instructions on [page 41](#).
The image window appears when acquisition is completed ([Figure 9.13 on page 164](#)). See [Table 3.2 on page 28](#) for more details on the Image window.

9.3 DyCE Analysis


Automatic or manual DyCE analysis is available. PerkinElmer recommends performing an automatic analysis first, followed by manual analysis to identify possible additional temporal components.

Automatic DyCE Analysis


1. Load a DyCE sequence. The  icon in the Living Image browser indicates a DyCE sequence.



NOTE: If the data is noisy, as is common for Cherenkov data, smooth all the images in the sequence using the Smoothing tools in the Image Adjust Tool Palette. This can be done in sequence view mode. See [Smoothing and Binning on page 62](#) for details.



2. Click the Analyze tab in the Spectral Unmixing/DyCE tools.
3. Select **Automatic** from the Methods drop-down list and click **Start Unmixing**.
The Auto Unmix Wizard appears and shows the purple data mask that specifies the analysis area (Figure 9.14). The data mask includes the entire subject by default.
4. If necessary, change the threshold level to adjust the purple mask so that it matches the underlying subject photograph as closely as possible without including any area outside the subject image.




5. If you do not want to analyze the entire subject, draw a mask on a particular area (Figure 9.15).
- a. Select **Draw Mask** and choose the **Rectangle** or **Ellipse** option.
 - b. Draw a mask over an area using the mouse. If necessary, click the mask to discard it, and redraw the mask.




Table 9.1 Data Mask Options

Option	Description
Photograph	If this option is chosen, the software automatically draws the data mask so that it includes the entire photograph.

Table 9.1 Data Mask Options (continued)

Option	Description
Threshold	If necessary use the threshold slider or  arrows to adjust the mask so that it matches the underlying subject photograph as closely as possible without including any area outside the subject image.
Draw Mask	Choose this option to manually draw a data mask on an area of the photograph.
Rectangle	Specifies a rectangular shape for the manual data mask.
Ellipse	Specifies an elliptical shape for the manual data mask.

6. Click the  button to add components to unmix.



NOTE: Two or three components are recommended for the initial automatic analysis. The DyCE results obtained from the automatic analysis can be manually analyzed to identify possible additional components (see [page 168](#) for details on manual analysis).


7. For Cherenkov data only:

- Choose the Decay Correction option to apply decay correction to the image data before analysis.



NOTE: If using Decay Correction, Cherenkov decay correction will be applied to every pixel in image, including pixels where the Cherenkov isotope is not present. Therefore, ensure that the data mask covers only the image region of interest, for example, only the mouse.

- If the radio-isotope used in the experiment was selected in the Imaging Wizard prior to acquisition, it will be displayed in the Isotope drop-down list. If the incorrect radio-isotope was selected at acquisition, choose a different radio-isotope from the drop-down list. If your radio-isotope is not available in the list, choose "Undefined" and enter the isotope half-life in minutes ([Figure 9.16](#)).

Figure 9.16 Auto Unmix Wizard – Cherenkov Data

8. Click **Finish**.
- The Unmixing window shows a plot of the temporal spectra, unmixed images, and a composite of the unmixed images (Figure 9.17).









Table 9.1 Spectrum List Toolbar

Item	Description
	Enables you to view and save the unmixed images as a sequence data set. The image adjust, corrections/filtering, image information, or ROI tools are available for the images.
	Enables you to subtract one spectrum from another (see page 174).
	Adds a temporal component to the spectrum list when performing a manual analysis. See page 168 for more details on manual analysis.
	Deletes the last component in the spectrum list. Click Unmix after deleting a spectrum to view updated DyCE results.

9. To save the results:
- a. Enter a name in the Results tab of the tool palette (Figure 9.18).

b. Click **Save**.




Manual DyCE Analysis

1. Load a DyCE image sequence. Alternatively, load DyCE results obtained from an automatic analysis (Figure 9.19).



NOTE: This section illustrates manual analysis of DyCE results obtained from an automatic analysis.



2. Click the Image Cube tab (Figure 9.20).

The image cube represents a “stack” of the DyCE sequence images. If the Overview option is selected, the image cube shows a composite of all images (Figure 9.20).

To view a particular image, remove the check mark next to Overview and move the slider or enter an image number (Figure 9.21).

Figure 9.20 Image Cube – Overview Mode






Figure 9.21 Image Cube – Single Image Mode



3. Move the mouse pointer over the image cube to see the temporal spectrum at a particular location. The temporal spectrum at the pointer location is updated as you move the pointer.




NOTE: If analyzing DyCE results, the Normalized option for the spectrum plot must be checked to see all of the temporal spectra when the mouse pointer is over the image cube.

4. To add another component to unmix:
 - a. Click the **+** button.
A new name appears in the spectrum list (Figure 9.22)
 - b. Specify the region by using the mouse to draw a mark on the image cube. If necessary, click the  button next to the spectrum name to select a different line thickness from the drop-down list.
 - c. If necessary, right-click the image cube to erase the mark.
5. Repeat [step 4](#) to specify additional temporal components.







NOTE: A maximum of 10 components can be unmixed.

Figure 9.22 Mark the Area of a Temporal Component on the Image Cube




6. Click **Unmix** after you finish marking the regions.
The software generates unmixed images for the new temporal spectra and updates the composite image with these components.

Table 9.2 Spectrum List Toolbar

Item	Description
	Enables you to view and save the unmixed images as a sequence data set. The image adjust, corrections/filtering, image information, or ROI tools are available for the images.
	Enables you to subtract one spectrum from another (see page 174).
	Adds a component to the spectrum list.
	Deletes the last spectrum in the spectrum list.

9.4 DyCE Results

The Unmixing window shows the DyCE results. The example in [Figure 9.23](#) shows three “temporal spectra” (signal as a function of time).



Viewing Unmixed Images




- An unmixed image shows the maximum signal of a temporal spectrum.
- Double-click an unmixed image to view it in an image window ([Figure 9.24](#)). The tool palette is available for viewing and analyzing the image.
 - Click the  button to view the unmixed images as a sequence([Figure 9.24](#)). The tool palette is available for viewing and analyzing the sequence. The software prompts you to save the sequence when closing the Sequence View window.

Figure 9.24 View an Unmixed image in an Image Window or View all Unmixed Images as a Sequence









Viewing the Composite Image

1. Double-click the composite thumbnail.
The Composite window opens.

Figure 9.25 Open the Composite Window


2. Add or remove the check mark next to an image to include or exclude the data from the composite image.
3. Use the image adjust tools at the bottom of the Composite window to adjust the appearance of the composite image.


Table 9.3 Composite Window

Item	Description
	<div><div>Sends the composite image to the “top” of the image cube.</div><div><div><div>Image Cube</div><div>Unmixed</div></div><div></div><div><div>Image Cube Viewer</div><div><input checked="" type="checkbox"/> Overview</div><div><div>1</div><div></div></div></div></div><div><div>Click the Image Cube tab in the Unmixing window to view the image cube. See Figure 9.20 on page 169 for more details on the image cube.</div><div><div><div>Image Cube</div><div>Unmixed</div></div><div></div><div><div>Image Cube Viewer</div><div><input checked="" type="checkbox"/> Overview</div><div><div>1</div><div></div></div></div><div>Composite image displayed on top of the image cube.</div></div></div></div>
	Copies the Composite window to the system clipboard.
	Opens a dialog box that enables you to export the composite image to a graphic file (for example, .png).
	Opens the print dialog box.
<div>Color: </div>	<div><div>Shows the color of the data for the highlighted image.</div><div><div>Image List:</div><div><div><input checked="" type="checkbox"/> Kidney</div><div><input checked="" type="checkbox"/> Brain</div><div><input checked="" type="checkbox"/> Bladder</div><div><input checked="" type="checkbox"/> Photograph</div></div><div></div></div><div>Click the color swatch to open the color palette which can be used choose a color for the selected image data.</div></div>
<div>Label: <input type="text" value="Bladder"/></div>	Data name for the highlighted image. Double-click the name to edit it.

Correcting Temporal Spectra

Temporal spectra can be corrected for overlapping spectra; for example, correcting fluorescence temporal spectra for tissue autofluorescence.

 **NOTE:** If correcting for tissue autofluorescence, one of the unmixed components of the data set should be tissue autofluorescence signal only.

- 1. Click the  button in the Unmixing window.
- 2. In the dialog box that appears, choose the spectra to subtract (Figure 9.26).





Table 9.4 Computed Spectrum

Item	Description
Normalized	Choose this option to normalize the spectra with respect to time zero.
Result: $C = A - x \cdot B$	The subtraction performed by the software where "x" is a factor that ensures the residual signal is positive.
Autoscaling	Choose this option to normalize spectra signal on a scale of zero to one.
Fit Offset	If this option is chosen, the software computes and removes an intensity baseline from the spectra.

Table 9.4 Computed Spectrum (continued)

Item	Description
Error Tolerance	The software computes a default error tolerance (the factor "x" for $A - x \cdot B$) such that signal B is maximally removed from signal A with no negative result. Moving the slider adjusts the error tolerance and automatically updates the computed spectrum.
	<p>Choose "New" to save computed spectrum with the specified name and color. Click Apply to add the computed spectrum to the line plot and spectrum list in the Unmixing window.</p> <p>Choose a spectrum number from the drop-down list to overwrite that spectrum with the computed spectrum when you click Apply.</p>

10 Reconstructing a 3D Surface


Generating a Surface

Managing Surfaces on page 181

Export or Import a Surface on page 182

A *surface* is a 3D reconstruction of the animal surface (topography) derived from a structured light image. A surface is a required input to DLIT or FLIT analyses ([Figure 10.1](#)).

You can also import a surface or export a surface for viewing in other 3D viewer applications



10.1 Generating a Surface


1. Load the image sequence for the reconstruction.
For example, a sequence that was acquired for DLIT analysis.
2. Select an orientation (dorsal or ventral) and subject in the surface topography tools.
3. Select a smoothing level.



NOTE: The default "Low" smoothing level is sufficient in most cases, but it may be necessary to modify this if there are tufts of hair on the animal which disrupt the surface smoothness. Furred mice are not recommended for DLIT or FLIT. If the Surface Topography tool appropriately generates the surface of a furred mouse, the optical data pattern can be grossly shifted by the fur

4. Click **Reconstruct**.


The entire subject is selected for reconstruction by default in the Single View Surface Tomography window that appears (Figure 10.2).




5. If you want to reconstruct only a particular region of the subject, resize the rectangle (drag a green handle ■) so that it includes only the area of interest.
6. Click **Next**.


The purple data mask appears. The mask is an overlay on the subject image that defines the area of interest for the surface topography reconstruction. The mask should match the underlying photograph of the subject as closely as possible without including any area outside the subject image.


Clicking the Show option also displays the pink selected data. The Surface Topography tool usually provides a surface which very closely resembles a surface that is "shrink-wrapped" around the pink mask.















Figure 10.3 Data Mask (purple)

7. If it is necessary, adjust the threshold value so that the mask fits the subject image as closely as possible. To change the threshold, do one of the following:
 - Press the left or right arrow keys on the keyboard.
 - Move the Threshold slider left or right.
 - Click the  arrows or enter a new value in the box.
8. Click **Finish**.

The surface and 3D tools appear in the Tool Palette. For more details on the Tool Palette, see *3D Optical Surface Tools* on page 218.


Figure 10.4 3D View and 3D Tools in the Toolbar and Tool Palette

**Table 10.1** 3D View Tools

Tool	Description
Info	Click to display or hide the image label. The image label includes information about the experiment that was entered in the Edit Image Labels dialog box and other image information automatically recorded by the software.
Image Tools	<p>A drop-down list of tools for viewing and working with the surface.</p> <p>Select  to:</p> <ul style="list-style-type: none"> Click and display measurement dimensions in the coronal, sagittal, or transaxial view (in the 3D view window). Drag a measurement cursor in the coronal, sagittal, or transaxial view and display measurement dimensions. (See page 68 for details on measurement cursors.) <p>Select  to zoom in or out on the image (use a click-and-drag operation).</p> <p>Select  to move the subject in the window (use a click-and-drag operation).</p> <p>Select  to rotate the subject around the x, y, or z axis (use a click-and-drag operation).</p>
	Click to hide or show the x, y, z-axis display in the 3D view window.
	Click to hide or show coronal, sagittal, and transaxial planes through the surface in the 3D view window.
	Click to show or hide a bounding box around the surface.
	Click to show or hide a grid under the surface.
	Select this tool from the drop-down list to change the view perspective (top, bottom, left, right, front, back, or perspective view). For examples of the views, see Figure 10.7 .
	Select this tool from the drop-down list to display the perspective view.
	Click to show or hide measurement cursors in the coronal, sagittal, or transaxial views. Click and drag the green handle (■) at either end of a measurement cursor to resize and reposition it.
	Use this tool to quantify source signal, measure source depth, determine source center of mass, and identify host organ (if an organ atlas is registered with the reconstruction). Load DLIT or FLIT results, click the  button, then draw a box over an area containing voxels in the 3D reconstruction. The source measurements appear in the 3D tools (see Figure 11.28 on page 208).
	Enables you to save the 3D view to a graphic file (for example, .jpg).

Changing the View Perspective

You can click and drag the surface to view it from different perspectives. Alternatively, do one of the following:

- Select  to change the view (Figure 10.6)
- Click the surface in the 3D View window, then press the **V** key to cycle through the different views of the surface
- Figure 10.7 shows examples of the available views. You can view the surface from different perspectives by doing one of the following:





Figure 10.7 Alternative Surface Views

Click the surface, then press the "V" key to change the view.



10.2 Managing Surfaces

After the surface is saved, it can be shared by the DLIT or FLIT tools.


Figure 10.8 Tool Palette – Surface Topography Tools

Table 10.2 Surface Topography Tools – Managing Surfaces

Item	Description
Name	Name of the selected surface.
Delete	Removes the selected surface from the system.
Load	Opens the selected surface.
Save	Saves a surface to the selected name.
Overwrite	Saves the surface and overwrites the previous surface results.

Export or Import a Surface

A surface can be shared with other users or viewed in other 3D viewer applications.



NOTE: Surface import capability is only available if “Show Advanced Options” is selected in the general preferences (see [page 275](#)).

1. Load a surface.
2. Select **File** → **Export (or Import)** → **3D Surface** on the menu bar.
3. In the dialog box that appears, select a folder, enter a file name, and select a file type (see [Table 10.3](#)).



NOTE: Importing a surface by this method is for viewing purposes only, not for registration with optical reconstructions in Living Image software. To import a surface or other organs for registration purposes, import an organ atlas. See [page 227](#) for more details.

Table 10.3 Surface File Types

Export Option	Description	Export	Import
Surface mesh (.xmh)	A native file format of the Living Image software that is used to exchange 3D surface information between Living Image software and other third party analysis tools. It is based on a basic indexed face set format which stores all of the vertex information first, then stores the triangle information in terms of indexes into the vertex list.	yes	yes
AutoCAD DXF (.dxf)	Drawing exchange format that is compatible with most DXF file viewers.	yes	yes
VRML 1.0 (.wrl)	VRML 1.0 (.wrl) - Virtual reality modeling language format that is compatible with most VRML viewers.	yes	no
Open Inventor (.iv)	The ASCII version of the IV file format which is supported by all IV viewers.	yes	yes
STL (.stl or ASCII format)	Stereo lithography binary format compatible with most STL viewers.	yes (binary)	yes

11 3D Reconstruction of Sources

Overview of Reconstructing Sources

Acquire a Sequence on page 186

Steps to Reconstruct Luminescent Sources Using DLIT on page 192

Steps to Reconstruct Fluorescent Sources Using FLIT on page 197

Including or Excluding Data for 3D Reconstruction on page 200

3D Reconstruction Results on page 203

Checking the Reconstruction Quality on page 205

Measuring Sources on page 207

Viewing Luminescent and Fluorescent Sources in One Surface on page 211

Comparing Reconstruction Results on page 211

Exporting a 3D Scene as DICOM on page 216

3D Optical Surface Tools on page 218

3D Optical Source Tools on page 220

3D Optical Registration Tools on page 222

3D Animation on page 229

DLIT/FLIT Troubleshooting on page 234

11.1 Overview of Reconstructing Sources

Living Image software provides algorithms which analyze 2-dimensional optical image data to render 3-dimensional (3D) reconstructions of luminescent or fluorescent sources located inside an animal (*tomographic analysis*). [Figure 11.1](#) shows an overview of the 3D reconstruction workflow.




TIP: See the technical note *DLIT and FLIT Reconstruction of Sources* for more details on the DLIT or FLIT algorithm (select **Help** → **Tech Notes** on the menu bar).


3D Reconstruction Algorithm	Description	See Page
Diffuse Tomography (DLIT)	DLIT provides a complete 3D reconstruction of the luminescent source distribution within the subject. DLIT places no constraints on the geometry or spatial variation of the source strength throughout the volume. DLIT is well-suited for analyzing complex and spatially extended luminescent sources. The 3D reconstruction is presented as voxels. If a luminescent quantification database is available, the number of cells per source can be determined in addition to source intensity (photons/sec).	192
Fluorescent Tomography (FLIT)	FLIT provides a complete 3D reconstruction of the fluorescent source distribution within the subject. The 3D reconstruction is presented as voxels. If a fluorescent quantification database is available, the number of fluorophore molecules or cells per source can be determined in addition to the total fluorescence yield.	197

Figure 11.1 Basic 3D Reconstruction Workflow



1. Set up a DLIT or FLIT sequence using the Imaging Wizard (see [page 186](#)). Acquire and load the sequence.




2. Generate or load a surface using the Surface Topography tools. See [Chapter 10](#) on [page 176](#) for more details.



3. In the DLIT or FLIT 3D Reconstruction tools, select the:
 - Wavelengths or excitation point images to analyze.
 - Tissue and source properties.
4. Reconstruct sources.
See [page 192](#) for detailed DLIT steps.
See [page 197](#) for detailed FLIT steps.



5. View source measurements (see [page 207](#)).



General Considerations

Animal Requirements

The best surface topography reconstruction is obtained from nude mice. It is possible to perform 3D imaging on white or light-colored furred mice if the fur is reasonably smooth over the mouse surface. Therefore it is recommended that you comb the fur before imaging to eliminate any "fluffy" areas that may alter the light emission pattern and/or trigger artifacts during the surface topography reconstruction. In this case, it is recommended that you shave the animals or apply a depilatory. 3D reconstructions using mice with black or dark-colored fur will give poor results.

Luminescent Exposure vs. Luciferin Kinetic Profile

It is important to consider the luciferin kinetic profile when you plan the image sequence acquisition. The DLIT algorithm currently assumes a stable luciferin kinetic profile. Therefore, to optimize the signal for DLIT 3D reconstruction, carefully plan the start and finish of image acquisition and ration the exposure time at each emission filter so that the sequence is acquired during the flattest region of the luciferin kinetic profile.

Reconstruction Inputs



NOTE: Use the Imaging Wizard to set up the DLIT or FLIT image sequence. See [Acquire a Sequence on page 186](#) for more details.

DLIT

Input data to the DLIT algorithm for a 3D reconstruction of luminescent light sources includes:

- A surface topography of the subject (generated from a structured light image).
- A sequence of two or more images of the light emission from the subject surface that is acquired at different filter bandpasses ([Table 11.1](#)).

Table 11.1 IVIS® Spectrum Filters for Luminescence or Fluorescence Tomography

Filters	Range	Bandwidth
10 excitation filters	415-760 nm	30 nm
18 emission filters	490-850 nm	20 nm

FLIT

Input data to the FLIT algorithm for 3D reconstruction of fluorescent light sources includes:

- A surface topography of the subject (generated from a structured light image).
- A sequence of images acquired at different transillumination excitation source positions using the same excitation and emission filter at each position ([Table 11.1](#)).

Quantification Database (Optional)

If a quantification database is available, it is possible to determine the number of cells in a DLIT source or the number of cells or dye molecules in a FLIT source. The database is derived from an analysis of images of known serial dilutions of luminescent or fluorescent cells, or dye molecules in a well plate.

See [Chapter 12 on page 235](#) for more details on generating a database. Using a quantification database is optional.

11.2 Acquire a Sequence

Use the Imaging Wizard to set up a DLIT (see below) or FLIT (see [page 188](#)) sequence. See [Manual Sequence Setup on page 191](#) if not using the wizard.


Bioluminescence Imaging




NOTE: The IVIS Spectrum should be initialized and the temperature locked before setting the imaging parameters. See [page 7](#) for more details.

1. Start the Imaging Wizard. See *Starting the Imaging Wizard on page 40* for instructions.
2. Double-click the Bioluminescence option. Double-click the DLIT option in the next screen (Figure 11.2).


Figure 11.2 Imaging Wizard – Choose Bioluminescence and DLIT Options



3. Select a probe from the Name drop-down list and click **Next** (Figure 11.3).




4. In the next screen (Figure 11.4):
 - a. Select the type of imaging subject.
 - b. Choose the Auto Settings option for the exposure parameters.
 - c. Select a field of view from the drop-down list.



5. Set the focus by doing either of the following:
 - Enter a subject height and choose the “use subject height” focus option.
 OR
 - Choose the “manual focus” option from the Focus drop-down list and set the focus parameters in the Manual Focus Window that appears.


6. If performing a time series study, choose the Time Series Study option. Enter the number of segments to acquire and the delay between segment acquisition. A segment is an image sequence.
7. Click **Next**.
The specified sequence appears in the sequence table (Figure 11.5).

Figure 11.5 Sequence Setup Complete



8. Acquire the sequence following the instructions on page 41.
9. The image window appears when acquisition is completed (Figure 11.6). See Table 3.2 on page 28 for more details on the Image window.

Figure 11.6 Image Window – Sequence View




Fluorescence Imaging




NOTE: The IVIS Spectrum should be initialized and the temperature locked before setting the imaging parameters. See page 7 for more details.


1. Start the Imaging Wizard. See *Starting the Imaging Wizard* on page 40 for instructions.
2. Double-click the Fluorescence option. Double-click the FLIT option in the next screen (Figure 11.7).

Figure 11.7 Imaging Wizard – Choose Fluorescence and Flit Options


3. Select a probe from the Name drop-down list and click **Next** (Figure 11.3).

Figure 11.8 Select a Probe


4. In the next screen (Figure 11.4):
 - a. Select the type of imaging subject.
 - b. Choose the Auto Settings option for the exposure parameters.
 - c. Select a field of view from the drop-down list.

Figure 11.9 Imaging Wizard – Fluorescence FLIT

5. Set the focus by doing either of the following:
 - Enter a subject height and choose the “use subject height” focus option.
 OR
 - Choose the “manual focus” option from the Focus drop-down list and set the focus parameters in the Manual Focus Window that appears.
6. If performing a time series study, choose the Time Series Study option. Enter the number of segments to acquire and the delay between segment acquisition. A segment is an image sequence.
7. Select the transillumination locations.
 - a. Click **Transillumination Setup**.
 - b. Choose the transillumination locations in the Transillumination Setup box that appears (Figure 11.10).
 See [Table 3.3 on page 37](#) for more details on Transillumination Setup.

Figure 11.10 Transillumination Setup Box

8. Click **Next**.
- The specified sequence appears in the sequence table (Figure 11.5).



9. Acquire the sequence following the instructions on page 41.
- The image window appears when acquisition is completed (Figure 11.6). See Table 3.2 on page 28 for more details on the Image window.

Manual Sequence Setup

This sections provides sequence requirements if you will not be using the Imaging Wizard and plan to manually set up the sequence.

Table 11.2 IVIS® Spectrum Filters for Luminescence or Fluorescence Tomography

Filters	Range	Bandwidth
10 excitation filters	415-760 nm	30 nm
18 emission filters	490-850 nm	20 nm

DLIT Sequence Requirements

A sequence must include:

- A structured light image
- Optical data from at least two different emission filters (560 - 660 nm), at a minimum:
 - Emission filter #1: Photographic, luminescent
 - Emission filter #2: Luminescent image

Analyzing more optical images usually produces more accurate results. Table 11.3 shows the recommended optical image sequence.

Table 11.3 Manual Sequence Setup – Recommended DLIT Optical Image Sequence

Image Type	Emission Filter Options					
	560	580	600	620	640	660
Photograph	✓	Select the Reuse option in the control panel.				
Luminescent	✓	✓	✓	✓	✓	✓




NOTE: It is recommended to set the binning level the same for all of the luminescent images.

FLIT Sequence Requirements

Use transillumination on the IVIS® Spectrum and the same excitation and emission filters from at least four source locations that form a rectangle. Acquire the following images:

- Fluorescent image and photograph at the first transillumination location
- Fluorescent image at the remaining transillumination locations
- A structured light image


Figure 11.12 shows an example image sequence.




11.3 Steps to Reconstruct Luminescent Sources Using DLIT

1. Load a DLIT image sequence.
2. Generate or load a surface using the Surface Topography tools. For details on generating the surface, see [Chapter 10 on page 176](#).
3. In the Tool Palette, choose **DLIT 3D Reconstruction**.

The Analyze tab shows the data that the algorithm automatically selects for the reconstruction ([Figure 11.13](#)). For more details about the Threshold %, see [page 200](#).




4. In the Properties tab, make a selection from the “Tissue Properties” and “Source Spectrum” drop-down lists (Figure 11.14).



5. To view the tissue properties (μ_a , μ_{eff} , μ'_s) for the tissue and source you selected, make a selection from the Plot drop-down.
6. Select a luminescent quantification database to compute the number of cells per source (optional).
For details on generating a luminescent quantification database, see [page 235](#).
7. In the Analyze tab, click **Start**.

The Data Preview window appears and displays the image data that will be included in the reconstruction. Usually, no data adjustment is required. However, it is possible to exclude or include user-selected pixel data from the analysis. See [Including or Excluding Data for 3D Reconstruction on page 200](#) for more details.

Figure 11.15 Data Preview Window



8. In the Data Preview window, click **Reconstruct.**

The reconstruction usually requires less than one minute, depending on the reconstruction volume, parameter settings, and computer performance. When the analysis is finished:

- The 3D View window displays the animal surface and the reconstructed sources.
- In the Tool Palette, the Results tab displays the results data and the algorithm parameter values. See [page 204](#) for details on managing results (for example, save, load, or delete).
- The 3D Optical Tools appear after a reconstruction is generated or loaded. Use these tools to modify the source display parameters. See [Table 11.4 on page 195](#) for an overview of the tools; more details are available on [page 218](#) to [page 229](#).

Table 11.4 3D Optical Tools






















Tool Palette	3D Optical Tools	Functions	See Page
	Surface Tools	Adjust the appearance of the reconstructed animal surface and photon density maps	218
	Source Tools	Adjust the appearance of reconstructed sources, make source measurements, export voxel measurements	220
	Registration Tools	Display organs on the reconstructed surface, adjust the location or scale of organs on the surface, import an organ atlas	222
	Animate Tools Select Tools → 3D Animation on the menu bar to view these tools.	Display preset animations of the 3D View scene. Enables you to create custom animations and record an animation to a movie file.	229

Figure 11.16 DLIT Reconstruction Results




**Table 11.5** 3D View Tools

Tool	Description
Info	Click to show or hide the image label. The image label includes information you enter in the Edit Image Labels dialog box (see Figure 3.5 on page 26) and other information automatically recorded by the software.
Image Tools 	<p>A drop-down list of tools for viewing and working with the surface or DLIT results.</p> <p> or  - Rotates or spins the surface in the x, y, or z-axis direction.</p> <p> - Moves the surface in the x or y-axis direction.</p> <p> - Zooms in or out on the image. To zoom in, right-click (Cmd key (apple key) +click for Macintosh users) and drag the  toward the bottom of the window. To zoom out, right-click and drag the  toward the top of the window.</p>
	Displays the x,y,z-axis display in the 3D view window.
	Displays coronal, sagittal, and transaxial cross-sections through the subject in the 3D view window.
	Displays a bounding box around the subject.
	Displays a grid under the subject.
	Select this tool from the drop-down list to change the view perspective (top, bottom, left, right, front, back, or perspective view). For examples of the views, see Figure 11.46 on page 227 .
	Select this tool from the drop-down list to display the perspective view.
	Rotates the 3D reconstruction results in the 3D view window (<i>3D scene</i>). Click the + or - key to increase or decrease the rotation speed. To stop the rotation, click the 3D scene or the  button.
	Displays measurement cursors in the coronal, sagittal, or transaxial views.
	Click this button, then select a source or a point in a source to obtain source measurements (total flux, volume, center of mass, host organ) in the 3D Optical Tools (Source tab). For more details, see page 207 .
	Copies or pastes voxels or a source surface so that DLIT and FLIT reconstructions can be displayed on one surface. For more details, see page 211 .
	Enables you to save the 3D view to a graphic file (for example, .jpg).


11.4 Steps to Reconstruct Fluorescent Sources Using FLIT

1. Load a FLIT image sequence.
2. Generate or load a surface in the Surface Topography tools. For details on generating the surface, see [Chapter 10 on page 176](#).
3. In the Tool Palette, choose **FLIT 3D Reconstruction**.

The Analyze tab shows the images that the algorithm automatically selects for the reconstruction based on an appropriate signal level ([Figure 11.13](#)). For more details about the Threshold %, see [page 200](#).




4. Select the type of image used in the reconstruction: Radiance or NTF Efficiency ([Figure 11.18](#)). NTF Efficiency data is the default because it affords higher sensitivity to the embedded fluorescence sources.
5. Make a selection from the “Tissue Properties” list in the Properties tab ([Figure 11.19](#)).



6. To view the tissue properties (μ_a , μ_{eff} , μ_s) for the tissue you selected, make a selection from the Plot drop-down.
7. Optional: Select a fluorescent quantification database to include the number of fluorescent molecules/source in the results.
For details on generating a fluorescent quantification database, see [page 235](#).
8. In the Analyze tab, click **Start**.
9. The Data Preview window appears and displays the image data that will be included in the reconstruction.

Usually, no data adjustment is required. However, it is possible to exclude or include user-selected pixel data from the analysis. See [page 200](#) for more details.


You can also include or exclude image data by adding or removing the check mark next to the images listed in the Analyze tab ([Figure 11.18](#)).



10. Click **Reconstruct.**

The reconstruction normally requires less than one minute, depending on the reconstruction volume, parameter settings, and computer performance. When the analysis is finished:

- The 3D View window displays the surface and the reconstructed sources. See [page 204](#) for details on managing results (for example, save, load, or delete).
- In the Tool Palette, the Results tab displays the results data and the algorithm parameter values ([Figure 11.24](#)). See [page 204](#) for details on managing results (for example, save, load, or delete).
- The 3D Optical Tools appear in the Tool Palette. Use these tools to modify the source display parameters. See [Table 11.4 on page 195](#) for an overview of the tools; more details are available on [page 218](#) to [page 229](#).



11.5 Including or Excluding Data for 3D Reconstruction

The Data Preview window shows the image data that are automatically selected for reconstruction ([Figure 11.22](#)). In special cases, you may want to include or exclude particular data from this default selection. There are two ways to do this:

- Change the Threshold % value (see below) – Applying a Threshold % value excludes or includes some pixels from the reconstruction. The software computes the minimum and maximum pixel values of an image based on an histogram of pixel intensities. If Threshold % = 0.5%, then pixels with intensity less than 0.5% of the maximum intensity value are excluded from the reconstruction. The Threshold % can be edited for individual images. The Data Preview window is updated when you change the Threshold % value. Min Counts translates the Threshold % to the minimum counts required for reconstruction. Keep the minimum counts > 200.
- Region selection (see [page 201](#)) – Use the pencil tool to mark particular regions to include in the reconstruction. This may be useful for noisy images with high intensity pixels where changing the Threshold % value is not helpful. You can also use this method to focus on particular sources to reconstruct and ignore others.


To change the Threshold % for a selected image:


1. Click **Start** in the Analyze tab ([Figure 11.22](#)).
The Data Preview window appears.
2. Click an image in the Data Preview window.



NOTE: Changes to Threshold % are applied to the selected image only. To apply the change to all images, choose the **Select All** option.

3. Click **Data Adjustment**.
4. In the window that appears, enter a new Threshold % value.
The new Threshold % appears in the Analyze tab.
5. To reset the Threshold % to the default value (for the selected images), click **Restore Threshold**.


Figure 11.22 Adjusting Threshold %**To select particular regions for reconstruction:**

1. Open the Data Preview window as shown in [Figure 11.22](#).
2. Click **Data Adjustment**.
3. In the window that appears, choose the **Draw** option and put the mouse pointer over the image so that the pencil tool  appears.
4. To automatically select all pixels in a source, right-click with the region with the pencil tool. Alternatively, put the pencil over the image and click the mouse key or press and hold the mouse key while moving the pencil over an area of the image.



NOTE: If the pencil tool markings are applied to the image, only the marked pixels are included in the analysis.

Figure 11.23 Selecting Regions to Include in Reconstruction



Use these tools to select particular image data to include in the analysis. See [Table 11.6](#) for details on the tools.



Choose the "Draw" option then mark the area to include in the reconstruction using the pencil tool . In this example, the red area marked with the pencil tool will be reconstructed.

Table 11.6 Region Selection Tools

Item	Description
Draw	Choose this option to display the pencil tool  when the mouse pointer is over the data adjustment image. Use this tool to apply markings that select regions to include in the reconstruction.
Erase	Choose this option to display the eraser tool. Use the eraser to remove pencil tool markings (exclude pixels from the image).
Painting size	Adjusts the width of the pencil tool mark or the eraser tool.
Segment	Colors available for the pencil tool.
Opacity	Adjusts the opacity of the pencil tool markings.
Reset	Removes all pencil tool markings.

11.6 3D Reconstruction Results

The Results tab displays information about the photon density, voxels, and algorithm parameters.

DLIT or FLIT Results



NOTE: For more details on DLIT, see the reference article *DLIT and FLIT Reconstruction of Sources* (select **Help** → **References** on the menu bar). Sometimes adjusting the DLIT algorithm parameters improves the fit of the simulated photon density to the measured photon density data.




Table 11.7 DLIT or FLIT 3D Reconstruction Results

Item	Description
Final voxel size (mm)	The voxel size (length of a side, mm) that produces the optimum solution to the DLIT or FLIT analysis.
Number of voxels	The number of voxels that describe the light source(s).
Reduced Chi2	A measure of the difference between the computed and measured photon density maps at the optimum solution. A smaller χ^2 value indicates a better quality of fit.
Index of Refraction	Refractive index of light for the imaged subject.
Angle Limit(deg)	Angle limit of surface normal to optical axis, above which data will not be used in the reconstruction.
Damping reduce	The damping parameter is calculated from this reduction factor, relative to the maximum singular value of the system matrix.
Data range	For multi-view data, the image views used in the reconstruction.
Mirror XOffset	For multi-view data, the mirror location from the x center line.
Starting voxel size	The voxel size at the start of the analysis. The length of the side of the voxel cube in mm units for the coarsest initial grid size in the adaptive gridding scheme.

Table 11.7 DLIT or FLIT 3D Reconstruction Results (continued)


Item	Description
Total # of data pts	The total number of data points used in the reconstruction.
Median Filter	Indicates whether or not a median filter was applied to the data.
Image Threshold	The percentage of the minimum radiance at each wavelength (DLIT) or source location (FLIT) is of the maximum radiance. This defines the minimum intensity included in the data.
Samples of Image	The data in each image is sampled. This parameter shows the number of pixels sampled from each image.
Tissue Properties	The tissue properties for modeling the photon propagation.
Source Spectrum	The emission spectrum of the type of luminescent source.
Quantification Selection	A user-selected quantification database used in the reconstruction to convert reconstruction voxel units to 'cells' or 'picomoles' units.
Sequence name	Image data sequence name.
Version	Living Image® software version

Managing 3D Reconstruction Results

Item in the DLIT/FLIT 3D Reconstruction Results Tab	Description
Name	The name for the active DLIT or FLIT results. Select results from this drop-down list.
Delete	Deletes the selected DLIT or FLIT results.
Load	Opens the selected reconstruction results in the 3D View.
Save	Saves the active DLIT or FLIT results to the selected name. The results are saved to the sequence click number folder and are available in the Name drop-down list.
Overwrite	If you reanalyze saved results, saves the new results and overwrites the previous results.
Export Results	Saves the results to a .csv file.

Copying Results to the System Clipboard

1. To copy all results:
 - a. Right-click the results and chose **Select All** from the shortcut menu.
 - b. Right-click the results again and select **Copy** from the shortcut menu.



2. To copy user-selected results:
 - a. Select the results.
 - b. Right-click the selection and choose **Copy** from the shortcut menu.

11.7 Checking the Reconstruction Quality

Comparing the measured and simulated photon density plots is a useful way to check the quality of a 3D reconstruction.

The photon density is closely related to the measured radiance. Photon density is the steady state measure of the number of photons in a cubic millimeter. Light sources inside the tissue contribute to photon density in other portions of the tissue.


The reconstruction algorithm first converts the luminescent or fluorescent image of surface radiance to photon density just inside the animal surface because this is what can be observed. The algorithm then solves for intensity values at locations inside the tissue which would produce the observed photon density near the surface.

For fluorescence reconstructions using NTF Efficiency data, the photon density of the fluorescence image is divided by the photon density of the transmission image, giving the NTF Efficiency. The NTF Efficiency values are the data just inside the animal surface for this type of data set.

Viewing Photon Density or NTF Efficiency Maps

1. After the reconstruction is finished or results are loaded, click **Photon Density or NTF Efficiency Maps** in the Results tab.
The photon density maps or NTF Efficiency maps for all image data are displayed (Figure 11.26).
2. To rotate the surface and view it from a different angle, move the thumb wheel to the left or right.


Figure 11.26 Photon Density Maps (or NTF Efficiency Maps for Fluorescence)



3. Select a wavelength from the drop-down list

The photon density or NTF Efficiency profiles at the crosshairs location are displayed. In a good reconstruction, the simulated photon density or NTF Efficiency curves (red) closely resemble the measured photon density or NTF Efficiency curves (blue).

Figure 11.27 Simulated (Red) and Measured (Blue) Photon Density (or NTF Efficiency) Plots



4. To view the photon density or NTF Efficiency profile at another location on the animal surface, drag the cross hairs or click a point on the photon density or NTF Efficiency map.

Table 11.8 Photon Density Maps Window

Item	Description
Image sources	A list of images used in the reconstruction. Select all images or a particular image number to display.
Angle of View	The thumb wheel position. Turn the thumb wheel to rotate the surface on the vertical axis.
Log Scale	Choose this option to display the photon density or NTF Efficiency using a log scale.
Simulated	The photon density or NTF Efficiency computed from DLIT or FLIT source solutions which best fit the measured photon density or NTF Efficiency.
Measured	The photon density or NTF Efficiency determined from the image measurements of surface radiance.
Horizontal Profile	The photon density or NTF Efficiency line profile at the horizontal plane through the subject at the crosshairs location.
Vertical Profile	The photon density or NTF Efficiency line profile at the vertical plane through the subject at the crosshairs location.
Position (mm)	Horizontal Profile: The y-axis position of the crosshairs horizontal line. Vertical Profile: The x-axis position of the crosshairs vertical line. The x-y positions are relative to the center of the FOV (where x = 0 and y = 0).

11.8 Measuring Sources

This section presents a convenient way to measure the source (voxels): total flux or total florescence yield, or if calibrated, the abundance in cells or picomoles after the reconstruction is finished or results are loaded.


The volume, center of mass, and depth at the center of mass are also reported in the 3D Optical Tools, Source tab.



NOTE: If the surface contains voxels pasted from other reconstruction results, choose a source in the 3D Source tools ([Figure 11.28](#)). For more details on pasting voxels, see [page 211](#).



Determining the Source Center of Mass


Follow the steps in [Figure 11.28](#) after reconstruction is finished or results are loaded to determine the source center of mass. Alternatively, use the 3D ROI tool for more precise measurements. See [page 114](#) for more details on 3D ROIs.



Measuring Source Depth

Follow the steps below after reconstruction is finished or results are loaded to measure source depth.

1. If the surface includes voxels pasted from other results, select a source from the drop-down list.
2. Confirm that “Display Voxels” is selected, not “Display Source Surface”.
3. Click the Measurement Cursor toolbar button .
The distance from the center of mass to the surface is measured in the three planes.
 - Coronal and transaxial planes display the shortest distance from the center of mass to the surface.
 - The sagittal plane displays the distance from the center of mass to the bottom of the subject.
4. Click the  button to display slice planes through the center of mass. See [page 210](#) for more information on planes.




Viewing Location Coordinates

Click a location in the reconstruction slice in the Coronal, Sagittal, or Transaxial windowpane.

The coordinates (mm) of the position are displayed ([Figure 11.30](#)). The coordinates are updated when you press and hold the mouse button while you drag the cursor.

Slice Plane	Displays...
Coronal	The x-y coordinates of a position.
Sagittal	The y-z coordinates of a position.
Transaxial	The x-z coordinates of a position.

Figure 11.30 Viewing Y,Z Coordinates in the Sagittal Plane

Displaying Slices Through a Reconstruction







1. Click a location on a source. Alternatively, click the  toolbar button, draw a box around a source, then click **Center of mass** in the 3D Source tools.
2. Click the  toolbar button.
The Coronal, Sagittal, and Transaxial windowpanes show a slice through the surface taken by the associated plane.

Figure 11.31 Planes Cutting a Reconstruction

3. To move a plane, put the mouse cursor over a line in the coronal, sagittal, or transaxial windowpane. When the cursor becomes a  or  arrow, drag the line.
The view is updated in the windowpanes as you move the line.





11.9 Viewing Luminescent and Fluorescent Sources in One Surface

When an experiment includes luminescent and fluorescent reporters, DLIT and FLIT reconstructions can be displayed in one surface if the luminescent and fluorescent imaging is done in the same imaging session, without moving the animal.



NOTE: If the DLIT and FLIT image sequences are acquired during the same session, the generated surfaces are nearly identical.

1. Load a DLIT reconstruction and a FLIT reconstruction.
2. Choose one of the reconstructions, click the  button and select **Copy source voxels**.
3. In the other reconstruction, click the  button and choose **Paste source voxels**.



NOTE: Pasted voxels can be measured. For more details on measuring sources, see [page 207](#).

11.10 Comparing Reconstruction Results

Multiple DLIT or FLIT reconstruction results can be viewed side-by-side in the Longitudinal Study window. Voxel intensity within the entire surface or a user-selected area can be measured in all results in the Longitudinal Study window.

The Longitudinal Study window provides a convenient way to compare different results, for example, results obtained at different time points or results from different types of reporters.




NOTE: The FLIT results selected for display in the Longitudinal Study window must have the same type of units. The DLIT results selected for display in the Longitudinal Study window must have the same type of units.


Viewing Results in the Longitudinal Study Window





NOTE: The Longitudinal Study window can display FLIT results or DLIT results, but not both at the same time. Only 3D reconstruction results with the same type of units can be loaded.

1. Load the DLIT or FLIT sequences with the results that you want to display. Select **Tools** → **Longitudinal Study** on the menu bar.
The Longitudinal Study window appears.






NOTE: After the Longitudinal Study window is open, more sequences can be added to the window by clicking the Open button  and selecting sequenceinfo.txt files (found in the sequence data folder).

2. To show particular results:
 - a. Select a sequence in the upper box.
 - b. Select one or more analysis results in the lower box. To choose multiple adjacent results, press and hold the Shift key while you click the first and last result. To choose non-adjacent results, press and hold the Ctrl key while you click the results.
 - c. Click **Load**.
3. To show more results, repeat step [step 2](#)
4. To remove results from the Longitudinal Study window, right-click a surface and select **Remove** on the shortcut menu. Alternatively, select a surface, click the Remove button  and choose **Selected Result**.
To remove all results, click the Remove button  and choose **All Results**.
5. To view a particular image in a sequence:
 - a. Click the surface.
 - b. For DLIT results, make a selection from the Wavelength drop-down list. For FLIT results, make a selection from the Image drop-down list.



**Table 11.10** Longitudinal Study Window


Item	Description
<p>Unit [DLIT]</p> <p><input checked="" type="radio"/> photons/sec</p> <p><input type="radio"/> cells</p> <p>Unit [FLIT]</p> <p><input checked="" type="radio"/> pmol M⁻¹ cm⁻¹</p> <p><input type="radio"/> pmol</p> <p><input type="radio"/> cells</p>	<p>DLIT – Select photons/sec or cells (results calibrated using a quantification database).</p> <p>FLIT – Select pmole M⁻¹ cm⁻¹ or pmoles (results calibrated using a quantification database). Only DLIT results of the selected unit will be enabled for loading in the Longitudinal Window.</p>
<p><input checked="" type="checkbox"/> Display Voxels [DLIT]</p> <p>Red</p> <p><input type="checkbox"/> Reverse <input type="checkbox"/> Log Scale <input checked="" type="checkbox"/> MIP</p>	<p>Voxel display controls:</p> <p>Display Voxels – Choose this option to show voxels within the surface.</p> <p>From the drop-down list, select a color scheme for the color scale. Move the sliders to adjust the color scale minimum and maximum values.</p> <p>Reverse – Choose this option to apply the colors of the selected color table in reverse order to the photon density scale. For example, the Red color table represents the source intensity (photons/sec) from low to high using a color scale from transparent to red. If Reverse is chosen, the source intensity (photons/sec) from low to high is represented using the color scale from red to transparent.</p> <p>Log Scale – Applies a log scale to the color scale.</p> <p>MIP – When this option is chosen, all maximum intensity voxels in the view are projected along the viewing direction into the viewing plane.</p>

Table 11.10 Longitudinal Study Window (continued)

Item	Description
	Copies the 3D View tab in the Longitudinal Study window to the system clipboard.
	Opens a dialog box that enables you to export the 3D View tab to a graphic file (for example, .png).
	Enables you to select voxels for measurement. Measurements are displayed in the Plots tab.

Measuring Intensity

1. Load 3D reconstruction results and click the  button.
By default, a selection box appears around each surface (Figure 11.35). This means that measurements for the entire surface will be computed.
2. To select a particular region of the surface for measurements, draw a box (by clicking and dragging the mouse) around the area.
The same box is applied to the other surfaces in the Longitudinal Study window.
3. To clear boxes, click the  button again.




Viewing Plots

To view a graph, make a selection from the Analysis Type and Plot drop-down lists in the Plots tab (Figure 11.36).


The following graphs are available in the Plots tab:

Plot Type	Description
Quantification Profile	Plots the measured intensity within the user-selected area on the surface. If no box was drawn on the surface, measures the total intensity for the entire surface.
Reduced Chi-Squared Profile	A measure of the difference between the computed and measured photon density maps at the optimum solution. A smaller χ^2 value indicates a better quality of fit.
Voxel Size	Plots the voxel size at the start of the 3D reconstruction and at the end of the 3D reconstruction.



11.11 Exporting a 3D Scene as DICOM


The items in the perspective 3D View are called a *3D scene*. For example, the 3D scene in [Figure 11.37](#) includes a surface and voxels. The 3D scene can be exported to DICOM format and viewed in the Living Image DICOM Viewer or third party software.



To export the 3D scene:

1. Load the results that you want to export.
2. Select **File** → **Export** → **3D Scene as DICOM** on the menu bar.
3. In the dialog box that appears, set the export options, and click **Export**.

For more details on the 3D Scene Exporter, see [Table 11.11](#).



4. In the Browse For Folder dialog box that appears, choose a folder for the DICOM files and click **OK**.

During the export operation, the 3D View window displays the each slice in the export. For example, if Transaxial Slice is selected for export, then the transaxial windowpane cycles through a display of each exported slice.

Table 11.11 3D Scene Exporter Dialog Box

Item	Description
Save DICOM as:	Single-Frame DICOMs - Exports multiple files that contain a single frame each. Multi-Frame DICOM - Exports a single file that contains multiple frames. Note: Choose the Single-Frame or Multi-Frame DICOM option, depending on the third party software you will use to import and view the 3D scene. Some applications cannot reconstruct multi-frame DICOM files.
Slice Orientation	Choose transaxial, coronal, or sagittal slices for the export.
Export voxels using original resolution	Choose this option to export source voxels without any smoothing or binning. The original resolution of the source voxels is the resolution obtained after DLIT or FLIT reconstruction (approximately 1mm resolution).
Slice Resolution	Sets the number of slices required to accommodate the slice orientation with good slice sampling/spacing.
Total Slices Slice spacing Pixel spacing	Parameters that determine the number and resolution of the slices to export.
Solid mesh	If this option is chosen, voxels generated inside the hollow mesh are assigned an intensity so that they are displayed as “tissue” when loaded into visualization software. If no intensity is associated with the voxels, they are considered noise or air and appear hollow.
Hollow mesh	The intensity of pixels inside the surface is set to zero so that the exported surface appears as a hollow empty structure.

Viewing DICOM Data

The 3D scenes exported to DICOM can be viewed in the Living Image® 3D Browser.

1. Select **File** → **Browse 3D Volumetric Data** on the menu bar.
2. In the dialog box that appears, select the DICOM data (.dcm or .dc3) and click **Open**.
The 3D Browser window appears.


Figure 11.39 Living Image 3D Browser

Table 11.12 Living Image 3D Browser DICOM Viewing Controls

Item	Description
Start Index	Specifies the first image (slice) for viewing.
Auto Preview	Select this option to automatically play back the images.
End Index	Specifies the last image (slice) for viewing.
Load	Opens the DICOM data in a 3D View window.
Load data in new window	If this option is selected, DICOM data are opened in a new 3D View window when you click Load . If this option is not selected, DICOM data are loaded in the active 3D View window.

11.12 3D Optical Surface Tools

Use the Surface tools to adjust the appearance of the reconstructed animal surface and photon density maps.

Figure 11.40 Surface Tools and Example DLIT Reconstruction with Photon Density or NTF Efficiency Maps























Table 11.13 3D Surface Tools

Item	Description
Display Subject Surface	Choose this option to display the surface in the 3D View window.
	<p>Drawing styles for the surface.</p> <div>     </div> <p>  Point cloud  Wire frame  Surface face  Wire frame & surface face </p>
	<p>Shading styles for the surface.</p> <div>     </div> <p>  Surface face  Smooth surface face  Reflect surface face  Reflect smooth surface face </p>
	Click to open the color palette from which you can select a display color for the surface and the cross section views.
Opacity	Adjusts the surface opacity.
Display Photon Density or NTF Efficiency Map	Choose this option to display the photon density or NTF Efficiency on the surface.
Apply	Choose measured or simulated photon density or NTF Efficiency maps for display.
Wavelengths (DLIT) Images (FLIT)	Choose the data to display in the photon density or NTF Efficiency map.
Intensity	Set the maximum intensity of the photon density or NTF Efficiency map using the slider or by entering a value.
Color Table	Color scheme for the photon density or NTF Efficiency map.
Reverse	Choose this option to apply the colors of the selected color table in reverse order. For example, the Red color table represents the mapped intensity from low to high using a color scale from transparent to red. If Reverse is chosen, the mapped intensity from low to high is represented using the color scale from red to transparent.
Log Scale	Choose this option to apply a logarithmic scale to the photon density or NTF Efficiency scale.

11.13 3D Optical Source Tools

Use the Source tools to:

- Adjust the appearance of sources in DLIT or FLIT reconstructions
- Make source measurements (page 207)
- Export voxel measurements (.csv)




Table 11.14 3D Source Tools



Item	Description
Select Source	A drop-down list of available sources. Original – Results saved with the data. <sequence name...SourceVoxels> – Pasted voxels. (Click the  button to remove pasted voxels from the surface.) See Viewing Luminescent and Fluorescent Sources in One Surface on page 211 for more details on copying and pasting sources from one sequence to another.
Display Source Surface	Choose this option to display the source surfaces reconstructed using DLIT or FLIT. A surface will be wrapped around the currently displayed voxels. Adjust the voxel display by moving the Threshold slider.
	Drawing styles for the source surface (see “ Display Source Surface ”).

Table 11.14 3D Source Tools (continued)








Item	Description
	Shading styles for the source surface (see " Display Source Surface ").
	Click to open the color palette from which you can select a display color for the source surface.
Opacity	Adjusts the source surface opacity.
Display Voxels	Choose this option to display the sources reconstructed using DLIT or FLIT.
Maximum Intensity Projection	Choose this option to project all maximum intensity voxels in the view along the viewing direction into the viewing plane.
Threshold (DLIT/FLIT)	Choose this option to apply a minimum threshold intensity to the voxel display.
Gradation (DLIT/FLIT)	Use this slider to set a threshold for the percentage voxel intensity above which voxels are opaque and below which voxels will gradually fade to transparent. The percentage voxel intensity is the percentage relative to the maximum intensity.
Voxel size	The 3D grid-spacing size for interpolation of the reconstructed source.
Smoothing	The smoothing box filter size.
Display voxels as	The voxel display mode (cubes, spheres, points, or texture).
Color Scale and Color Table	 <p>Voxel Min: Use the slider, up/down arrows, or enter a value to set the minimum value of the source color scale. Voxels with intensities less than the color scale minimum are not displayed in the reconstruction.</p> <p>Color Table – Color scheme for voxel display. Use the left and right sliders, up/down arrows, or enter values to set the minimum and maximum colors. The Color Table Max is allowed up to 100 times the voxel maximum.</p> <p>Reverse – Choose this option to apply the colors of the selected color table in reverse order to the source voxel scale. For example, the Red color table represents the source intensity from low to high using a color scale from transparent to red. If Reverse is chosen, the source intensity from low to high is represented using the color scale from red to transparent.</p> <p>Log scale – Choose this option to apply a logarithmic scale to the color table.</p>

Table 11.14 3D Source Tools (continued)

Item	Description
Measured Sources	<p>Quantification (DLIT) – For uncalibrated sources, the total flux measured for the sources selected using the Measure Source tool . For calibrated sources, this unit will be in [cell] units. For details on using this tool, see page 207.</p> <p>Quantification (FLIT) – For uncalibrated sources, the fluorescence yield measured for the voxels selected using the Measure Source tool . Fluorescence yield is expressed in units of [pmol M⁻¹cm⁻¹] here for uncalibrated sources. For calibrated sources, this unit will be in either [cells] or [pmol]. For details using this tool, see page 207.</p> <p>Volume – Volume of the selected source (mm³).</p> <p>Center of Mass (DLIT or FLIT) – The weighted average x, y, and z-coordinates of the selected voxels, where the weights are the flux of each highlighted voxel.</p> <p>Host Organ – The location of the selected source can be referenced to an organ atlas, and the organ from the atlas that is closest to the source will be reported here. This information is available if you select and register an organ atlas with the reconstruction. For more details, see page 227.</p>
Export Voxels	Enables you to export the voxel measurements in their x-, y-, and z-coordinates and source intensities (.csv file).
 Center of mass	Click to compute the center of mass for the source selected with the Measure Source tool  . For details using this tool, see page 207 .

11.14 3D Optical Registration Tools

Mouse anatomy reference atlases are available for registration with 3D reconstructions. A mouse anatomy reference atlas is used when volumetric data from another imaging modality is not available. A reference atlas provides guidance for the bioluminescent or fluorescent source anatomical location.




Use the Registration tools to:

- Display organs in the surface ([page 224](#))
- Manually adjust the location or scale of organs in the surface ([page 225](#))
- Check the organ fit ([page 226](#))
- Import an organ atlas ([page 227](#))


You can check the organ fit in the 3D View window ([page 226](#))

Figure 11.42 3D Registration Tools and Surface with Fitted Organs (“Skin” not displayed)

Table 11.15 3D Registration Tools (continued)

Item	Description
	Shading styles for the organs (see “ Display Organs ”).
Opacity	Adjusts the opacity of the organ display.
Organ Atlas	Choose a type of organ atlas.
	Click to select all organs in the database and display them on the surface.
	Click to clear the selected organs and remove all organ diagrams from the surface.

Displaying Organs With a Reconstruction

1. Load reconstruction results and confirm that the surface is in the perspective view (click the  toolbar button in the 3D View window or press the **R** key).
2. In the 3D registration tools, choose the Display Organs option and select an organ atlas. The organs in the selected atlas appear on the surface.
3. To fit the organs to the surface, click a registration tool:





Rigid registration: Performs linear transformation, but keeps the shape of the atlas surface.






Full registration: Performs linear transformation and volume deformation.




NOTE: For an optimum fit when there is a large difference between the orientation or size of the atlas organs and surface, first use the transformation tool to manually register the surface and atlas organs, then click a registration tool to automatically fit the organs. (See [Manually Adjusting the Scale or Location of Organs on page 225](#) for more details.)

4. If necessary, adjust the opacity of the organs using the slider or enter a number in the box. The organs are easier to view if you uncheck Skin in the Organs list.
5. To clear all organs from the surface, click the **Deselect All** button . To hide a particular organ, remove the check mark next to the organ name.
6. To display a specific organ(s), choose the organ name. To display all organs on the surface, click the **Select All** button .




NOTE: After fitting organs to the surface using the  or  tool, if necessary, you can click Reset button  to restore the default fit.

Manually Adjusting the Scale or Location of Organs

1. Load reconstruction results and confirm that the surface is in the perspective view (click the  toolbar button in the 3D View window or press the **R** key).



NOTE: It may be helpful to view the 3D image from different perspectives to check the organ position and size. To turn and rotate the 3D image, press and hold the left mouse key, then drag the mouse when the hand  appears.



2. In the 3D registration tools, choose the Display Organs option and select an organ atlas. The organs in the selected atlas appear on the surface. In [Figure 11.43](#), only “Skin” is selected.
3. Click the **Transform tool** button . The transform tool appears. [Figure 11.44](#) explains the tool functions.

Figure 11.43 Displaying the Transform Tool





Pink = skin surface
from the organ atlas

Gray = surface of the
imaged mouse

Use the transform tools to match
atlas organs (pink skin in this
example) as closely as possible to
the mouse surface (gray)



Figure 11.44 Transform Tools

- Press the Tab key to switch between the transform tools.
The position of the organ(s) is updated in the slice windowpanes (coronal, sagittal, and transaxial views) after each adjustment.
- Turn off the transform tool when you are done adjusting the position of the organ(s) (click the  button).

Checking the Organ Fit


- Check the fit in the coronal, sagittal, and transaxial windowpanes.
- Click the **Change view** toolbar button  .
The Top view is displayed.

Figure 11.45 Skin (pink) Fitted to Surface (gray)



3. Press the **V** key or the  button to display alternative views of the surface.

Figure 11.46 Alternative Surface Views

In this example, “skin” is selected from the organ atlas (pink surface). The mouse surface is gray.




Importing an Organ Atlas

An organ atlas (.iv, .dxf, or .stl, one organ per file) consisting of segmented organ surfaces derived from an MRI or CT scan can be imported into the Living Image software for registration with the animal surfaces derived from IVIS data. Organ files must be segmented from MRI or CT 3D volumetric data in third party medical imaging analysis software.



NOTE: The imported atlas must include a surface (skin) file which delineates the animal surface. The file name must include the word “skin”, for example *rat skin.iv*.

1. Load a DLIT or FLIT image sequence that is associated with the mouse comprising the organ files in *.iv, *.dxf or *.stl format.
2. Select **File** → **Import** → **Organ Atlas** on the menu bar.
3. In the dialog box that appears, click **Add Organ Files** (Figure 11.47).

Figure 11.47 Importing an Organ Atlas

4. In the next dialog box that appears, select all of the files (.iv, .dxf, .stl) that you want to include in the atlas (one file per organ) and click **Open**.
5. In the Select Skin Mesh drop-down list, select the skin organ file, which must include 'skin' in the file name.
6. Click **Generate Mesh Coefficients**.
7. Enter a name for the atlas and click **Save Organ Atlas**.

The organ atlas (.atlas) is created and is added to the Organ Atlas drop-down list (in the 3D Optical Tools, Registration tab).

11.15 3D Animation

Living Image software can create an animation from a sequence of 3D views (*key frames*). For example, an animation can depict a rotating 3D scene (Figure 11.48). The animation (series of key frames) can be recorded to a movie file (.mov, .mp4, or .avi).

Use the animation tools to:

- View a preset animation (generated from a factory-loaded animation setup) (page 231)
- Create a custom animation (created from your custom animation setup) (page 234)
- Save an animation setup (page 233)
- Record an animation to a movie file (page 233)
- Edit an animation setup (page 233)

Figure 11.48 Individual 3D Views (key frames) in the Preset Animation “Spin CW on Y-Axis”









Table 11.16 3D Animation Tools

Item	Description
Time Scale%	The time stamp of a key frame in the animation on a time scale of 0-100. For example, if the animation is 10 sec long and includes five key frames: Key frame 1: Time stamp= 0; first frame of the animation. Key frame 2: Time stamp = 25%; frame occurs 2.5 seconds after the start of animation. Key frame 3: Time stamp = 50%; frame occurs 5.0 seconds after the start of animation. Key frame 4: Time stamp = 75%; frame occurs 7.5 seconds after the start of animation. Key frame 5: Time stamp = 100%; last frame of the animation.
Presets	A drop-down list of predefined animation setups.
Key frame	A 3D view. The software interpolates the key frames to create intermediate frames in real time, then generates an animated sequence from all of the frames. Each successive key frame in a sequence should differ slightly from the preceding one, so that motion is smoothly depicted when the frames are shown at a proper frame rate (frames/second). The Living Image software provides preset key frames or you can specify the 3D views for the key frames.
Preset Key Frame Factor	Determines how many key frames are used to generate one revolution in a spinning animation (No. of frames = $(4 \times \text{Key Frame Factor}) + 1$). Increasing the key frame factor reduces the time period between key frames and creates the appearance of finer movement. Decreasing the key frame factor increases the time period between key frames and creates the appearance of coarser movement.
FPS	Frames displayed per second in the animation sequence.
	Creates a new key frame from the current 3D view.
	Updates the selected key frame to the current 3D view.
	Deletes a selected or all key frames from the key frame box.
	Moves a selected key frame up in the key frame box.
	Moves the selected key frame down in the key frame box.
Total Duration	The total time of the animation sequence.
Play	Click to view the animation sequence defined by the current key frames and animation parameters.
Record	Displays a dialog box that enables you to save the current animation to a movie (.mov, .mp4, or .avi, .mpg).
Animation Setup	
Load	Displays a dialog box that enables you to open an animation setup (.xml).
Save	Displays a dialog box that enables you to save the current key frames and animation parameters to an animation setup (.xkf).

Viewing a Preset Animation

Preset animations are factory-loaded animation setups. They include predefined key frames which are used to generate the animation.

To view a preset animation:


1. Open an image sequence and load 3D reconstruction results.
2. Select properties to display in the 3D View window (for example, organs, sources, surface, or photon density maps).
3. Select **View** → **3D Animation** on the menu bar.
4. In the 3D Animation tools that appear:
 - a. Clear the key frame box if necessary (click the **X** button and select **Delete All**).
 - b. Make a selection from the Presets drop-down list. See [Table 11.16, page 230](#) for a description of the preset animations.

After a preset animation is selected, a list of the key frames appears.



NOTE: You can view multiple animations sequentially. For example, if you select Spin CW on X-Axis and Spin CW on Y-axis from the Presets drop-down list, the animation shows the 3D reconstruction spinning clockwise on the x-axis, then spinning clockwise on the y-axis.

5. Click **Play** to view the animation.





Creating a Custom Animation

To create an animation, specify a custom animation setup or edit an existing setup.

1. Open an image sequence and load 3D reconstruction results.
2. Select properties to display in the 3D View window (for example, organs, sources, surface, or photon density maps).
3. Select **View** → **3D Animation** on the menu bar.

The 3D Animation tools appear (Figure 11.50).

4. Clear the key frame box if necessary (click the  button and select **Delete All**).








5. To capture the first key frame, click the  button.
The first key frame is added to the key frame box.
6. Adjust the position of the reconstruction in the 3D View using an image tool (for example,  , or ). For more details on the image tools, see [page 196](#).
7. Click the  button.
The second key frame is added to the key frame box.

Figure 11.51 Example Key Frames for a Custom Animation



8. Repeat step 6 to step 7 until all of the key frames are captured. For details on how to edit the key frame sequence, see [page 233](#).

Click a key frame to display the associated 3D view and the time stamp (position in the time scale (0-100) at which the frame occurs in the animated sequence).

9. Confirm the defaults for FPS (frames per second) and Total Duration (length of animation) or enter new values.

$\text{FPS} \times \text{Total Duration} = \text{No. of frames generated to create the animation}$. The number of generated frames should be \geq to the number of key frames. Otherwise, the frames may not be properly animated.

10. To view the animation, click **Play**. To stop the animation, click **Stop**.

An animation setup (series of key frames) can be saved (.xkf) or recorded to a movie (.mov, mp4, .avi, mpg).

Managing Animation Setups

To save an animation setup:

1. Click **Save**.
2. Select a directory and enter a file name (.xkf) in the dialog box that appears.


To record the animation to a movie:






1. Click **Record**.
2. Choose a directory, enter a file name (.mov, mp4, .avi), and click **Save** in the dialog box that appears.

To edit an animation setup:


1. Open an image sequence and load a reconstruction.
2. Open an animation setup:
 - To select a predefined setup, make a selection from the Preset drop-down list.
 - To select a saved user-defined setup:
 - a. Click **Load**.
 - b. Select an animation setup (.xkf) in the dialog box that appears.

Figure 11.52 List of Key Frames in the Selected Animation




3. To add a key frame:
- a. Adjust the position of the reconstruction in the 3D view using an image tool (for example,  , or ). For more details on the image tools, see [page 196](#).
 - b. Click the  button.
 - c. To reorder a key frame in the sequence, select the key frame and click the  or  arrow.


To update a key frame:

- a. Select the key frame and adjust the 3D view.
- b. Click the  button.

To delete a key frame:

- a. Select the key frame that you want to remove.
- b. Click the  button and select **Delete Current**.

11.16 DLIT/FLIT Troubleshooting

Issue	Solution
No sources in solution	This can occur in DLIT or FLIT if the surface is not correct. For example, if a surface is imported into the 3D View from another source other than a Surface Topography analysis.
Surface has spikes	<p>The most common source of spiky surfaces are folds in the animal skin or fur, which corrupt the desired smooth lines projected on the animal from the laser galvanometer.</p> <ul style="list-style-type: none">■ Choose the 'Fur Mouse' option for 'Subject'.■ Smoothing the surface by using the 'Smooth' feature in the Surface Topography tools can help improve the surface. <div></div>
Bad Photon Density or NTF Efficiency fit	The optical properties or source spectrum may have been incorrectly chosen. For example, 'Mouse Tissue' optical property is appropriate for mice, but 'XPM-2/XFM-2' is only appropriate for the mouse phantom.

12 Quantification Database

Preparing and Imaging the Samples

Creating a Quantification Database on page 236

Managing Quantification Results on page 240

It is possible to determine the number of cells in a DLIT source or the number of dye molecules or cells in a FLIT source if a quantification database is available. The database is derived from an analysis of images of known serial dilutions of luminescent cells or fluorescent cells or dye molecules.

12.1 Preparing and Imaging the Samples


1. Prepare a well plate (4 × 6, 6 × 4, 8 × 12, or 12 × 8 well format) that contains a dilution series of luminescent cells or fluorescent dye at four or more concentrations.
2. Include at least four background wells that contain diluent only.
3. Place the well plate on the IVIS stage, positioning it so that it is centered and square in the field of view.



NOTE: All of the wells must be within view in the image. For wells containing fluorophores, FOV D is recommended to reduce shadows from well walls and ensure more uniform excitation of the wells.

4. Acquire the images:
 - Bioluminescent samples – Acquire one 'Open' filter image of the well plate.
 - Fluorescent samples – Acquire reflectance-illumination Filter Scan images using the appropriate excitation and emission bandpass filters.

The well plate in Figure 14.1 contains a dilution series of a sample at four concentrations. The image sequence is a filter scan set of images with the excitation filter centered at 465 nm for all the images, and emission filter images centered at 520 nm, 540 nm, 560 nm, and 580 nm.




12.2 Creating a Quantification Database


1. Load the well plate image sequence.
2. Select **Tools** → **Well Plate Quantification for “<name>_SEQ”** on the menu bar.

The Well Plate Quantification window appears.

3. For fluorescent samples, choose the Dye molecules or Cells option.



4. Select the well plate dimensions from the Well Plate Type drop-down list.
The first image in the sequence opens and a grid ROI appears on the image.






5. Adjust the grid ROI to closely fit the plate wells.
6. In the well plate table, select the table cells for the samples, and click **Set** (Figure 12.4). Clicking a row or column header selects the entire row or column.
7. To remove the “sample” designations from table cells, select the table cells and click the  button.
8. To apply a color to table cells:
 - a. Select the table cells and click the  button. Alternatively, right-click the selected table cells and choose Background Color on the shortcut menu.
 - b. Choose a color from the color palette that appears.

Figure 12.4 Select Sample Wells and Enter Cell Numbers or Molecule Concentrations



Well Plate Quantification Window



For Sequence: EL20090414101005_SEQ Image: EL20090414101005_001

Well Plate Type:  Well Plate Type

Fluorophore Type: ☒ Dye molecules ☐ Cells

Measurement

Sample Wells: 3D : 3A  

☐ Background Wells:  


Well Volume [μ L]: 1

☐ Apply to Sequence

Well Plate Quantification Plots Results

Set position and enter dilution values in nanoMolar units

	D	C	B	A
1				
2				
3	500000	250000	125000	62500
4				
5				
6				

For Image: EL20090414101005_001 




9. Enter the concentration values in the table cells in nanomolar units, if calibrating fluorescent dyes. Enter the cell values in dimensionless units if calibrating cells.
10. To delete a concentration or cell value, select the table cell and press the Delete key. Alternatively, right-click a selected value to view a shortcut menu of edit commands (for example, cut, copy, paste).
11. If calibrating fluorescent molecules, enter the fluid volume (microliters) for the highlighted wells. The highlighted well volumes must be equal.
12. Choose the Apply to Sequence option.
13. Choose the Background Wells option.
14. In the well plate table, select the background wells and click **Set**.
Clicking a row or column header selects the entire row or column. To remove the “background” well designations, click the  button.

Figure 12.5 Set Background Wells




15. Click **Quantify**.
The results are displayed



16. Check the linear fit of the data for each image in the quantification plot.

A good fit to the straight line gives confidence to the results values. Large deviations of individual points from a straight line could indicate possible issues with the dilution series or errors when entering sample dilution values.

17. To export the quantification plot values:

- Click the  button.
- In the dialog box that appears, select a folder for the file (.csv) and click **Save**.


18. To copy the quantification plot values to the system clipboard, click the  button.

Table 12.1 Quantification Results

Item	Description
Fluorescence	
Excitation (nm) Emission (nm)	The excitation and emission filter wavelengths for the image. 'Excitation' and 'Emission' filters will be specified for fluorescent images, and the 'Open' filter for 'Emission' will be specified for bioluminescent images.
Extinction Coeff	A measure of excitation photon absorption interaction with the well plate samples based on a base-10 logarithmic derivation. The quantum efficiency factor of the conversion of the absorbed photon to the emission wavelength is also included.
Cross Section	A measure of excitation photon absorption interaction with the well plate samples based on a natural logarithmic derivation. The quantum efficiency factor of the conversion of the absorbed photon to the emission wavelength is also included.
Bioluminescence	
Total Flux/cell	A measure of total flux (photon/sec) emitted from a single cell. This number can be used to estimate the number of cells from the total flux in the 3D quantification.

12.3 Managing Quantification Results

Quantification results can be saved with the image sequence and as a calibration database that is available in the DLIT or FLIT 3D reconstruction tools (Properties tab). If a calibration database is selected when defining the properties for performing 3D reconstruction, the 3D reconstruction results will be displayed in calibrated units for cell numbers or molecule quantities in picomole units.







Table 12.2 Managing Quantification Results

Item	Description
 <p>Sequence Name : PC3M_GFP Delete Load Overwrite</p>	<p>Delete - Removes the active quantification results from the image sequence.</p> <p>Load - Opens quantification results from the sequence path.</p> <p>Save - Saves the quantification results with the selected image sequence.</p> <p>Overwrite - Saves the results with the selected image sequence and overwrites previous results.</p>
 <p>Database Name : WPQUANT_1 Delete Load Overwrite</p>	<p>Delete - Deletes the database from the system.</p> <p>Load - Opens quantification results from the system path.</p> <p>Save - Saves the quantification results to a system database that is available for DLIT or FLIT reconstruction.</p> <p>Overwrite - Saves the results to the selected database name and overwrites previous results.</p>

Exporting Quantification Results

Right-click the results table to view copy and export options.

- Copy – Copies the selected rows to the system clipboard
- Select All – Selects all rows in the results table
- Export Results – Opens a dialog box that enables you to export the selected results to a text file



13 3D Multi-Modality Tools

About the 3D Multi-Modality Tools

Classifying 3D Volumetric Data on page 243

Volume Display Options on page 247

Smoothing a Volume on page 252

Viewing and Rendering Slices on page 253

Volume Information and Results on page 256

Registering Optical and Volumetric Data on page 257

Volume Data Viewer on page 265

Viewing RAW Volumetric Data on page 266

13.1 About the 3D Multi-Modality Tools

Use the 3D Multi-Modality tools to:

- Classify volumetric data (3D image data).
- View slices.
- Refine the appearance of the volume (*volume processing*).
- Register optical and imported volumetric data (for example, CT, MRI, or PET data).

3D Multi-Modality Tool Requirements

The Living Image 3D Multi-Modality tools require a separate license. Additionally, the graphics processing unit (GPU) must meet the minimum specifications shown in [Table 13.1](#).

If the appropriate license is not installed or the GPU does not meet these specifications, the 3D Multi-Modality tools will not appear in the Tool Palette.

Table 13.1 Minimum Graphics Card Specifications

Specification	Description
OpenGL Version Requirement*	OpenGL 2.0 and above
OpenGL Extension Requirement*	GL-EXT-texture3D
Graphics Card Memory	Minimum: 256MB (Dedicated + Shared) Recommended: 1GB (Dedicated)
Consumer Graphics Cards (Desktop/ Mobile, Windows/Mac)	Supported: <ul style="list-style-type: none"> ■ NVIDIA® GeForce® 8 Series and above (8, 9, 100, 200, 300 and 400 series) ■ ATI Radeon™ HD 4000 Series and above (4000 and 5000 series) Recommended: <ul style="list-style-type: none"> ■ Desktop - NVIDIA GeForce GT 240 and above ■ Mobile - NVIDIA GeForce GT 230M and above
Workstation Graphics Cards (Desktop/ Mobile, Windows/Mac)	Supported: <ul style="list-style-type: none"> ■ NVIDIA® Quadro® NVS Series and Above (NVS and FX series) ■ ATI FireGL™ V5600 and Above (FireGL, FirePro and CrossFire series) Recommended: <ul style="list-style-type: none"> ■ Desktop - Quadro FX 1800 and above ■ Mobile - Quadro FX 880M and above


*If these specifications are not met, the 3D Multi-Modality tools do not appear in the Tool Palette.

13.2 Classifying 3D Volumetric Data

3D Multi-Modality tools provide a histogram-based method for classifying 3D volumetric data. The histogram represents the distribution of voxel intensities in the 3D volumetric data and their color-opacity values. The goal of classification is to set color and opacity values for different intensity ranges so that the color-opacity map shows the volume regions that you are interested in (opaque in the map) and hides unimportant regions (transparent in the map).

For example, [Figure 13.1](#) shows how the histogram tool designed a color-opacity map that shows both the skin and bone. The histogram tool enables you to easily re-design the color-opacity map to show only the skin or only bone.

3D Multi-Modality tools also enable you to classify the volumetric data by specifying color and opacity values for different intensity ranges so that you can view or hide certain parts of the data as needed. A color-opacity map can be saved.



Specifying a Color-Opacity Map

- After the surface and volume data are loaded, confirm that the Display Volume option is selected.

Histogram Display Options



- To change the color table for the color-opacity map, make a selection from the Color table-Opacity Map drop-down list. To apply the reverse color table, select the Reverse option.
- To view the histogram in a separate window, click the  button (Figure 13.2).
- Select units for histogram display (Hounsfield, counts, or absorption).
- If the histogram intensity range appears narrow or suppressed, choose the Logarithmic Histogram option.
- This will enhance histogram display by magnifying the smaller regions of interest in the histogram while keeping noise and air-related intensity peaks high. It helps bring out hidden regions of the histogram for easier identification of interesting intensity ranges.

Figure 13.2 3D Multi-Modality Tools



Managing Control Points


Edit the 3D volumetric data color-opacity map using "control points" (Figure 13.3). During volume rendering, the color-opacity map is used to map color and opacity to the corresponding intensity value as well as interpolate color and opacity for all data between adjacent control points.

1. Place a control point on the histogram by clicking anywhere on the histogram between the ● point (represents the lowest intensity in the volume) and ○ point (represents the highest intensity in the volume).
2. Drag any control point up or down to set the opacity level that is associated with the intensity value represented by the point. Drag a user-added control point left or right to change the intensity associated with the opacity specified by the point.

When you add, delete, or modify a control point, the color-opacity map and the rendering of the volume data are updated in real-time.



NOTE: The minimum and maximum intensity levels associated with the ● and ○ control points cannot be changed. The opacity level associated with these points can be changed.



3. To select a color for particular data, double-click a control point. In the color palette that appears, choose a color and click **OK**. The software interpolates the color range between adjacent control points.
4. To delete a control point, right-click the point. To delete all control points, click the button.



NOTE: The ● and ○ control points cannot be deleted from the histogram.


Saving a Color-Opacity Map

A color-opacity map can be saved and applied to any volumetric data set.

1. Click the Save button (Figure 13.4).
2. In the dialog box that appears, select a folder for the file (.tfn) and enter a file name.
3. Click **Save**.

Loading a Color-Opacity Map

1. Click the Open button (Figure 13.4).
2. In the dialog box that appears, navigate to the map file (.tfn), and click **Open**.



13.3 Volume Display Options

Adjusting Image Quality

By default, the color-opacity map displays the volumetric data at original (1×) resolution. This means, for example, if the volume comprises 512 slices, then all of the 512 slices are displayed. You can increase or decrease the resolution of the data display from 0.5× to 3.0× resolution (see [Table 13.2](#) for examples).

If the resolution is increased, the software interpolates the data and adds slices to the volume. If the processing performance is impacted at the original resolution, you may want to reduce the resolution to improve performance. Reducing the resolution down-samples the data and fewer slices are displayed.

To adjust the image resolution:

1. Move the “Level of Detail Slider” to the left or right ([Figure 13.5](#)).
The color-opacity map is updated.
2. To return the resolution to 1×, click the Reset button





Table 13.2 Example Volume With 512 Slices at 1x Resolution

Volume Resolution	No. of Slices Displayed
0.5X	256
1X (original resolution)	512
1.5X	768
2X	1024
2.5X	1280
3X	1536

Adjusting Volume Opacity

Adjust the volume opacity using the slider in the 3D Multi-Modality tools.



Maximum Intensity Projection

A maximum intensity projection (MIP) projects all maximum intensity voxels in the view along the viewing direction into the viewing plane (Figure 13.7). Living Image automatically extracts a 2D image from the MIP. The extracted image is similar in appearance to an X-ray image and is available in the image window. For example, Figure 13.8 shows an overlay of a luminescent image on an extracted 2D "x-ray" image.




NOTE: If you change the volume opacity (see Figure 13.6 on page 248) and want to extract new 2D images from the MIP, click the  button (Figure 13.7).

Figure 13.7 Volume Display – Maximum Intensity Projection





Figure 13.8 Image Window – Luminescent Image on X-ray Image

Gradient Illumination

Gradient illumination is based on the idea that light is reflected at boundaries between different voxel intensities, but is not affected when passing through homogeneous regions. Choosing this option illuminates the voxels at boundaries more than voxels within a homogeneous region. The boundaries are based on the gradient magnitude between heterogeneous regions or the change in intensities between neighboring voxels in heterogeneous regions. Using this option enhances the variation in tissue properties and may be helpful for visualizing the boundaries of different tissues.


Two options are available for Gradient Illumination:

- **Quality** – This option will be automatically selected if your system has an appropriate graphics card. If this option is selected, the volume is displayed with more detail (Figure 13.9).



NOTE: If the system graphics card does not meet the recommended specifications (Table 13.1 on page 243), choosing the Quality option causes slow performance of actions such as rotating the volume.

- **Denoise** – Filters out image noise in the volume rendering. The raw data are not modified by this filter.

Figure 13.9 Volume Display – Gradient Illumination

Modifying Volume Resolution

Changing the pixel or slice spacing modifies the volume resolution. Increasing the pixel or slice spacing reduces resolution, while reducing either increases resolution.



1. In the Volume tab, click the Edit Space and Orientation button .
2. In the dialog box that appears (Figure 13.10), edit the pixel or slice spacing.

Figure 13.10 Volume Information Dialog Box

13.4 Smoothing a Volume

Smoothing a volume reduces noise in a CT, MRI, or PET image such as excessive variation in voxel grayscale values.

To apply smoothing:



1. Load the volumetric data.
2. Choose the type of smoothing and group size in the Process tab of the 3D Multi-Modality tools (Figure 13.11).
 - Mean – Applies the average grayscale value of a group of voxels (for example, a 3x3 group) to the central voxel of the group.
 - Median – Applies the median grayscale value of a group of voxels to the central voxel of the group.
 - Gaussian – Applies the weighted mean to the central voxel of the group. The weight distribution is similar to a normalized Gaussian shape with the highest value at the central voxel of the group.
3. Click the  button to remove the smoothing.

Figure 13.11 3D Multi-Modality Tools – Process Tab



13.5 Viewing and Rendering Slices

Viewing Slices

View volume slices by double-clicking the Coronal, Sagittal, or Transaxial windowpane (Figure 13.12).











Table 13.3 Volume Slice Viewer

Item	Description
Orientation	Select a slice orientation from the drop-down list.
Slice Spacing	The distance between each slice in the Volume Slice Viewer. Enter a smaller value to increase the number of slices in the viewer or a larger value to decrease the number of slices in the viewer.

Table 13.3 Volume Slice Viewer (continued)


Item	Description
Total Slices	The number of slices shown in the viewer.
Render Thick Slice	This option is used to create a sequence of 3D or maximum intensity projection (MIP) renderings from the image stack. When this option is selected, "Slice Spacing" changes to "Slice Thickness". Increasing the slice thickness causes more slices to be extracted from the volume before creating the rendering.
Gradient Illumination	Gradient Illumination is based on the idea that light is reflected at boundaries between different voxel intensities, but is not affected when passing through homogeneous regions. Choosing this option illuminates the voxels at boundaries more than voxels within a homogeneous region. The boundaries are based on the gradient magnitude between heterogeneous regions or the change in intensities between neighboring voxels in heterogeneous regions. Using this option enhances the variation in tissue properties and may be helpful for visualizing the boundaries of different tissues.
Maximum Intensity Projection (MIP)	Projects all maximum intensity voxels in the view along the viewing direction into the viewing plane.
 <p>Min: -15.00 mm Max: 15.00 mm</p> <p>Min – The slice coordinate of the first slice being viewed. Zero is defined as the center plane of the image. Max – The slice coordinate of the last slice being viewed. Specify the position range to include in the viewer using the Min and Max sliders or enter values.</p> <p>-7.50 Slice position</p>	
	Click to show the single view of the active slice in the multi view. Alternatively, double-click a slice in the multi view to show the single view.
	Click to show the multi view.
	If the single view has been magnified, click this button to zoom out incrementally.
	Magnifies the single view.
	Resets the single view to the default magnification.
	Click to export the slice view as a graphic file (for example, .bmp)

Rendering Slices

The Slice tab in the 3D MM tools contains rendering and viewing options for slices.




Table 13.4 3D Multi-Modality Tools for Rendering Slices

Item	Description
Slice Color Table	Choose this option to apply the color table selected from the Color Table drop-down list.
Volume Color Table	Choose this option to apply the volume color table of the volume color-opacity map that was selected in the Volume tab.
<div><div>Color Table</div><div>Red</div><div>Reverse</div></div>	Color table options. Choose the Reverse option to apply the inverse color table.
Opacity	Move the slider to adjust the color opacity.
Color Scale	Min – Sets the intensity level associated with the lowest color scale value. Max – Sets the intensity level associated with the maximum color scale value.
<div><div>Slice Viewer</div><div>Orientation: Coronal</div><div></div></div>	<p>The Tool Palette provides an alternative way to view slices and access the Slice Viewer.</p> <p>Choose a slice orientation from the drop-down list. Use the slider to move through the slices. Double-click the slice view or click the  button to open the Slice Viewer. The selected slice is highlighted in the Slice Viewer.</p>

13.6 Volume Information and Results

The Results tab displays information about the loaded data taken from the DICOM file header (Figure 13.14).



Saving the registered and classified data provides a convenient way to share data. The software saves the following:

- Level of detail setting
- Color tables for the opacity map and slices
- Histogram tool control settings and the resulting color-opacity map
- Multi-modal registration settings
- Crop settings

Managing Results

Saving Registered Results

1. In the Results tab, confirm the default name in the Name drop-down list or enter a name.
2. Click **Save**.

The registered 3D volumetric data, along with the color-opacity settings, appear in the 3D View window.



NOTE: The results are saved in XML format in the XML folder within the sequence folder. The results can only be accessed from the same optical data set.

Loading Results

1. Select the results from the Name drop-down list.
2. Click **Load**.

Deleting Results

1. Select the results from the Name drop-down list.
2. Click **Delete**.
3. Click **Yes** in the confirmation message that appears.

13.7 Registering Optical and Volumetric Data


Registering *multi-modal* data (optical and volumetric data) provides an anatomical context for interpreting biological (functional) information. Two registration methods are available:

- Automatic fiducial registration – For experiments in which the optical data are acquired on the IVIS Spectrum and the CT data are acquired on the Quantum FX μ CT instrument. The subject must be contained in the Mouse Imaging Shuttle during both optical and CT imaging, and the CT data must be exported to DICOM format. See [page 261](#) for more details.
- Manual registration – Use the 3D Multi-Modality tools to register a 3D surface reconstruction with 3D volumetric data acquired on a third party instrument. See [page 263](#) for more details.

[Figure 13.15](#) shows an overview of the steps to register these types of multi-modal data. After registration, classify the 3D volumetric data to help identify and separate objects (see [page 243](#)).

Figure 13.15 Registering Multi-Modal Data

1. Load the optical data:
 - Bioluminescence or fluorescence image sequence and structured light surface
 - 3D source reconstruction (DLIT or FLIT results) ([page 192](#) or [page 197](#))
2. Load 3D volumetric data (CT or MRI) ([page 259](#)).
3. Register the 3D source reconstruction and the 3D volumetric data by performing either:
 - Automatic fiducial registration–Available for data acquired on the Quantum FX µCT instrument using the Mouse Imaging Shuttle ([page 257](#))
 - or
 - Manual registration–Match animal surface representations using the Manual Registration tool ([page 263](#))
4. Classify the 3D volumetric data to help identify and separate objects ([page 243](#)). Save the color-opacity map (optional).
5. Save the registered 3D multi-modality results ([page 256](#)).



Loading Data for Registration

1. Load a DLIT or FLIT image sequence and the 3D reconstruction results.



NOTE: The 3D Multi-Modality tools appear in the Tool Palette after you load optical image data. If the 3D Multi-Modality tools do not appear in the Tool Palette, confirm that the 3D Multi-Modality Tools license is installed and that the workstation graphics card meets the specifications in [Table 13.1 on page 243](#).


2. Select the DICOM or TIFF volumetric data

- a. Select **File** → **Browse 3D Volumetric Data** on the menu bar.
- b. Select a data folder in the Browse For Folder box that appears and click **OK**.

The Living Image 3D Volumetric Browser appears ([Figure 13.16](#)).



NOTE: Only DICOM or TIFF data can be added to the 3D Volumetric browser. For details on loading other data types (.raw or .vox files) see [page 265](#).



NOTE: The next time you start the Living Image software and open the Browse For Folder box, the software automatically returns to the last folder visited.

The 3D Volumetric Browser automatically previews a playback of the data along with other information about the data ([Figure 13.17](#)).

 - DICOM file

 - TIFF file

3. Load the volumetric data with the optical data:

- a. Confirm that the “Load in a new window” option is not selected. (If this option is selected, the volumetric data are loaded in a new window.)
- b. Double-click the data row in the browser. Alternatively, select the data row and click **Load**.

The 3D volumetric data appears in the 3D View window of the optical data ([Figure 13.18](#)). The software converts loaded volumetric data into an 8-bit representation to reduce memory overhead and for easier color mapping. The 3D Multi-Modality tools provide an 8-bit color-opacity map for volume visualization which maps each voxel to an RGB color, or a color and opacity value.

A histogram of voxel intensities appears in the Multi-Modality tools and the software sets a default air/noise boundary.





Table 13.5 3D Volumetric Data Browser

Item	Description
Add to List	If this option is chosen, the data selected in the Browse for Folder box is added to the 3D Volumetric Data Browser. If this option is not chosen, the data selected in the Browse for Folder box replaces the contents of the 3D Volumetric Data Browser, except for loaded data.
Browse	Opens the Browse For Folder box.
Load in a new window	If this option is chosen, multiple data sets can be loaded, each in a separate window. If this option is not chosen, only one data set can loaded at a time.
Load	Click to open the data selected in the 3D Volumetric Data Browser.



Registering Multi-Modal Data


Automatic Fiducial Registration

About the Mouse Imaging Shuttle

The Mouse Imaging Shuttle (PN 127744) contains the subject during imaging and enables the subject to be transferred between an IVIS® Imaging System and the Quantum FX µCT instrument without disrupting the subject's position.

The Mouse Imaging Shuttle must be correctly docked to the docking station in the IVIS Imaging System and the Quantum FX µCT instrument. The docking station in the Quantum FX µCT system is marked with a triangle-shaped fiducial pattern under the plane where the Mouse Imaging Shuttle docks. Automatic fiducial registration is available if both sides of the triangle fiducial pattern are included in the CT images. For more details on using the Mouse Imaging Shuttle, see the *Mouse Imaging Shuttle Instructions* (PN 127820_RevA).

To perform automatic fiducial registration:

1. Load the data that you want to register (see [page 259](#)).
2. Click the Fiducial Registration button .


The multi-modal data are automatically registered and cropped ([Figure 13.19](#)).
3. To undo the registration, click the Reset Registration button .
4. To save the registration information:
 - a. Confirm the default name or enter a name for the results in the Results tab.
 - b. Click **Save**.

Figure 13.19 Registered 3D Optical and 3D Volumetric Data





Figure 13.20 3D Multi-Modality Tools – Results




NOTE: Registration information is saved with the results for the volumetric data and is specific for a particular optical data set.

Manual Registration

To manually register data, use the 3D Multi-Modality tools to translate, scale, or rotate the 3D volumetric surface so that features common to both surfaces are matched and aligned in the x, y, and z planes. Examine the matched surfaces in the 3D slice views to help you fine tune the registration.


Figure 13.21 Example Surfaces Before and After Registration



To manually register data:

1. Load the data that you want to register (see [page 259](#) for more details).
The software determines a default air/noise boundary for the 3D volumetric data ([Figure 13.22](#)).
2. If you need to remove noise from the 3D volumetric data, move the air/noise boundary to the right in the histogram tool.

Figure 13.22 Adjusting Air/Noise Boundary Using Histogram Tool



3. If the volumetric data needs cropping (for example, to remove structures such as the stage from the CT view), follow step a to step c below. If cropping is not needed, proceed to step 4.

To crop the data:

- a. Click the crop tool button .

The crop tool appears and has six control points:





-  – Crops the data along the x-axis.
-  – Crops data along the y-axis.
-  – Crops data along the z-axis.

Figure 13.23 Crop Data Along X,Y, or Z-Axis







- b. Click and hold a control point while you move the crop plane. As you move the crop plane, the slice views are updated. Release the mouse button to crop the data.
- c. To reset the crop planes, click the  button. When finished cropping, press the Tab key to turn off the crop tool.
4. Click the Manual Registration button .
- The transformation tool appears (Figure 13.24). The tool has three modes that enable you to translate, scale, or rotate the 3D volumetric data (press the Tab key to change the tool mode). The slice views are automatically updated when you use the tool.

Figure 13.24 Manual Registration Tool – Transformation Modes

NOTE: Make sure that you click the transformation tool so that it is highlighted before you use it. Otherwise the dragging operation is applied to the optical data (structured light surface).

5. To return the 3D volumetric data to the default position and size, click the Reset Registration button .
6. Save the registration information (see [page 261](#)).



NOTE: Registration information is saved with the results for the volumetric data and is specific for a particular optical data set.

13.8 Volume Data Viewer

The Living Image software provides a viewer for volumetric data. The 3D Multi-Modality tools are not required to view DICOM or TIFF data.




1. Select **View** → **Volume Data Viewer** on the menu bar.
The Volume Data Viewer appears.
2. Select volume data by doing either of the following:
 - Drag the data file (DICOM, TIFF) from Windows Explorer to the Volume Data Viewer window
or
 - In the Volume Data Viewer, click the Open button , and in the dialog box that appears, select a DICOM or TIFF file, and click **Open**.
3. To clear the Volume Data Viewer, click the  button.


Figure 13.25 Drag Volume Data from Windows Explorer to Volume Data Viewer

13.9 Viewing RAW Volumetric Data

1. Drag a single RAW file (*.raw or *.vox) from Windows Explorer to the 3D Multi-Modality tools (Figure 13.26).




NOTE: Only single *.raw or *.vox files consisting of multiple slices of a 3D volume can be loaded into Living Image.

Figure 13.26 Opening RAW Volumetric Data

2. In the Volume Information dialog box that appears (Figure 13.27), enter the:
 - Data width, height, and the number of slices.
 - Slice row, column pixel size, and the slice spacing in millimeters.

Figure 13.27 Volume Information




Select a data type.

Enter the:

- Width, height, and number of slices
- Slice row, column pixel spacing, and the slice spacing in millimeters

3. If loading the data will cause low memory, you are prompted to down-sample the data (Figure 13.28) to improve memory and performance. Decrease the slice resolution by moving the Slice Resolution slider to the left until the Memory Requirement Status is “Good”.

Figure 13.28 Down-sample 3D Volumetric Data



Changing the Orientation of RAW Volumetric Data

Occasionally, RAW files (*.raw or *.vox) may be loaded with the orientation “flipped” or reversed along the x, y, or z-axis. As a result, the slice views (transaxial, coronal, sagittal) may be flipped or rotated so that the actual view that is displayed does not match the 3D View windowpane name (for example, the Sagittal windowpane does not display a sagittal slice), or the data appears flipped with respect to the surface derived from the IVIS® Spectrum.

In such cases, you can:

- Invert the data along the x, y, or z-axis
- Manually rotate the data using the Transformation tool (see page 265 for more details).

To invert the subject orientation:



1. Click the Edit Spacing & Orientation button .
2. In the dialog box that appears, choose a “Subject Orientation” option and click **OK**.

Figure 13.29 Volume Information


Appendix A IVIS Acquisition Control Panel

Control Panel

Manually Setting the Focus on page 272

A.1 Control Panel

The control panel provides the image acquisition functions (Figure A.1).




 **NOTE:** The options available in the IVIS acquisition control panel depend on the selected imaging mode, the imaging system, and the filter wheel or lens option that are installed.

Table A.1 IVIS acquisition control panel




Item	Description
 Luminescent	Choose this option to acquire a luminescent image.
 Fluorescent	Choose this option to acquire a fluorescent image. If the Fluorescent option is selected, the following options also appear in the control panel: Transillumination - Choose this option to acquire a fluorescent image using transillumination (excitation light located below the stage). Normalized - This option is selected by default when the Fluorescent and Transillumination options are chosen so that NTF Efficiency images can be produced.
 Photograph	Choose this option to automatically acquire a photograph. The illumination lights at the top of the imaging chamber are on during a photographic image so that the system can acquire a black and white photograph of the sample(s). Note: You can adjust the appearance of the photographic image using the Bright and Contrast controls (see Viewing Intensity Data on page 65).

Table A.1 IVIS acquisition control panel (continued)


Item	Description
 Structure	<p>Choose this option to take a structured light image (an image of parallel laser lines scanned across the subject) when you click Acquire. The structured light image is used to reconstruct the surface topography of the subject which is an input to the Diffuse Luminescence Imaging Tomography (DLIT™) algorithm that computes the 3D location and brightness of luminescent sources.</p> <p>When this option is chosen, the f/stop and exposure time are automatically set to defaults for the structured light image (f/8 and 0.2 sec, respectively). The spatial resolution of the computed surface depends on the line spacing of the structured light lines. The line spacing and binning are automatically set to the optimal values determined by the FOV (stage position) and are not user-modifiable.</p>
Overlay	If this option is chosen, the system automatically displays the overlay after acquisition is completed (for example, luminescent image on photograph).
Exposure time	<p>The length of time that the shutter is open during acquisition of an image. The luminescent or fluorescent signal level is directly proportional to the exposure time. The goal is to adjust the exposure time to produce a signal that is well above the noise (>600 counts recommended), but less than the CCD camera saturation of ~60,000 counts.</p> <p>Luminescent exposure time is measured in seconds or minutes. The minimum calibrated exposure time is 0.5 seconds. The exposure time for fluorescent images is limited to 60 seconds to prevent saturation of the CCD. There is no limit on the maximum exposure time for luminescent images; however, there is little benefit to exposure times greater than five minutes. The signal is linear with respect to exposure time over the range from 0.5 sec to 10 minutes. Integration times less than 0.5 seconds are not recommended due to the finite time required to open and close the lens shutter.</p>
Binning	<p>Controls the pixel size on the CCD camera. Increasing the binning increases the pixel size and the sensitivity, but reduces spatial resolution. Binning a luminescent image can significantly improve the signal-to-noise ratio. The loss of spatial resolution at high binning is often acceptable for <i>in vivo</i> images where light emission is diffuse. For more details on binning, see the reference article <i>Detection Sensitivity</i> (select Help → References on the menu bar).</p> <p>Recommended binning: 1-4 for imaging of cells or tissue sections, 4-8 for <i>in vivo</i> imaging of subjects, and 8-16 for <i>in vivo</i> imaging of subjects with very dim sources.</p>
F/stop	<p>Sets the size of the camera lens aperture. The aperture size controls the amount of light detected and the depth of field. A larger f/stop number corresponds to a smaller aperture size and results in lower sensitivity because less light is collected for the image. However, a smaller aperture usually results in better image sharpness and depth of field.</p> <p>A photographic image is taken with a small aperture (f/8 or f/16) to produce the sharpest image and a luminescent image is taken with a large aperture (f/1) to maximize sensitivity. For more details on f/stop, see the reference article <i>Detection Sensitivity</i> (select Help → References on the menu bar).</p>
Excitation Filter	<p>A drop-down list of fluorescence excitation filters. For fluorescent imaging, choose the appropriate filter for your application. For luminescent imaging, Block is selected by default. If you select Open, no filter is present. For systems equipped with spectral imaging capability, choose the appropriate emission filter for your application.</p> <p>Note: On some models with standard filter sets, the excitation filter selection automatically sets the emission filter.</p>
Emission Filter	<p>A drop-down list of fluorescence emission filters located in front of the CCD lens. The emission filter wheel is equipped with filters for fluorescence or spectral imaging applications. The number of filter positions (6 to 24) depends on the system. For luminescent imaging, the Open position (no filter) is automatically selected by default.</p>

Table A.1 IVIS acquisition control panel (continued)




Item	Description
Lamp Level	<p>Sets the illumination intensity level of the excitation lamp used in fluorescent imaging (Off, Low, High, and Inspect). The Low setting is approximately 18% of the High setting. Inspect turns on the illumination lamp so that you can manually inspect the excitation lamp.</p> <p>Note: Make sure that the filters of interest are selected in the filter drop-down lists before you select Inspect. The Inspect operation automatically positions the selected filters in the system before turning on the lamp. Subsequent changes to the filter popup menus will have no effect until another Inspect operation is performed.</p>
Lights	Turns on the lights located at the top of the imaging chamber.
Alignment Grid	Choose this option to illuminate an alignment grid on the stage when the imaging chamber door is opened. The alignment grid shows the sizes and positions of the possible fields of view. If subject alignment is not completed in two minutes, place a check mark next to Alignment Grid to turn on the grid.
Field of View	Sets the size of the stage area to be imaged by adjusting the position of the stage and lens. The FOV is the width of the square area (cm) to be imaged. A smaller FOV gives a higher sensitivity measurement, so it is best to set the FOV no larger than necessary to accommodate the subject or area of interest. The FOV also affects the depth of field (range in which the subject is in focus). A smaller FOV results in a narrower depth of field. Select the FOV by choosing a setting from the drop-down list. See Table A.2 for more details on the calibrated FOV positions.
Service	Moves the stage to a position for cleaning the imaging chamber below the stage.
Load	Moves the stage from the cleaning position back to the home position.
MIS	Choose this option if the subject will be contained in the Mouse Imaging Shuttle during image acquisition.
Subject height (cm)	<p>Sets the position of the focal plane of the lens/CCD system by adjusting the stage position. The subject height is the distance above the stage that you are interested in imaging. For example, to image a mouse leg joint, set the subject height to a few mm. To image the uppermost dorsal side of a mouse, set the subject height to the 1.5 - 2.0 cm. The default subject height is 1.5 cm.</p> <p>IMPORTANT! The IVIS® instrument has a protection system to prevent instrument damage, however always pay close attention to subject height. For example, it is possible for a large subject (10 cm ventral-dorsal height) to contact the top of the imaging chamber if you set the subject height = 0 and choose a small FOV.</p>
Focus	<p>Drop-down list of focusing methods available:</p> <p>Use subject height – Choose this option to set the focal plane at the specified subject height.</p> <p>Manual – Choose this option to open the Focus Image window so that you can manually adjust the stage position. For more details on manual focusing, see page 272.</p>
Batch Sequences	Choose this option if you want to specify multiple, separate image sequences for batch acquisition (multiple image sequences are automatically acquired, one after another, without user intervention). See page 45 for more details.
Temperature   	<p>The temperature box color indicates the temperature and status of the system:</p> <ul style="list-style-type: none"> White box – System not initialized. Red box – System initialized, but the CCD temperature is out of range. Green box – System is initialized and the CCD temperature is at or within acceptable range of the demand temperature and locked. The system is ready for imaging. <p>Click the temperature box to display the actual and demand temperature of the CCD and stage. See page 8 for more details.</p>
Acquire	Click to acquire an image using the settings and options selected in the control panel or to acquire an image sequence specified in the Sequential Setup table.

Table A.1 IVIS acquisition control panel (continued)

Item	Description
Sequence Setup	Click to display the sequence table so that you can specify and manage sequence acquisition parameters, or open sequence acquisition parameters (xsq). See page 47 for more details on manually setting up an image sequence.
Imaging Wizard	Click to start the Imaging Wizard
Sequence Setup	Click to open the sequence table.
Image Setup	Click to close the sequence table.
Initialize	Click to initialize the IVIS® Spectrum. See page 7 for more details on initializing the system.

Table A.2 Field of view (FOV) settings


FOV Setting	FOV (cm)
A	4
B	6.5
C	13
D	22.5

A.2 Manually Setting the Focus


The IVIS Imaging System automatically focuses the image based on subject height. If you do not want to use the automatic focus feature, you can manually set the focus.

1. In the control panel, choose **Manual Focus** in the Focus drop-down list.

The Manual Focus window appears.



2. To mark the center of the camera in the window, put a check mark next to Display CCD Center.
3. Select the size of the step increment that the stage moves: **Coarse**, **Normal**, or **Fine**.

4. Click **Up** or **Down** to move the stage and change the focus.
5. If necessary, select another F/stop setting from the drop-down list and adjust the light level using the  arrows.
6. Click **Update** to apply the settings.
The resulting focal plane (cm above the stage) is automatically entered in the Subject height box.
7. Click **OK** when the image is focused.

Appendix B Preferences

General Preferences

Options on page 276


Acquisition on page 277

Theme on page 278

Optical Properties on page 281

You can manage user IDs and specify defaults for some parameters that are associated with the user ID selected at the start of a new session.

After you log on, select **Edit** → **Preferences** on the menu bar to view the user-modifiable preferences.

 **NOTE:** Any changes made to the Preferences are implemented at the start of the next session. The Acquisition tab is only available in the Living Image software that controls the IVIS Imaging System.

B.1 General Preferences






Table B.1 General Preferences

Item	Description
Start Up Defaults	Dock Tool Palette - Choose this option to set the position of the Tool Palette in the application window. Choose left or right. Note: To undock the Tool Palette, click on the palette title bar and drag it a distance greater than its width.
Window Size	Specifies the dimensions of the main application window. Width, Height - Sets the dimensions of the image window. Restore Defaults - Click to apply the default settings.

Table B.1 General Preferences (continued)

Item	Description
Apply Individual Color Scale for Sequences	Choose this option to apply a separate color scale to each thumbnail of a sequence. If this option is not chosen, all of the thumbnails are displayed using the same color scale.
Show Transillumination Locations	Choose this option to display a cross hair at each transillumination location when you load transillumination data. When you mouse over a cross hair, a tool tip displays the transillumination coordinates. If this option is not chosen, you can choose the Transillumination Location option in the sequence view window to display the transillumination locations.
Show Advanced Options	If this option is selected, advanced features are available in the menu bar and Tool Palette, including: <ul style="list-style-type: none"> ■ Additional ROI functionality for Auto ROI parameters. ■ Additional export and import option for 3D surfaces and voxels. ■ Planar Spectral Imaging tools in the Tool Palette.
Show Activity Window on:	A drop-down list of options for when to display the activity log (Figure B.2).
Save Settings	<p>Save float-corrected image - Saves an image after all corrections are applied (read bias subtraction, flat field correction, cosmic correction).</p> <p>Color Selections - Applies the color settings of the active image data to subsequently opened image data.</p> <p>Folder Locations - Sets the default folder path to the current folder path setting. Click the Export button  in the image window to view the current folder path setting (Figure B.2).</p> <p>Window Size & Position - Applies the active image window size and position settings to subsequently opened image data.</p> <p>Most Recently Used Dataset History - Defines the number of recently opened data sets to remember and display when you select File → Recent Files → Menu.</p>
Display ROI Label As Measurement	Sets the type of measurement in counts, radiance (photons), or efficiency to show in the ROI label

Some of the general preferences specify how the main application window is organized. To undock the Tool Palette, click on the palette title bar and drag it a distance greater than its width. To dock the Tool Palette in the main window, drag the palette to the right or left side of the window and release.



B.2 Options




Table B.2 User Preferences

Item	Description
Edit label Choices	Opens a dialog box that enables you to edit the Living Image Universal label set.
Default Units	Choose counts or radiance (photons) for image display.

B.3 Acquisition






Table B.3 Auto Exposure Settings

Item	Description
Luminescent/Fluorescent Auto Exposure Preferences	
First Preference Second Preference Third Preference	During auto exposure, the software acquires a luminescent or fluorescent image so that the brightest pixel is approximately equal to the user-specified Target Count (Minimum). If the target minimum count cannot be closely approximated by adjusting the first preference (for example, exposure time), the software uses the first and second or first, second and third preferences to attempt to reach the target max count during image acquisition.
Target Count (Minimum)	A user-specified intensity.
Range Values Exp Time (sec) Binning F/Stop	The minimum and maximum values define the range of values for exposure time, F/Stop, or binning that the software can use to attempt to reach the target max count during image acquisition.
Restore Defaults	Click to apply default settings.

**Table B.4** Camera Settings

Item	Description
Default Image Exposure	Sets the default exposure settings that appear in the IVIS acquisition control panel.
Default Image Binning	Standard - Binning choices include Small, Medium and Large. These are predetermined, factory-loaded binning values that depend on the imaging system camera. Manual - Allows the user to choose a binning value (1, 2, 4, 8 or 16)
Auto Save	Specifies the folder where images are automatically saved. Click the  button to select a folder.
Restore Defaults	Click to apply the default settings.

B.4 Theme




Table B.5 Image View Preferences




Item	Description
Color Palette	Use these controls to select a color table for luminescent and fluorescent image data. Choose the Reverse option to reverse the min/max colors of the selected color table.
Use saved color palette while loading datasets	If this option is chosen, data are displayed using a user-specified color palette. For example, after you load data, specify a color table in the Image Adjust tools, and save the data. The user-specified color table is automatically applied whenever the data are loaded.
Background & Text Color	<p>Sets the color of the:</p> <ul style="list-style-type: none">■ Background in the image window (shown below)■ Text for the color bar <p>To change a color, click the  button that opens the color palette.</p> 
ROI Color	<p>Sets the colors for the ROI outline. To change a color, click the  button that opens the color palette.</p> <p>Luminescent - Color of the ROI outline on a luminescent image.</p> <p>Fluorescent - Color of the ROI outline on a fluorescent image.</p>
Restore Defaults	Click to apply the default settings.

Figure B.7 3D View Preferences






Table B.6 3D View Preferences

Item	Description
Color Theme	<p>Predefined color schemes available for the 3D View window shown here. Click the  button to restore the defaults for the selected color theme.</p> 
Background Color	<p>Settings that modify the appearance of the background in the 3D View window.</p> <p>Solid Color - Choose this option to apply a non-gradient background color to the 3D view in the image window.</p> <p>Gradient Color - Choose this option to apply a gradient background color to the 3D view in the image window. Top = the color at the top of the window; Bottom = the color at the bottom of the window.</p>
Surface & Text Color	Settings that modify the display of the surface and text in the 3D View window.
Color Palette	<p>Source voxels - Choose a color table for voxel display.</p> <p>Reverse - Choose this option to reverse the min/max colors of the selected color table.</p>
Restore Defaults	Click to apply the default settings.

B.5 Optical Properties

Figure B.8 Set Default Optical Properties Preferences (left) for the Properties Tab in the Planar Spectral Imaging, DLIT, or FLIT Tools




Table B.7 Preferences – Optical Properties

Item	Description
Tissue Properties	Choose a default tissue type that is most representative of the area of interest. This tissue type will be used if a Subject Type is not selected in the Imaging Wizard and saved during acquisition.
Source Spectrum	Choose the default luminescent source spectrum. This Source Spectrum will be used if a Subject Type is not selected in the Imaging Wizard and saved during acquisition for DLIT sequences.
Plot	<p>Tissue Properties - Choose this option to display a graph of the absorption coefficient (μ_a), effective attenuation coefficient (μ_{eff}), and reduced scattering coefficient (μ'_s or μ_{sp}).</p> <p>Source Spectrum - Choose this option to display the source spectrum for DLIT reconstructions.</p> <p>Bioluminescent Spectrum - Choose this option to display the spectrum of the bioluminescent source (available for DLIT reconstructions only).</p> <p>Fluorescent Spectrum - Choose this option to display the spectrum of the fluorescent source (available for FLIT reconstructions only).</p>
Restore Defaults	Click to restore the defaults in the Optical Properties tab.

Appendix C Menu Commands, Toolbars, and Shortcuts

Figure C.1 Living Image toolbar




Table C.1 Menu bar commands and toolbar buttons







Menu Bar Command	Toolbar Button	Description
File → Open		Displays the Open box so that you can select and open an image data file. Double-click a SequenceInfo.txt file or ClickInfo.txt file to open the image data file (see page 56).
File → Browse		Displays the Browse For Folder box so that you can select and an image data folder. The selected folder is displayed in the Living Image Browser.
File → Browse 3D Volumetric Data		Displays the Browse For Folder box so that you can select a volumetric data folder (for example, DICOM format, TIF data). The selected folder is displayed in the 3D Browser.
File → Save		Saves (overwrites) the AnalyzedClickInfo text file to update the analysis parameters, but the original image data files are not altered.
File → Save As		Displays the Browse For Folder box so that you can specify a folder in which to save the image data. The original data is not overwritten.
File → Import → 3D Surface		Opens a dialog box that enables you to import a surface. Note: This menu item is only available if “Show Advanced Options” is selected in the Preferences (see page 274).
File → Import → 3D Voxels		Opens a dialog box that enables you to import a source volume. Note: This menu item is only available if “Show Advanced Options” is selected in the Preferences (see page 274).
File → Import → Atlas		Opens a dialog box that enables you to import an organ atlas (.iv, .dxf, .stl).
File → Export → Image/Sequence as DICOM		Opens the Browse for Folder dialog box that enables you to export the active image data to DICOM format (.dcm).
File → Export → 3D Surface		Opens a dialog box that enables you to save the 3D surface of the active data to a file such as Open Inventor format (.iv).
File → Export → 3D Voxels		Opens a dialog box that enables you to save the voxel information from the active data.
File → Export → 3D Scene as DICOM		Opens a dialog box that enables you to save a 3D reconstruction and/or surface in DICOM format. The Multi-Frame DICOM option supports 3D CT reconstruction in third party software.
File → Print		Displays the Print box.
File → Print Preview		Displays the Print Preview box that shows what will be printed.

Table C.1 Menu bar commands and toolbar buttons (continued)



Menu Bar Command	Toolbar Button	Description
File → Recent Files		Shows recently opened data sets. Note: The number of files displayed can be set in the Preferences box (select Edit → Preferences and click the General tab).
File → Logout		Opens the Select/Add User ID dialog box so that another user can logon or a new user ID can be added to the system.
File → Exit		Closes the Living Image software.
Edit → Copy		Copies the active image window to the system clipboard.
Edit → Image Labels		Opens the Edit Image Labels dialog box that enables you to edit the label set information for the active data (see page Figure 4.16 on page 69).
Edit → Preferences		Opens the Preferences box (see page 274).
View → Tool Bar		Choose this option to display the toolbar.
View → Status Bar		Choose this option to display the status bar at the bottom of the main window.
View → Tool Palette		Choose this option to display the Tool Palette.
View → Activity Window		Displays the Activity window at the bottom of the main application window. The Activity window shows a log of the system activity.
View → Image Information		Displays the Image Information box that shows the label set and image acquisition information for the active data.
View → ROI Properties		Displays the ROI Properties dialog box (see page 100).
View → 3D ROI Properties		Displays the 3D ROI Properties dialog box (see page 119).
View → ROI Measurements		Displays the ROI Measurements table.
View → Volume Data Viewer		Enables you to open and view DICOM data.
View → Image Layout Window		Opens the Image Layout window that enables you to paste an image of the active data in the window.
Tools → 3D Animation		Opens the 3D Animation window that enables you to view a preset animation or create an animation.
Tools → Longitudinal Study		Opens the Longitudinal Study window for side-by-side comparisons of DLIT or FLIT results.
Tools → Well Plate Quantification for ...		Opens the Well Plate Quantification window.
Tools → Image Overlay for...		Opens the Image Overlay window for the active data.
Tools → Colorize		Opens the Colorized View tab for the active sequence.
Tools → Image Math for...		Opens the Image Math window for the active data.
Acquisition → Background → Measure Dark Charge		Opens a dialog box that enables you to acquire a dark charge measurement.
Acquisition → Background → Add or Replace Dark Charge		Opens a dialog box that enables you to select an instrument luminescent background. This background measurement is subtracted from luminescent images.

Table C.1 Menu bar commands and toolbar buttons (continued)

Menu Bar Command	Toolbar Button	Description
Acquisition → Background → Measure and Replace Dark Charge		Measures the dark charge under the same conditions as the currently selected image. When the measurement is complete, the newly acquired dark charge image will be included in the dataset of the current image, replacing any existing dark charge image that may be present in the dataset.
Acquisition → Background → View Available Dark Charge		Opens a dialog box that enables you to view the dark charge measurements for the system.
Acquisition → Background → Clear Available Dark Charge		Clears all dark charge images from the system.
Acquisition → Background → Auto Background Setup		Opens a dialog box that enables you to acquire background images, or schedule or disable automatic background acquisition.
Acquisition → Fluorescent Background → Measure Fluorescent Background		Starts a measurement of the instrument fluorescent background.
Acquisition → Fluorescent Background → Add or Replace Fluorescent Background		Opens a dialog box that enables you to select an instrument fluorescent background measurement for the active image data. If the "Fluorescent Background" Subtraction option is chosen in the Corrections/Filtering tool palette, the background measurement is subtracted from the image data.
Acquisition → Fluorescent Background → Measure and Replace Fluorescent Background		Measures fluorescent background under the same conditions as the currently selected image. When the measurement is complete, the newly acquired background image will be included in the data set of the current image, replacing any existing background image that may be present in the data set.
Acquisition → Fluorescent Background → View Available Fluorescent Background		Opens a dialog box that displays the fluorescent background measurements for the system. If a fluorescent background is selected, the "Fluorescent Background Subtraction" option appears in the Corrections/Filtering tool palette. Choose the "Fluorescent Background Subtraction" option to subtract the user-specified background measurement from the image data.
Acquisition → Fluorescent Background → Clear Available Fluorescent Background		Opens a dialog box that enables you to remove the fluorescent background measurements from the system.
Acquisition → Auto-Save		If Auto-Save is selected, all images are automatically saved to a user-selected folder.
Acquisition → CT → Generate Alignment data		Acquires images of the Rotation Stage Alignment tool that are used to generate alignment data for the IVIS® Spectrum CT.
Acquisition → CT → Acquire Reference images		Acquires dark and bright reference images that are used to determine corrections that are applied to the raw projection images during the CT reconstruction process.
Acquisition → Auto-Save To		Opens a dialog box that enables you to select a folder where images will be saved to automatically.
Window → Close		Closes the active image window.
Window → Close All		Closes all image windows.
Window → Cascade		Organizes the open image windows in a cascade arrangement (see Figure 4.7 on page 60).
Window → Tile		Organizes the open image windows in a tiled arrangement (see Figure 4.7 on page 60).

Table C.1 Menu bar commands and toolbar buttons (continued)



Menu Bar Command	Toolbar Button	Description
Window → 1. <Image or Sequence name> Window → 2. <Image or Sequence name> Window → etc.		A list of the open image windows. Click a window in the list to make it the active window (indicated by a check mark).
Window → Other Windows → <window name>		Lists other windows that are open. For example, If the Living Image Browser is open, use these commands to make the browser the active window and display it on top of all other open windows.
Help → User Guide		Displays the Living Image Software Manual.
Help → Tech Notes		Displays a folder of technical notes. Note: For the most recent collection of technical notes, please see the IVIS University download page .
Help → License information		Displays the license information.
Help → Plug-in Information		Displays a list of tool plug-ins and Tool Palette plug-ins.
Help → IVIS Reagents		Opens the PerkinElmer web page for In Vivo Imaging Reagents.
Help → About Living Image		Displays information about the Living Image software and PerkinElmer technical support contact information.
		Click this button, then click an item in the user interface to display information about the item.

Table C.2 Keyboard shortcuts

Keys	Shortcut Description
Ctrl + B	Opens the Living Image Browser.
Ctrl + C	Copies the active image to the system clipboard.
Ctrl + D	Arranges open windows in a cascade.
Ctrl + O	Displays a dialog box that enables you to open data.
Ctrl + P	Open the Print dialog box.
Ctrl + S	Saves the active file or window.
Ctrl + T	Tiles the open windows.
Ctrl + W	Closes the active window.
Shift + F1	Changes the mouse pointer to the “What’s This” tool  . Click this button, then click an item in the user-interface to display information about the item.



NOTE: Macintosh users use the Cmd key (apple key) instead of the Ctrl key.

Index

Symbols

3D Multi-Modality tools

- adjusting image resolution 247
- classifying 3D volumetric data 243
 - color-opacity map 244
 - control points 245
- fiducial registration 257
- gradient illumination 250
- loading data 259–261
- manual registration 263–265
- maximum intensity projection 249
- requirements 242–243
- volume display options 247

3D reconstruction

- checking reconstruction quality 205–207
- DLIT 183
 - inputs 185
 - luminescent sources 192–202
 - sequence
 - acquisition (Imaging Wizard) 186–188
 - requirements (for manual setup) 191
- FLIT 183
 - fluorescent sources 197–200
 - inputs 185
 - sequence
 - acquisition (Imaging Wizard) 188–191
 - requirements (for manual setup) 192
- include/exclude data 200–202
- quick guide 184
- results 203–204
- troubleshooting 234

3D ROI (volumetric data)

- draw 115–119
- measurements 117–119
 - configure the 3D ROI Measurements table 120–122
- properties 119–120

3D tools

- Animate tab 229–234
- Registration 222–228
- Source tab 220–222
- Surface tab 218–219

3D Volumetric Browser 259

A

- adaptive fluorescent background subtraction 65
- animation 229–234
 - custom 231–233
 - edit an animation setup 233–234
 - preset 231

auto exposure 10, 277

autofluorescence (subtract)

- average background ROI 97–99
- image math 125–127
- spectral unmixing 128

B

- background-corrected ROI measurement 97–99
- batch mode (sequence acquisition) 45–46
- binning 62
- browse
 - optical image data 52–56
 - volumetric data 259

C

- cascade images 60
- CCD temperature 8
- center of mass 207–208
- Cherenkov imaging 28
 - DyCE 162–163
 - spectral unmixing 134–136
- classifying 3D volumetric data 243
 - color-opacity map 244
 - control points 245
- colorize data 77–78
- color-opacity map 244
- comments 72
- composite image
 - DyCE 172–173
 - image math 123–125
 - spectral unmixing 139
- control panel 269–272
- corrections tools
 - adaptive fluorescent background subtraction 65
 - binning 62
 - cosmic correction 65
 - dark background subtraction 64
 - flat field correction 64
 - smoothing 62
- cosmic correction 65
- crop box 68

D

- dark background subtraction 64
- Data Preview window 200–202
- DICOM data
 - export 3D scene as DICOM 216–217
 - view (Living Image 3D Browser) 217–218
 - view (Volume Data Viewer) 265–266
- distance measurements 68–69
- DLIT 183, 192–202

- inputs 185
- results 203–204
- sequence
 - acquisition (Imaging Wizard) 186–188
 - requirements (for manual setup) 192
- troubleshooting 234
- DyCE 156–175
 - analysis
 - automatic 164–168
 - manual 168–171
 - imaging
 - bioluminescence 157–159
 - Cherenkov 162–163
 - fluorescence 159–162
 - results
 - composite image 172–173
 - temporal spectra 171, 174–175
 - unmixed images 171–172

E

- edit
 - image label 70
 - sequence 49–50, 81
- export
 - 3D ROI measurements 122
 - images 51, 79–80
 - ROI measurements (optical data) 113
 - surface 182

F

- fiducial registration 257, 263–265
- flat field correction 64
- FLIT 183, 197–200
 - inputs 185
 - results 203–204
 - sequence
 - acquisition (Imaging Wizard) 188–191
 - requirements (for manual setup) 192
 - troubleshooting 234
- fluorescent imaging
 - DyCE 159–162
 - epi-illumination 29–34
 - quick guide 29
 - spectral unmixing 131–134
 - transillumination 35–39
- focus manually 272–273
- FOV settings 25

G

- gradient illumination 250

H

- host organ 208

I

- image

- add comments 72
- adjust appearance 60–61
- cascade 60
- distance measurements 68–69
- export 51, 79–80
- label 26, 27, 69–72
- pan 64
- print 79–80
- tags 73
- tile 60
- zoom 63
- image data
 - colorize 77–78
 - correcting 64–65
 - load (open)
 - from menu bar or toolbar 56–57
 - using the Living Image browser 52–56
- image layout window 79–80
- image math 123–125
 - subtract autofluorescence 125–127
- image overlay tool 74–77
- image window 57
- imaging
 - Cherenkov 28
 - DyCE 162–163
 - spectral unmixing 134–136
 - fluorescent
 - DyCE 159–162
 - epi-illumination 29–34
 - quick guide 29
 - spectral unmixing 131–134
 - transillumination 35–39
 - luminescent 23–28
 - DyCE 157–159
 - quick guide 23
 - spectral unmixing 129–131
 - sequence
 - batch mode 45–46
 - manual 47–50
 - using Imaging Wizard 39–45
- imaging modes 10–11
- Imaging Wizard 39–45
- import
 - organ atlas 227
 - surface 182
- initialize system 7–8

L

- line profile 66
- Living Image 3D browser 217–218
- Living Image browser 52–56
- Living Image software
 - starting 5–7
- load (open) image data
 - from menu bar or toolbar 56–57

- using the Living Image browser 52–56
- Longitudinal Study window 212, 213
- luminescent imaging 23–28
 - DyCE 157–159
 - quick guide 23
 - spectral unmixing 129–131

M

- manual focusing 272–273
- maximum intensity projection 249
- measure
 - distance 68–69
 - source center of mass 207–208
 - source depth 209
 - source intensity 208
 - source volume 208
- measurement ROI (optical data) 89–96
- menu commands 282–285
- mirror ROI 94–96
- Mouse Imaging Shuttle 261
- Multi-Modality Tools. See 3D Multi-Modality Tools.
- multiple reporters per image 74–77

N

- NTF Efficiency 110, 269

O

- open (load) image data
 - from menu or toolbar 56–57
 - using the Living Image browser 52–56
- organ atlas - import 227
- organ display 224–227
- organ registration tools 222–228
- overlying images 74–77

P

- pan an image 64
- passwords 20
- PerkinElmer contact information 4
- photon density 205
- photon density maps 206
 - measured 207
 - simulated 207
- preferences
 - acquisition 277–278
 - general 274–276
 - optical properties 281
 - options 276
 - theme 278–280
- print images 79–80

Q

- quantification database 235–241
 - create 236–240
 - manage results 240
 - samples 235

R

- RAW volumetric data 266–268
- reconstruction. See 3D reconstruction.
- reduced Chi2 203
- registering multi-modal data
 - fiducial registration 257
 - loading data 259–261
 - manual registration 263–265

ROI

- Measurements table 87
- tools 87–88
- ROI (optical data)
 - background-corrected measurement 97–99
 - delete
 - from image 108
 - from system 109
 - draw
 - automatic 90–91
 - free draw 93
 - manual 92–93
 - quick guide 86
 - edit dimensions 104
 - label 107
 - load 108
 - measurements 109–113
 - configure ROI Measurements table 111–112
 - copy or export 113
 - move 103
 - properties 100–103
 - ROI line 106
 - save 107
 - tools 87–89
 - types 83–85
 - average background 84
 - measurement 89–96
 - mirror ROI 94–96
 - subject 97

ROI (volumetric data). See 3D ROI (volumetric data).

S

- save image data (manual) 51
- segment 48
- sequence
 - acquire using Imaging Wizard 39–45
 - batch mode acquisition 45–46
 - create from individual images 82
 - edit 49–50, 81
 - manual setup 47–50
- slices
 - rendering 255
 - viewing 253–254
- smoothing 62
- source
 - center of mass 207–208
 - depth 209

- quantification 208
- volume 208
- spectral unmixing 128–155
 - analysis
 - automatic 142–145
 - guided 137–140
 - library 140–141
 - manual 145–148
 - methods overview 137
 - correcting spectra 149–150
 - imaging
 - bioluminescence 129–131
 - Cherenkov 134–136
 - fluorescence 131–134
 - results 150–155
 - analyzing images 154
 - composite image 139, 153
 - spectra plot 151–152
 - sequence requirements 128
- stage temperature 9
- subject ROI (optical data) 97
- surface
 - export or import 182
 - generate 177–178
 - manage 181
- system activity 22

T

- tags 73
- Tech Notes 3–4
- Technical Support 4
- temperature
 - CCD 8
 - stage 9
- temporal spectra 171, 174–175
- tile images 60
- tool palette 59
- toolbar 282–285
- tools overview 12–19
- transillumination overview 73
- troubleshooting DLIT or FLIT 234

U

- user accounts
 - adding users 19
 - changing or adding passwords 20
 - deleting users 20
 - locking or unlocking accounts 21
- user activity 22
- user preferences 274–281

V

- Volume Data Viewer 265–266
- volume slices
 - information and results 256
- volumetric data

- classify 243, 245
- color-opacity map 244
- display options 247
- information and results 256
- rendering slices 255
- smoothing 252
- viewing slices 253–254
- vsize, starting 203

X

- x,y coordinates 65, 73

Z

- zoom 63