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The screenshot shows the Bioray software interface. At the top, there is a menu bar with 'File' and 'Configuration'. Below it is a toolbar with icons for 'Open Protocol', 'Save', 'Run Protocol', 'Stop', 'Protocol Info', 'Undo', 'Redo', 'Delete', 'Connect', 'Deck In', 'Deck Out', 'Exchange Solo Syringe', and 'Mix LCP'. On the left side, there is a vertical toolbar with icons for 'Loop', 'Pause', 'Aspirate', 'Dispense', 'Mix', 'Exchange', 'Wash', 'Solo Aspirate', and 'Solo Dispense'. The main area is titled 'Protocol Start' and contains a warning message: 'No tray defined Unknown liquid class' and 'Aspirate 0  $\mu$ l'. Below the warning, there is a large red-bordered box representing the 'Aspirate' configuration window. This window has a central area with a purple circle and the letter 'A' next to it, labeled '1'. Below this area are several input fields: 'Volume ( $\mu$ l)' set to 40, 'Air Gap' with 'Volume ( $\mu$ l)' set to 0, 'Pre-dispense' checkbox, 'Number' set to 0, 'Above Bottom (mm)' set to 0, and 'Liquid Class' set to 'Water'. A red-bordered box highlights the 'Volume ( $\mu$ l)' field. A red-bordered box highlights the 'Pre-dispense' checkbox and 'Number' field, with a text box containing the Chinese text: '图示两个位置的参数必须选择。' (The parameters at these two positions must be selected). On the right side of the interface, there is a 'Deck' area with icons for 'Wash Fluoro', 'Intelli-Plate 3-Well', 'Pre Dispense Slide Trimm', and 'Deep Well - 1 ml'. A red-bordered box highlights the 'Deep Well - 1 ml' icon, with a text box containing the Chinese text: '点击 1 号托盘，则 96 针机头会从 1 号托盘取液' (Clicking the 1st tray, the 96-pin head will take liquid from the 1st tray). Below the 'Deck' area, there is another configuration window for the 'Deep Well - 1 ml' tray, with fields for 'Volume ( $\mu$ l)' set to 0, 'Air Gap' with 'Volume ( $\mu$ l)' set to 0, 'Pre-dispense' checkbox, 'Volume ( $\mu$ l)' set to 0, 'Number' set to 0, 'Above Bottom (mm)' set to 0, and 'Liquid Class' set to 'Water'. A red-bordered box highlights this entire configuration window.

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The screenshot displays the Bioray software interface. The 'Protocol Start' list includes the following steps:

- Aspirate 40  $\mu\text{l}$  from Tray 1
- Dispense 20  $\mu\text{l}$  in Tray 2 in location A,1
- Solo Syringe Dispense 0.2  $\mu\text{l}$  drop in Tray 2
- Dispense 0.2  $\mu\text{l}$  in Tray 2 in location A,1
- Empty in Tray 2 in location A,1
- Wash 2 times 45  $\mu\text{l}$  in Tray 3

The 'Wash' step is highlighted. A red box highlights the configuration parameters for this step:

- Number of Cycles: 2
- Volume ( $\mu\text{l}$ ): 45
- Distance Above Bottom (mm): 0
- Liquid Class: Water
- Fill Time (s): 5
- Empty Time (s): 5

A red callout box points to these parameters with the text: 清洗体积要大于吸液体积 5 $\mu\text{L}$ ，一般推荐循环次数为 5

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