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


The screenshot displays the BD FACSDiva software interface with several key components highlighted by red boxes and annotated with Chinese text:

- Browser, 实验管理 (Browser, Experiment Management):** Located on the left side, showing a hierarchical tree view of folders and experiments.
- Cytometer, 参数设置 (Cytometer, Parameter Settings):** A central panel showing a table of parameters for the Cytometer. The table includes columns for Parameter, Voltage, Log, A, H, and W.
- Worksheet, 实验结果展示 (Worksheet, Experiment Results Display):** A large central area displaying a flow cytometry plot (Sheet1) with a histogram showing 'Count' vs 'FSC-A'. The plot is titled 'Specimen_001-Tube_001'.
- Acquisition dashboard, 实验操作控制面板 (Acquisition Dashboard, Experiment Operation Control Panel):** A panel below the Cytometer settings, showing acquisition status, elapsed time, and various control buttons like 'Record' and 'Restart'.
- 侧液流窗口, 调节侧液流偏转角度, dropdelay确定等 (Side Flow Window, Adjust Side Flow Deflection Angle, Drop Delay Confirmation, etc.):** A panel at the bottom left showing various flow parameters and controls.
- Sort Layout, 分选设置 (Sort Layout, Sorting Settings):** A panel at the bottom right showing sorting parameters for 'Tube_001: Sort Layout_012', including 'Sort Rate', 'Confl. Cnt.', and 'Efficiency'.
- 断点窗口 (Breakpoint Window):** A vertical panel on the far right showing a series of droplet images and a 'Sweet Spot' indicator.

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





3.5、测量drop delay (手动)

The screenshot displays the BD FACSDiva software interface with several key components:

- Browser - Accudrop Drop Delay:** A tree view on the left showing the experimental setup hierarchy. The 'Accudrop Drop Delay' folder is highlighted with a red box.
- Cytometer - FACSArial (1):** A central panel showing acquisition parameters. The 'Parameters' tab is active, displaying a table of settings:

Parameter	Voltage	Log	A	H	W
FSC	200	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SSC	300	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
FITC	500	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

- Acquisition Dashboard:** A floating window showing real-time acquisition data for 'Tube_001':
 - Current Activity: Tube_001
 - Threshold Rate: 1292 evt/s
 - Stopping Gate Events: 0 evt
 - Elapsed Time: 00:01:55
- 85 micron:** A floating window for the nozzle. The 'Optical Filter' is highlighted with a red box. It shows a scale from 74 to 78 with a red box around the value 28. Other parameters include Voltage Center (4,500), Plate Voltage (4,500), and Drop Delay (28.26).
- Tube_001: Sort Layout_001:** A floating window for the sort layout. The 'Precision' dropdown is highlighted with a red box and set to 'Fine Tune'. It shows a sort rate of 1292 evt/s and various performance metrics.
- Normal Worksheet - Sheet1:** A plot window showing a histogram of 'Specimen_001-Tube_001' with a peak labeled 'P1'.

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3.5、测量drop delay (自动)

The screenshot displays the BD FACSDiva software interface with several key components:

- Browser - Accudrop Drop Delay:** A tree view on the left showing a folder structure. The 'Accudrop Drop Delay' folder is highlighted with a red box.
- Cytometer - FACSArial (1):** A table showing parameters for the instrument. The 'Parameters' tab is active, displaying a table with columns for Parameter, Voltage, Log, A, H, and W.
- Normal Worksheet - Sheet1:** A plot titled 'Specimen_001-Tube_001' showing a histogram of 'Count' vs 'FSC-A (x 1,000)'. A gate labeled 'P1' is visible.
- Acquisition Dashboard:** A central panel showing 'Current Activity' with 'Active Tube/Well' set to 'Tube_001', 'Threshold Rate' at '1292 evt/s', 'Stopping Gate Events' at '0 evt', and 'Elapsed Time' at '00:01:55'. It includes 'Basic Controls' (Next Tube, Unload, Stop Acquir..., Record Data, Restart) and 'Acquisition Setup' (Stopping Gate, Storage Gate, Events To Record, Events To Display, Stopping Time, Flow Rate).
- 85 micron:** A control panel at the bottom showing 'Drop Delay' set to '28.26' with an 'Auto Delay' button highlighted in red. It also displays 'Voltage Center' and 'Plate Voltage'.
- Tube_001: Sort Layout_001:** A panel on the right showing 'Sort Rate' at '1292 evt/s' and 'Efficiency' at '100%'. It includes 'Sort', 'Pause', and 'View Counters' buttons.



Parameter	Voltage	Log	A	H	W
FSC	200	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SSC	300	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
FITC	500	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Sort Rate	Conf. Cnt	Conf. Rate	Efficiency
1292 evt/s	0 evt	0 evt/s	100%

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Global Sheet1: Sort Layout_001

Device: 2 Tube Precision: Purity Target Events: Continuous Save Sort Reports: Ask User Save Conflicts: Index Sorting:

	Left	Right
Sort Rate:	NA	13 evt/s
Confl.Cnt:	NA	13876 evt
Confl. Rate:	NA	5 evt/s
Efficiency:	NA	71%

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