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

Initiating Instrument Connection...

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About

MicroCal PEAQ-ITC Control Software

Malvern MicroCal PEAQ-ITC Control Software

Run Experiment Maintenance Design Experiment

Start Experiment Start Experiment Load Run Clean

Methods Experiments

Browse... C:\ > Users > Public > Documents > Malvern Instruments > MicroCal PEAQ-ITC > Methods

Name	Modified Date
13 Injections.itcm (MicroCal Method)	
19 Injections.itcm (MicroCal Method)	

Open

Preview - 13 Injections.itcm

This experimental method will make a single 0.4 μL injection, followed by 12, 3 μL injections.

Temperature ($^{\circ}\text{C}$)	25.0
Reference Power (μW)	41.9
Feedback	High
Stir Speed (rpm)	750
Initial Delay (s)	60.0
Injection Spacing (s)	150
Injection Duration (s)	6

Instrument Offline

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The screenshot shows the MicroCal PEAQ-ITC Control Software interface. A red box highlights the main navigation bar at the top, which includes the following elements:

- Buttons: Run Experiment (highlighted), Maintenance, Design Experiment
- System icons: Malvern logo, a red warning icon with the number '2', and HELP/ABOUT links.
- Settings: A gear icon and a zoom control set to 100%.

Below the navigation bar, the interface is divided into several sections:

- Run Section:** Contains 'Start Experiment', 'Load', 'Run' (highlighted), and 'Clean' buttons.
- Experiment Parameters:** Includes fields for [Syr] (M), [Cell] (M) (set to 0), and a Comment text area.
- Instrument Settings:** A list of parameters including Temperature (°C) set to 25.0, Reference Power (10.0), Feedback (High), Stir Speed (750), Initial Delay (60), and Injection Spacing (150).
- Process Flow:** A horizontal timeline with four stages: Idle (active), Setting Temperature, Equilibrating, and Injecting (0/13).
- Graph:** A plot of DP (μcal/s) vs Time (min) with a y-axis from -0.04 to 0.04 and an x-axis from -0.4 to 0.5.
- Bottom Bar:** Shows 'Instrument Offline' and buttons for 'Start', 'Analyze', and 'Save As Method'.

Callouts with arrows point to specific features:

- '实验参数设置' points to the 'Run Experiment' button.
- '实验设计' points to the 'Design Experiment' button.
- '设置' points to the gear icon in the navigation bar.
- '仪器维护视频' points to the 'Start Experiment' button.
- '字体大小设置' points to the zoom control in the navigation bar.

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Instrument Settings

Temperature (°C) 25.0

This experimental method will make a single 0.4 μL injection, followed by 12, 3 μL injections.

The following settings will be used:

Reference Power ($\mu\text{cal/s}$)	10.0
Feedback	High
Stir Speed (rpm)	750
Initial Delay (s)	60
Injection Spacing (s)	150

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Injection	Volume (μL)	Duration (s)	Spacing (s)
1	0.4	0.8	150
2	3.0	6.0	150
3	3.0	6.0	150
4	3.0	6.0	150
5	3.0	6.0	150
6	3.0	6.0	150
7	3.0	6.0	150
8	3.0	6.0	150

DP ($\mu\text{cal/s}$)

Graph showing DP ($\mu\text{cal/s}$) vs. Injection number.

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4003384toTD2-2.itc - MicroCal PEAQ-ITC Control Software


Malvern MicroCal PEAQ-ITC Control Software

Run Experiment Maintenance Design Experiment

Clean Start Experiment Load Run Clean

0 1 2 3 4 5 6

0 Introduction

 To view a video depicting all the cleaning steps, click the Play button.
To enter the Clean Instrument workflow, click Next.

Back Next

- 1 Choose Cleaning Method(s)
- 2 Insert Cell Cleaning Tool
- 3 Attach Fill Port Adapter
- 4 Move Pipette to Clean Location
- 5 Detach Fill Port Adapter
- 6 Remove Cell Cleaning Tool

Show entire video in external player

Instrument Online DP (μcal/s): 4.95 Temperature (°C): 22.7

Stop Idle Temperature (°C): 25 Set Syringe: Load Plunger Down Cell Clean: None Syringe Clean: None Clean

Windows taskbar: ITC 4003384toTD2-2.itc ... Untitled* - Paint 3D 10:36 AM 4/5/2020

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MicroCal PEAQ-ITC Control Software

Run Experiment Maintenance Design Experiment

Run Start Experiment Load Run Clean

Experiment Information

[Syr] (M) 500e-6
[Cell] (M) 33.0e-6
Comment

Instrument Settings

Temperature (°C) 25.0

This experimental method will make a single 1 μ l injection, followed by 18, 2 μ l injections.

The following settings will be used:

Reference Power (μ cal/s) 5.00
Feedback High
Stir Speed (rpm) 750
Initial Delay (s) 60
Injection Spacing (s) 120
Injection Duration (s) 4.0

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
New Analysis - MicroCal PEAQ-ITC Analysis Software

Malvern MicroCal PEAQ-ITC Analysis Software

Analyze Experiment(s) Design Experiment

Assign Controls

Start Analysis Overview Adjust Baseline **Assign Controls** Adjust Fit Presentation

Experiments

Sort by Bin

CaEDTAgGetStart *
Binding >

CaEDTAgGetStart_ctrl *
Control

Control Parameters

Type Single

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Titrant -> Buffer (-) CaEDTAgGetStart_ctrl Method Line

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
Reset Show Excluded Injections Show Legends

Experiment Information


Filename	CaEDTAgGetStart
Temperature (°C)	25.0
[Syr] (M)	1.00e-3
[Cell] (M)	100e-6
Ref. Power (μcal/s)	5.00 (5.35)
Comment	
Control Type	Single
Control Experiment	CaEDTAgGetStart_ctrl
Results	
Bin Comment	Binding
Competitive Model	No
[Syr] (M)	1.00e-3
[Cell] (M)	100e-6
Ligand In Cell	No
N (sites)	0.969 ± 1.5e-3

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New Analysis - MicroCal PEAQ-ITC Analysis Software

Malvern MicroCal PEAQ-ITC Analysis Software

Analyze Experiment(s) Design Experiment

Adjust Fit

Start Analysis Overview Adjust Baseline Assign Controls **Adjust Fit** Presentation

Experiments

Sort by Bin

CaEDTAGetStart *
Binding >
CaEDTAGetStart_ctrl *
Control

CaEDTAGetStart

Experiment Information

Filename: CaEDTAGetStart
Temperature (°C): 25.0
[Syr] (M): 1.00e-3
[Cell] (M): 100e-6
Ref. Power (μcal/s): 5.00 (5.35)
Comment:

Control Type: Single
Control Experiment: CaEDTAGetStart_ctrl

Results

Bin Comment: Binding
Competitive Model: No
[Syr] (M): 1.00e-3
[Cell] (M): 100e-6
Ligand In Cell: No
N (sites): 0.969 ± 1.5e-3

Fitting Parameters

Ligand is in Cell

Parameter	Vary	Initial Value	Lower Bound	Upper Bound
N (sites)	<input checked="" type="checkbox"/>	0.969	1.0e-3	10.0
Offset (kcal/mol)	<input type="checkbox"/>	0	-80.0	80.0
[Syr] (M)	<input type="checkbox"/>	1.00e-3	0	1.00

Use Competitive Model

Unknown Binder: Strong

[Weak] (M): 100e-6

Known Weak Parameters: Enter Manually

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New Analysis - MicroCal PEAQ-ITC Analysis Software

Malvern MicroCal PEAQ-ITC Analysis Software

Analyze Experiment(s) Design Experiment

Presentation

Start Analysis Overview Adjust Baseline Assign Controls Adjust Fit Presentation

Result Table Final Figure Scatter Plot Injection Table Statistics Plot Signature Plot Raw Plot Integrated Heat Plot

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



29.5 μ M Protein titrated with
1.1mM Compound

11.5 μ M Compound titrated with
179 μ M Protein


Parameter	Ligand in syringe	Ligand in cell
n	0.99	0.97
K_d	104 nM	105 nM
ΔG°	-9.4 kcal/mol	-9.4 kcal/mol
ΔH°	-11.9 kcal/mol	-12.4 kcal/mol
$T\Delta S^\circ$	-2.5 kcal/mol	-3.0 kcal/mol

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
Peptide binding to Batch 1 Protein

N = 1.05
K_D = 97 nM

100% of Batch 1
protein Active based
on Stoichiometry

Peptide binding to Batch 2 Protein

50 μM Peptide
10 μM Protein X

N = 0.235
K_D = 135 nM


Re-analyzed as
2.3 μM Protein X

N = 1.02
K_D = 135 nM

23% of Batch 2
protein Active based
on Stoichiometry

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